



**US
DENTAL
DEPOT**

US DENTAL DEPOT

6555 NW 9th Avenue suite 212
Fort Lauderdale 33309
Florida USA

TEL 1 1866 625 5558 – 1 954 874 6325
www.USDENTALDEPOT.com

ART-M3 II (UPGRADED SERIES) Magnetostrictive Ultrasonic Scaler

NEW MODEL!! An improved features and functions of the ART-M3 Magnetostrictive Ultrasonic Scaler, this system has our unique "AGC" (Auto Gain Control), and an LCD digital display screen that shows scaling mode and power setting.



ART-M3 II (UPGRADED SERIES) Magnetostrictive Ultrasonic Scaler

System is available in either 25KHz or 30KHz stand-alone system. Up to a maximum power setting of 50 with auto tune system that locks in the **frequency and stabilize power setting**.

The well-equipped feature called **"AGC" (Auto Gain Control)** stabilizes the power setting and prevents stalling during scaling operation. System comes with a **built in LCD digital screen** that displays the mode and power setting, and allows the user to choose from **Scaling, Perio and Endo mode** with one turn of the knob, and an easy to access **water control knob at the front control panel**. System includes 1 plastic grip internal flow insert and 1 metal grip external flow insert.

Available Now! Call us for your authorized dealer.

ISO certified; CE marked; 510K registered.

5-year limited warranty on the internal electronic module

for this model only.

WWW.USDENTALDEPOT.COM