



# US DENTAL DEPOT

## US DENTAL DEPOT

6555 NW 9<sup>th</sup> Avenue suite 212  
Fort Lauderdale 33309  
Florida USA

TEL 1 866 625 5558 – 1 954 874 6325  
[www.USDENTALDEPOT.com](http://www.USDENTALDEPOT.com)

### ART-P3 II (UPGRADED SERIES) Piezo Electric Scaler

NEW MODEL!! An improved features and functions of the ART-P3 Piezelectric Scaler, this system features an "AGC" (Auto Gain Control) capability, and a built in LCD digital display screen that shows scaling mode and power setting.



### ART-P3 II (UPGRADED SERIES) Piezo Electric Scaler

An improved model that almost resembles the look of the old ART-P3 is equipped with more features and functions. Up to a **maximum power setting of 50** with auto tune system that locks the frequency. The well-equipped feature called **"AGC" (Auto Gain Control)** stabilizes the power setting and prevents stalling during scaling operation. System comes with **a built in LCD digital screen** that displays the mode and power setting, allows the user to choose from **Scaling, Perio and Endo mode** with one turn of the knob, and has an easy to access **water control knob at the front control panel**. System includes 4 various tips (2 pointed tips & 2 flat tips), 4 individual plastic autoclavable wrenches, 1 aluminum sterilization box.

Available Now!! Call us for your authorized dealer.

**ISO certified; CE marked; 510k registered.**

**5-year limited warranty on the internal electronic module;  
and 2 year warranty on the hand piece. (This Model Only)**

[WWW.USDENTALALDEPOT.COM](http://WWW.USDENTALALDEPOT.COM)