



**US
DENTAL
DEPOT**

US DENTAL DEPOT

6555 NW 9th Avenue suite 212
Fort Lauderdale 33309
Florida USA

TEL 1 1866 625 5558 – 1 954 874 6325
www.USDENTALDEPOT.com

ART-P4 Piezo Electric Scaler w/ AGC Feature (Blue)

NEW MODEL!! Bonart's newly designed ART-P4 defines the art of modern innovative dental scaling technology. The unique features of this system include the "AGC" Auto Gain Control, capability and a built in LCD digital display.



ART-P4 Piezo Electric Scaler w/ AGC Feature (Blue)

Up to a maximum power setting of 50 with auto tune system that locks the frequency. The well-equipped feature called **"AGC" (Auto Gain Control)** stabilizes the power setting and prevents stalling during scaling operation. System comes with **a built in LCD digital screen** that displays the mode and power setting, allows the user to choose from **Scaling, Perio and Endo mode** with one turn of the knob, and has an easy to access **water control knob at the front end of the hand piece cable.** System includes 4 various tips (2 pointed tips & 2 flat tips), 4 individual plastic autoclavable wrench, 1 aluminum sterilization box.

Available Now!! Call us for your authorized dealer.

ISO certified; CE marked; 510k registered!

5-year limited warranty on the internal electronic module;

and 2 year warranty on the hand piece. (This Model Only)

WWW.USDENTALALDEPOT.COM