



US DENTAL DEPOT

US DENTAL DEPOT

6555 NW 9th Avenue suite 212
Fort Lauderdale 33309
Florida USA

TEL 1 1866 625 5558 – 1 954 874 6325
www.USDENTALDEPOT.com

ART-M4 Magnetostrictive Ultrasonic Scaler w/ AGC Feature (Blue)

NEW MODEL!! Introducing its newest model of magnet scaler system-ART-M4. The new system brings out the true definition of modern innovative technology designed by Bonart's finest engineering team.



ART-M4 Magnetostrictive Ultrasonic Scaler w/ AGC Feature (Blue)

System is available in either 25KHz or 30KHz stand-alone system. Up to a maximum power setting of 50 with auto tune system that locks the frequency and stabilize power setting

The well-equipped feature called **"AGC" (Auto Gain Control)** stabilizes the power setting and prevents stalling during scaling operation. System comes with **a built in LCD digital screen** that displays the mode and power setting, and allows the user to choose from **Scaling, Perio and Endo mode** with one turn of the knob, and an easy to access **water control knob at the front end of the hand piece cable**. System includes 1 plastic grip internal flow insert and 1 grip handle external flow insert.

Available Now! Call us for your authorized dealer.

Note: Blue is the only color available for the , Canadian,

and South and Central American customers.

ISO certified; CE marked; 510K registered.

**5-year limited warranty on the internal electronic module
for this model only.**

WWW.USDENTALALDEPOT.COM