



US DENTAL DEPOT

US DENTAL DEPOT

6555 NW 9th Avenue suite 212
Fort Lauderdale 33309
Florida USA

TEL 1 1866 625 5558 – 1 954 874 6325
www.USDENTALDEPOT.com

ART-E1 Electrosurgery / Cutting Unit With 7 Pc. Electrodes

State of the art radio frequency electrosurgical unit proudly manufactured and offered by Bonart. Unit is comparable to major brands. Great quality & economical. One of the best selling and most popular products of Bonart.



ART-E1 Electrosurgery / Cutting Unit With 7 Pc. Electrodes

This is an electrosurgery system and not a cautery unit.

BONART ART-E1 Electrosurgery/cutting unit comes with 7 piece variety electrodes set, 1 hand piece, 1 indifferent or grounding plate, 1 power cord, AND 1 foot switch. Operation mode of cut, coagulation, cut-coagulation. Operation frequency is at 1.7 MHZ. 10 different power settings. Unit is a monopolar mode but bipolar converter is optional and available upon request. Electrodes are autoclavable, 1/16" diameter shaft. Made universal to other major electrosurgical brands. Ideal for human (dental and dermatology use) and veterinary field.

ISO Certified; CE Certified; 510k registered