



US DENTAL DEPOT

US DENTAL DEPOT

6555 NW 9th Avenue suite 212
Fort Lauderdale 33309
Florida USA

TEL 1 1866 625 5558 – 1 954 874 6325
www.USDENTALDEPOT.com

ART-SP1 Magnetostrictive Scaler Polisher Combo Unit

Magnetostrictive type combo unit provides the ease of use in scaling performance and polishing function. No need for 2 individual units to perform each operation.



ART-SP1 Magnetostrictive Scaler Polisher Combo Unit

With the simple press of the button, BONART ART-SP1 Magnetostrictive scaler and polisher combo unit can function in either the scaling or polishing mode. This model uses magnetostrictive technology for scaling, and a low speed system for polishing. This 2 and 1 system is available in both the 115V and 235volts. Output power for this unit ranges from 10W TO 30W with working frequency of 24.5KHZ while water supply capacity is at 25-60 PSIG (172-414KPa). Polisher speed is from 2,000-30,000 RPM. INCLUDED ACCESSORIES ARE E-TYPE LOW SPEED MICROMOTOR WITH STAND, STRAIGHT HAND PIECE ATTACHMENT, STAINLESS STEEL PROPHY ANGLE(SNAP-ON OR SCREW ON) AND 25K INSERT. PROPHY CUPS AND PASTE ARE NOT INCLUDED BUT CAN BE PURCHASED FROM ANY DENTAL SUPPLY DEALERS OR DISTRIBUTORS.

UNIT IS FOR VETERINARY USE ONLY. NOTE:

Accessories included with the unit might vary among regions. All optional accessories can be purchased at

additional costs. NOTE: Prophy cups and paste not included.

www.usdentaldepot.com