



**US
DENTAL
DEPOT**

US DENTAL DEPOT

6555 NW 9th Avenue suite 212
Fort Lauderdale 33309
Florida USA

TEL 1 1866 625 5558 – 1 954 874 6325
www.USDENTALDEPOT.com

ART-SP2 Piezo Electric Scaler Polisher Combo Unit

With the simple press of the button, BONART ART-SP2 Piezoelectric scaler and polisher combo unit can function in either the scaling or polishing mode. This model uses piezo technology for scaling, and a low speed system for polishing. This 2 and 1 system is available in both the 115V and 235volts. Output power for this unit ranges from 10W TO 30W with working frequency of 29.5-32KHZ while water supply capacity is at 25-60 PSIG (172-414KPa). Polisher speed is from 2,000-30,000 RPM. INCLUDED ACCESSORIES ARE E-TYPE LOW SPEED MICROMOTOR WITH STAND, STRAIGHT HAND PIECE ATTACHMENT, STAINLESS STEEL PROPHY ANGLE(SNAP-ON OR SCREW ON), 3 VARIOUS TIPS (1 POINTED & 2 FLATS) AND 1 PLASTIC AUTOCLAVABLE WRENCH. PROPHY CUPS AND PASTE ARE NOT INCLUDED BUT CAN BE PURCHASED FROM ANY DENTAL SUPPLY DEALERS OR DISTRIBUTORS.

UNIT IS FOR VETERINARY USE ONLY. NOTE:
Accessories included with the unit might vary among regions. All optional accessories can be purchased at additional costs. NOTE: Prophy cups and paste not included.