



**US
DENTAL
DEPOT**

US DENTAL DEPOT

6555 NW 9th Avenue suite 212
Fort Lauderdale 33309
Florida USA

TEL 1 1866 625 5558 – 1 954 874 6325
www.USDENTALDEPOT.com

ART-L2-B Blue Halogen Lamp With Built In Radiometer Curing Light

Bonart introduced its conventional halogen lamp curing light in 2002 to meet its customers' demand and usability needs.



ART-L2-B Blue Halogen Lamp With Built In Radiometer Curing Light

This system features an accurate **built-in radiometer**, 1 set of the **9.5MM light probe**, and **eye protector**. It uses an original **Philips 75 WATT** halogen light bulb.

Built with 20sec, 40sec, 60sec and 120 sec internal timer and beeper. High output light intensity within the range of **500-550mW/cm²**. It also features a **10 SECOND BEEPER REMINDER**.

AVAILABLE IN EITHER 110 OR 220 VOLTS WITH A DIMENSION OF 25CM X 12.5CM X21.5CM.

ISO Certified; 510k registered