



**US
DENTAL
DEPOT**

US DENTAL DEPOT

6555 NW 9th Avenue suite 212

Fort Lauderdale 33309


Florida USA

TEL 1 1866 625 5558 – 1 954 874 6325

www.USDENTALDEPOT.com

BONART (P-10) Series Straight Angle Ultrasonic Inserts

Quality crafted inserts comparable to major brands. An external flow metal handle insert built to outlast its competitors. U.S. made materials.



BONART (P-10) Series Straight Angle Ultrasonic Inserts

P-10 series ultrasonic inserts are available in both 25khz and 30khz. Universal design fits major brands' scalers like Dentsply, Henry Schein, Parkell, Autoscaler, scalers etc. This insert is made of well-crafted stainless steel tips, and nickel strip stem. External water flow.

510k Registered