



**US
DENTAL
DEPOT**

US DENTAL DEPOT

6555 NW 9th Avenue suite 212
Fort Lauderdale 33309
Florida USA

TEL 1 1866 625 5558 – 1 954 874 6325
www.USDENTALDEPOT.com



BONART (IF-50) Streamline #10 Ultrasonic Inserts-30KHZ

Quality 30khz internal flow plastic handle insert designed to be universally compatible with the major brands. Inserts are color coded for khz classifications. Compatible with major brands like Dentsply Cavitron, Parkell, Hu-Friedy, AutoScaler among the few.

BONART (IF-50) Streamline #10 Ultrasonic Inserts-30KHZ

IF50 Streamline #10 Ultrasonic Insert-30khz is designed with a plastic handle for easier and stable grip, and water flows from the tip internally. Made to be comparable with the Parkell, Dentsply, Hu-Friedy, Premier etc. 30khz insert is color coded in orange.

FDA 510K Registered!!