



**US
DENTAL
DEPOT**

US DENTAL DEPOT

6555 NW 9th Avenue suite 212

Fort Lauderdale 33309

Florida USA

TEL 1 1866 625 5558 – 1 954 874 6325

www.USDENTALDEPOT.com



BONART(IF5) Beavertail #3 Ultrasonic Inserts-25KHZ

Quality 25khz internal flow plastic handle insert is made to be universally compatible with major brands. Quality made product to be cost economical. U.S. made materials.

BONART(IF5) Beavertail #3 Ultrasonic Inserts-25KHZ

IF5 Beavertail Ultrasonic Insert-25khz is a flat wide spatula insert designed to get around the hard calculus in the hidden mouth areas more easily. Compatible with the markets' leading brands. 25khz insert is color coded in yellow.

FDA 510k Registered!!

www.USDENTALDEPOT.com