



**US
DENTAL
DEPOT**

US DENTAL DEPOT

6555 NW 9th Avenue suite 212
Fort Lauderdale 33309
Florida USA

TEL 1 1866 625 5558 – 1 954 874 6325
www.USDENTALDEPOT.com

BONART (IF-5) Beavertail #3 Ultrasonic Inserts-30KHZ

A 30khz internal flow plastic handle beavertail insert designed to be universally compatible with the major brands. Inserts are color coded for khz classifications.



BONART (IF-5) Beavertail #3 Ultrasonic Inserts-30KHZ

BONART (IF-5) Beavertail #3 Ultrasonic Inserts-30KHZ

IF5 Beavertail Ultrasonic Insert-30khz is a flat wide spatula insert to get around the calculus in the hidden areas more easily. Compatible with the markets' leading brands. 30khz insert is color coded in orange.

FDA 510K Registered!!