



**US
DENTAL
DEPOT**

US DENTAL DEPOT

6555 NW 9th Avenue suite 212
Fort Lauderdale 33309
Florida USA

TEL 1 1866 625 5558 – 1 954 874 6325
www.USDENTALDEPOT.com

ART-M1 Magnetostrictive Ultrasonic Scaler

The most basic Bonart Magnet Scaler system, this system offers powerful easy to use scaling performance that outlasts its competitors.



ART-M1 Magnetostrictive Ultrasonic Scaler

This system gives you excellent quality at an affordable price.

Specifications are as follows:

Operation Frequency: AVAILABLE IN EITHER A 25KHZ OR 30 KHZ.

Output Power: From 10 to 30 Watts of output power during normal scaling speed.

Turbo Function: Increase power during scaling with one touch of a button. C25W to 35W output power.

Water Supply: 25-60 PSIG (172-414 Kpa)

Optional accessories: 25KHz or 30KHz inserts

ISO Certified; CE Marked; 510K registered.

18 month limited warranty on the internal electronic module

WWW.USDENTALDEPOT.COM