

EQUIVALENCIAS ALEACIONES SEMIPRECIOSAS

Category name	ARGEN	JELENKO	DENTSPLY	GENERIC
High Gold for Ceramics	ARGEDENT 87	JEL O		
Palladium Base for Ceramics	ARGELITE 76SF+	LIBERTY	Option	
Reduced Gold for Crown and Bridg	ARGESOL YSF			
Reduced Gold for Crown and Bridg	ARGESOL LO			
Palladium Base for Crown and Bri	ARGENCO Y+		Ney-Oro 41	
Palladium Base for Crown and Bri	ARGENCO W+	SYMPHONY	Ney-Oro 66 CF	
High Gold for Crown and Bridge	ARGENCO 50	SUNCAST	JEL BIOS 99+	
Reduced Gold for Ceramics	ARGEDENT 62		Miracast	
Reduced Gold for Ceramics	ARGEDENT 500		Ney-Oro 60	
High Gold for Ceramics	ARGEDENT 90	JELBIOS 90	Bioney	
High Gold for Ceramics	ARGEDENT 75	JELENKO #75		
High Gold for Crown and Bridge	ARGENCO 83	BENCHMARK I		
High Gold for Crown and Bridge	ARGENCO 77	MODULAY		
High Gold for Crown and Bridge	ARGENCO 68	JEL 7	Ney-Oro G-3 Y Primalor G	
High Gold for Crown and Bridge	ARGENCO 62	JEL 62		
High Gold for Crown and Bridge	ARGENCO 56	LAB 44 LABORATORY 44	Argedent Yellow 2	
High Gold for Crown and Bridge	ARGENCO 52 HN	SUNCAST D	Argedent Yellow 4	
High Gold for Crown and Bridge	ARGENCO 5	JEL 5		
Reduced Gold for Crown and Bridg	ARGESOL WSF		Degulor C	
High Gold for Ceramics	ARGEDENT BIO 88			
High Gold for Ceramics	ARGEDENT BIO 89PF	DIAMOND		
High Gold for Ceramics	ARGEDENT YELLOW 2		Degudent GS	
High Gold for Ceramics	ARGEDENT YELLOW 4		Degudent H	
Reduced Gold for Ceramics	ARGEDENT 52	CAMEO	Stabilor H	
Reduced Gold for Ceramics	ARGENCO 40 HN	JEL 40	Stabilor G	
High Gold for Crown and Bridge	ARGENCO BIO 1			
Palladium Base for Crown and Bri	ARGEDENT 52SF (SPECIAL)	OLYMPIA	Ney-Oro #5	
Palladium Base for Crown and Bri	ARGEDENT 52SF			
Palladium Base for Crown and Bri	ARGEDENT 52SF			
High Gold for Ceramics	GOLDTECH BIO 2000	JELBIOS 99+		
Multi-Purpose	ARGISTAR BIO 75PF	ALLROUND 1		
Multi-Purpose	ARGISTAR BIO 74			
Multi-Purpose	ARGISTAR BIO 72	ALLROUND 4		
Multi-Purpose	ARGISTAR 73			

Reduced Gold for Crown and Bridg	ARGESOL 650		Ney-Oro B-2	
Reduced Gold for Ceramics	ARGENCO 46	MIDAS	Ney-Oro CB	
Reduced Gold for Ceramics	ARGENCO 4		Densilay	
Reduced Gold for Ceramics	ARGENCO 20	PROSPECTOR	Ney-Oro B-20	
Solders			B-Inlay	
High Gold for Ceramics	ARGEDENT BIO 88PF		Degudent G	
High Gold for Ceramics	ARGEDENT BIO YELLOW PF			
Solders			Ney-Oro-A-1	
High Gold for Ceramics	ARGEDENT YELLOW SPECIAL		BiOcclus 4	
Solders			Ney 76	
Palladium Base for Ceramics	ARGELITE 55		Applause	
Palladium Base for Ceramics	ARGELITE 75+6	ACCU-STAR		
Solders			Palliag M	
Palladium Base for Ceramics	ARGELITE 76+	LEGACY XT	Degubond Ultra	
Palladium Base for Ceramics	ARGELITE 80SF	FREEDOM		
Palladium Base for Ceramics	ARGELITE 80SF+			
Palladium Base for Ceramics	ARGELITE A			
Palladium Base for Crown and Bri	ARGEDENT 52SF	OLYMPIA	Duallor R	
Palladium Base for Ceramics	ARGELITE PF+		Ultima Gold	
High Gold for Ceramics	Argedent Bio 78PF	ALLAROUND 1	Excel	DEGUNOM
High Gold for Ceramics	Argedent 40		Veritas	
High Gold for Ceramics	Argedent 45			
High Gold for Ceramics	Argedent 45SF		Encore	
High Gold for Ceramics	Argedent 52		Ceramco White	
High Gold for Ceramics	Argedent Yellow 4		Image 2	
Palladium Base for Ceramics	Argelite 78SF		Ultima Lite	
	Argelite 58		Pors-on 4	
	Argelite 55			
	Argelite 60		Degubond 60	
	Argelite 60		Lunar	
	Argelite 61		Pors-on Lite	
	Argelite 70+		Ovation	
	Argelite 49		Acclaim	
	Argelite 50		Acclaim Lite	
	Argelite 50		Degustar F	
	Argelite 52+		Pors-on Plus	
	Argelite 54		Pd-Pro	

	Argelite 54+		Degustar GA-2	
	ARGENCO W-LITE			
	ARGENCO W	ALBACAST		
	ARGENCO BIO 1			
	ARGENCO WR+	MAESTRO	Quantum 2	
	ARGEDENT 65SF	JEL SPAN	G-Cast	
	ARGEDENT 54	COLLEGIATE	C Special	
			Ney-Oro LGY	
			Degudent LS	
			Degunorm	
			Multi-PF	
			Majority	
			Ney-Oro PF-CB	
			Stabilor LS	
			Argedent 3	
			Argedent 3	
	Argedent 690		SMG-CF	
	Argedent 70SF		Degudent SF	