

Comparing Flexite Plus to Valplast

Flexite

Valplast

Color Stability, Retains Original Color Pigment	Yes	No
Stains Easily Removed By Pumicing	Yes	No
No Preheating of Plastic in Cartridges	Yes	No
FDA Grade Plastics	Yes	No
Easy to Add to and Reline	Yes	No
Unbreakable Plastic	Yes	Yes
Finishing and Polishing	Easier to Finish and Polish	More Difficult to Finish and Polish
Pink Shades Available	Light Pink, Dark Pink, and Translucent Pink	Light Pink and Standard Pink
Ethnic Shades Available	Ethnic Brown, Ethnic Purple, Coe Ethnic Moderate	Standard Meharry and Light Meharry

Flexite Certified Labs must take a two day course of instruction which emphasizes surveying and designing
 Flexite rarely uses gum tissue for retention and explores all types of clasps to maximize esthetics.
 Our ultimate goal is to minimize unnecessary chair time spent with adjustments by the dentist.

Flexite Plus and Valplast plastics are nylons and will stain over time. Valplast will change color while Flexite Plus can be pumiced back to its original color.

Kleenite is our choice for cleaning and maintenance. Kleenite can be purchased in most drug stores.