

# EURO WORLDENT



*Spin Cycle*  
SLIM LINE

# Spin Cycle

SLIM LINE

## Wall Mount Vacuum Mixing System

- **Produces ultra-smooth mixing for all investments and gypsums**  
*Three dimensional mixing is all done simultaneously*
  1. Horizontal mixing
  2. Vertical mixing
  3. Diagonal mixing
- **Heavy duty mixing bowls and "Inox" stainless steel spatulas**  
*Shows Max Fill Line for liquid + powder*
- **"Hands Free" operation, vacuum automatically holds bowl in place**  
*Increases productivity and saves time*
- **Includes built-in automatic programmable digital timer with "Auto-Stop"**  
*Select the amount of time and Spin Cycle will stop mixing automatically*
- **Eliminates 100% of all vacuum pump maintenance**  
*No oil or oil vapor contamination and less moving parts*
- **Built-in air pressure regulator provides optimum vacuum levels**  
*Automatically converts 50-55 psi air pressure to "optimum vacuum levels"*
- **Overflow Filter Protection System**
- **Includes stainless steel wall mounting plate**  
*Slim-line space saving design*



### *Spin Cycle* SLIM LINE, p/n 80SC

*Includes 1 ea. "Large" mixing bowl (500 ml) with Stainless Steel Spatula,  
Built-in Programmable Digital Timer with Audio Tone,  
Auto Mixing Stop Feature,  
Stainless Steel Wall Mounting Plate  
Optional Bowls/Spatulas Available: small 200 ml and XL 900 ml*

For additional information contact an authorized dealer or Euro-Worldent, Inc.



9600 W. Sample Road, Suite 403  
Coral Springs, FL 33065  
Tel. (561) 470-9287  
e-mail: [WORLDENT1@aol.com](mailto:WORLDENT1@aol.com)  
website: [www.worldentgroup.com](http://www.worldentgroup.com)