



**US
DENTAL
DEPOT**

US DENTAL DEPOT

6555 NW 9th Avenue suite 212
Fort Lauderdale 33309
Florida USA

USA - 1 1866 625 5558 – 1 954 874 6325

Colombia 57 1 489 9051

Peru 51 1 707 3444

Venezuela 58 212 729 1639

Areentina 54 11 5258 7148 Cba 54 351 569 1409

www.USDENTALDEPOT.com

G.F.M.

TEMPERATURE CHART
INSTRUCTIONS FOR INJECTION

G.F.M.

V- plas

530 Fareng

270 Celius

20 Minutes

Clear

545 Fareng

280 Celius

20 Minutes