



**US
DENTAL
DEPOT**

US DENTAL DEPOT

6555 NW 9th Avenue suite 212
Fort Lauderdale 33309
Florida USA

TEL 1 1866 625 5558 – 1 954 874 6325
www.USDENTALDEPOT.com



TWISTER evolution, 120 V

Reliable, optimal, and secure mixing results on demand, highest level of user comfort, time savings, makes work easier.

Item Number
100885

New, user-friendly operating concept with programmable parameters where all functions can be intuitively selected via symbol keys.

Large potential time and work savings with a significant improvement in quality as a result of the pre-mixing function, recallable individual programming of all mixing parameters, even during a mixing process. The one unit for all the most common applications.

Comfortable, intuitive, and therefore very simple operating concept

All mixing parameters can be individually stored in up to 10 programmes

Automatic start-up

Pre-mixing function

Soft-start

Pre- and post-vacuum

Interval (inversion of mixing direction)

Special paddle geometry

The mixing process can be adjusted in 5 sec. intervals

Inlet filter with overflow protection
 Automatic bowl coupling
 Adjustable speed
 7-character LED display
 Practically maintenance-free
 Can be used as a benchtop or wall-mounted unit
 500ml bowl included in the standard delivery
 incl. service package: 3-year warranty, 10-year spare parts availability, loaner unit

Technical details	
Mains voltage:	230 V, 50-60 Hz 120 V, 50-60 Hz
Power consumption:	150 W
Speed:	150 to 450 1/min.
Unit fuse:	T 1.6 A (230 V) T 3.15 A (120 V)
Motor fuse:	T 8 A
Vacuum pump performance:	15 l/min.
Vacuum, max.:	approx. -870 mbar, corresponds to a bowl pressure of approx. 100 mbar
Dimensions (WxHxD):	140 x 325 x 230 mm [5,51 x 12,80 x 9,06 inches]
Weight:	8,1 kg, w./o. bowl

INCLUDED IN BOX

- 1 *Twister evolution* vacuum mixing unit
- 1 500 ml bowl, incl. paddle
- 1 Operating instructions
- 1 Spare parts list
- 1 Quick reference card
- 1 Drilling template
- 1 Fastening set