

*Piezoelectric*  
*Ultrasonic Scaler*  
**ART-P3**

BONART  
BONART  
BONART  
BONART

**Challenge Against Perfection**

**PRODUCT SPECIFICATION**

- Power Supply  
115V $\pm$ 10%-50/60Hz 0.3A 35VA  
230V $\pm$ 10%-50/60Hz 0.15A 35VA
- Output Power  
Perio model : 5W $\sim$ 10W $\pm$ 10%  
Scaling model : 10W $\sim$ 35W $\pm$ 10%  
Turbo model : 15W $\sim$ 40W $\pm$ 10%
- Working Frequency : 26KHz $\sim$ 32KHz
- Water Supply : 25 $\sim$ 60 PSIG  
(172 $\sim$ 414KPa)

