

WE MAKE PORCELAIN BEAUTIFUL

THE NEW VIP UNIVERSAL X-PRESS™ FROM JELRUS TRANSFORMS ALL PORCELAINS INTO EXTRAORDINARY WORKS OF ART.

Introducing the new VIP Universal X-Press multi-function porcelain/pressing furnace. The VIP Universal X-Press features triple-axis precision temperature control, large wide screen displays, an up front pressure gauge for easy viewing, plus an active display of the plunger position during the press cycle. The VIP Universal X-Press will press up to 15 grams per 300 gram ring and has the capability to press aluminum oxide, zirconium and metal. The rapid entry transport system keeps pressing rings hot. Specifically engineered for high temperature ceramics, the VIP X-Press is suitable for use with all low and high fusing porcelains.

Learn more about the new VIP Universal X-Press furnace. Contact your authorized Jelrus dealer or visit our web site at www.jelrus.com.



Jelrus International • Dental Lab Division of Air Techniques, Inc.
800.Jelrus-1 • Fax: 516.433.7684 • worldent1@aol.com • www.jelrus.com