



US DENTAL DEPOT

US DENTAL DEPOT

6555 NW 9th Avenue suite 212
Fort Lauderdale 33309
Florida USA

USA - 1 1866 625 5558 – 1 954 874 6325

Colombia 57 1 489 9051

Peru 51 1 707 3444

Venezuela 58 212 729 1639

Argentina 54 11 5258 7148 Cba 54 351 569 1409

www.USDENTALDEPOT.com

Cavity Preparation Tip



Ball diamond tip
(dia size: \varnothing 1.2m)





US DENTAL DEPOT

US DENTAL DEPOT

6555 NW 9th Avenue suite 212
Fort Lauderdale 33309
Florida USA

USA - 1 1866 625 5558 – 1 954 874 6325













Colombia 57 1 489 9051

Peru 51 1 707 3444

Venezuela 58 212 729 1639

Argentina 54 11 5258 7148 Cba 54 351 569 1409

www.USDENTALDEPOT.com

	 
SB1	Ball diamond tip(dia: \varnothing 1.2m). Preparation of the occlusal surface and the neck of the teeth.
	 
SB3	Half ball diamond tip(dia: \varnothing 1.2m) .Preparation of the interdental surface without damaging the adjacent teeth.
	 
SBL	Half ball diamond tip. Set 45° to the left. Allows get rid of the carious teeth without damaging the adjacent teeth.
	 
SBR	Half ball diamond tip. Set 45° to the right. Allows get rid of the carious teeth without damaging the adjacent teeth.