



**US
DENTAL
DEPOT**

US DENTAL DEPOT

6555 NW 9th Avenue suite 212
Fort Lauderdale 33309
Florida USA

USA - 1 1866 625 5558 – 1 954 874 6325

Colombia 57 1 489 9051

Peru 51 1 707 3444

Venezuela 58 212 729 1639

Argentina 54 11 5258 7148 Cba 54 351 569 1409

www.USDENTALDEPOT.com

Teeth Whitening Unit ECCO TW-A



Teeth Whitening Unit ECCO TW-A

Remark:

1. WHITE PRINCIPLE

Model: ECCO TW Teeth Whitening Unit is a new type whitening technology, based on the condition of strong blue LED light and whitening kit accelerant. By the character of nanometer pellet on the surface and Light decomposition, activate the most H₂O₂ in short order on low or normal temperature. Produce mass oxygen atomy permeate teeth and decompose color molecule, to meet the effect of teeth whitening.



**US
DENTAL
DEPOT**

US DENTAL DEPOT

6555 NW 9th Avenue suite 212
Fort Lauderdale 33309
Florida USA

USA - 1 1866 625 5558 – 1 954 874 6325

Colombia 57 1 489 9051

Peru 51 1 707 3444

Venezuela 58 212 729 1639

Argentina 54 11 5258 7148 Cba 54 351 569 1409

www.USDENTALDEPOT.com

2. CHARACTER

- 1. The specific wavelength and high density light speed catalyze whitening gel.**
- 2. 6 powerful LED tubes emit a high intensive cold blue light with an output at 6000 MW/CM²**
- 3. Special design, innocuous silica gel muzzle, make sure stopping intersectional infection.**
- 4. The new model of power supply make sure steady power output and long-time continuous working even in some cities and places with mutative voltage.**
- 5. Time display could shows treatment time to help dentist master treating time and progress and reduce the waiting time of patient.**

3. TECHNICAL PARAMETERS

- 1. Power supply: AV110-240V,50HZ**
- 2. Power: 60W**
- 3. Light size: 30*80mm²**
- 4. Power density: >6000MW/CM²**
- 5. Light timing: 0-20 mins**
- 6. Power instability <1%**
- 7. Working method: time output & continuous working**
- 8. Coupling: polyester monowave**