



**US
DENTAL
DEPOT**

US DENTAL DEPOT

6555 NW 9th Avenue suite 212
Fort Lauderdale 33309
Florida USA

USA - 1 1866 625 5558 – 1 954 874 6325

Colombia 57 1 489 9051

Peru 51 1 707 3444

Venezuela 58 212 729 1639

Argentina 54 11 5258 7148 Cba 54 351 569 1409

www.USDENTALDEPOT.com

ECCO S-B-I Led Curing light



LED CURING LIGHT ECCO S-B-I

Remark:

1. Product introduction

This product is used in fast solidification by lighting dental resin.

It is a new generation of standard developed LED curing light,

with humanity intelligent chip, the world advanced high intensity



**US
DENTAL
DEPOT**

US DENTAL DEPOT

6555 NW 9th Avenue suite 212
Fort Lauderdale 33309
Florida USA

USA - 1 1866 625 5558 – 1 954 874 6325

Colombia 57 1 489 9051

Peru 51 1 707 3444

Venezuela 58 212 729 1639

Argentina 54 11 5258 7148 Cba 54 351 569 1409

www.USDENTALDEPOT.com

LED light and Li-on rechargeable battery, which is based on YY0055-2002 □□Curing Light□□. It can cure all the dental resin in the market

nowadays, Its performance is further surpass the halogen curing light's,

which can solidify 4-5mm deepness in 20 seconds.

2. Standard accessories

- 1.Curing light + connect line (each one set)
- 2.Special power supply
- 3.Operation manual and maintenance card

3. Technical parameters

Power voltage AC100-240V 50-60Hz

Input power 5W

Light intensity □□600mw/cm²

□□1200mw/cm²

Wavelength 400-500nm