



US DENTAL DEPOT

US DENTAL DEPOT

6555 NW 9th Avenue suite 212
Fort Lauderdale 33309
Florida USA

USA - 1 1866 625 5558 – 1 954 874 6325

Colombia 57 1 489 9051

Peru 51 1 707 3444

Venezuela 58 212 729 1639

Argentina 54 11 5258 7148 Cba 54 351 569 1409

www.USDENTALDEPOT.com



ECCO - LED Curing Light S-IV

Item # 102285

Ecco LED Curing Light is a perfect combination of the highly – efficient LED tube, the micro-control electric circuit and the rechargeable lithium battery. So, it is featured with high-efficiency, outstanding function, convenient operation, and longevity, making it become the latest top product in the global market.

Features:

1. Light intensity: **1400-2200mW/cms**, 5 watt LED
2. BEEP at each 5-second to indicate remaining curing time
3. Various Curing Time: seven preset curing time 5/10/15/20/25/30/40 seconds



US DENTAL DEPOT

US DENTAL DEPOT

6555 NW 9th Avenue suite 212
Fort Lauderdale 33309
Florida USA

USA - 1 1866 625 5558 – 1 954 874 6325

Colombia 57 1 489 9051

Peru 51 1 707 3444

Venezuela 58 212 729 1639

Areentina 54 11 5258 7148 Cba 54 351 569 1409

www.USDENTALDEPOT.com

4. Three preset programs
 - a. Fast mode: Full power all the time
 - b. Soft 1 mode: Generally increasing power in the first 5 seconds, til the end
 - c. Soft 2 mode: Generally increasing power in the first 5 seconds, then half power til the end
5. Equipped with radiating fan to ensure the durable performance
6. Inner Build-in light meter within the base
7. Cord/cordless inchargeable: when running out of battery, the unit can be operated directly with the switching power adapter
8. High output: Full supply of power can satisfy 10-second curing time for 200times

Standard accessories:

- *Handpiece
- *Charging Base
- *Charger
- *Eye-protecting cover x 3pcs
- *Optic light guide
- *User Manual

