



**US  
DENTAL  
DEPOT**

**US DENTAL DEPOT**

6555 NW 9<sup>th</sup> Avenue suite 212  
Fort Lauderdale 33309  
Florida USA

USA - 1 1866 625 5558 – 1 954 874 6325

Colombia 57 1 489 9051

Peru 51 1 707 3444

Venezuela 58 212 729 1639

Areentina 54 11 5258 7148 Cba 54 351 569 1409

[www.USDENTALDEPOT.com](http://www.USDENTALDEPOT.com)

**ECCO IntraOral Camera - 1.3Mp Comps + USB**

**1.3 MEGA PIXELS COMS USB INTRAORAL CAMERA**

Item # 105800



Super-wide automatically  
adjusting-focus lens



Handpiece holder with  
automatic power-off function



Quick connector



Freeze button  
light on/off button



USB2.0 high-speed transmission



# US DENTAL DEPOT

## US DENTAL DEPOT

6555 NW 9<sup>th</sup> Avenue suite 212  
Fort Lauderdale 33309  
Florida USA

USA - 1 1866 625 5558 – 1 954 874 6325

Colombia 57 1 489 9051

Peru 51 1 707 3444

Venezuela 58 212 729 1639

Argentina 54 11 5258 7148 Cba 54 351 569 1409

[www.USDENTALDEPOT.com](http://www.USDENTALDEPOT.com)

### System Configurations



### Specification:

Image sensor	1/4" CMOS
Image resolution	1280X960 pixels, 640X480 pixels, 320X240 pixels
Video output	USB2.0(All digitize)
Lamp	6 pcs white LED(5600K)
Accumulation point	5mm~50mm
photographing visual angle	70°
Handpiece weight	34g
Power input	DC 5V 300mA (Provided by computer directly)
Cable length	2m