



Bleaching Units

Bleaching Units



ECCO-101380 Led Bleach. Unit YS- TW-A 6000

US\$ 497.65

Teeth Whitening Unit ECCO TW-A

1. Bleaching Mode:

ECCO TW Teeth Whitening Unit is a new type whitening technology, based on the condition of strong blue LED light and whitening kit accelerant.

By the character of nanometer pellet on the surface and Light decomposition, activate the most H₂O₂ in short order on low or normal temperature. This produce mass oxygen atomy permeate teeth and decompose color molecule, to meet the effect of teeth whitening.

2. CHARACTER

1. The specifical wave length and high density light speed catalyze whitening gel.
2. 6 powerful LED tubes emit a high intensive cold blue light with an output at 6000 MW/CM²
3. Special design, innocuous silica gel muzzle, make sure stopping intersectional infection.
4. The new model of power supply make sure steady power output and longtimecontinuous working even in some cities and places with mutative voltage.
5. Time display could shows treatment time to help dentist master treating time and progress and reduce the waiting time of patient.

3. TECHNICAL PARAMETERS

1. Power supply: AV110-240volts ,50HZ
2. Power: 60W
3. Light size: 30*80mm²
4. Power density: >6000MW/CM²
5. Light timing: 0-20 minutes non Stop
6. Power instability + - 1%
7. Working method: time output & continuous working
8. Mouth Coupling: polyester monowave

Bleaching products used with this unit :

- LUMAWHITE - 101484 - Lumawhite Patients Kit - bleaching kit 12 patients 35% Perox
- SHOFU- 111727 - Shofu - Niveous Whitening Kit 25 % stabilized hydrogen peroxide
- IBRITE- 101578 - IBRITE- Chairside Tooth Whitening System - 10 patient kit - 35% H₂O₂ - powder/liquid
- IBRITE- 101575 - IBRITE- Chairside Tooth Whitening System - 5 patient kit - 35% H₂O₂ - gel type -





Catalog Index / Índice del Catálogo

Bleaching Units 1