



Camera Intraoral Intraoral Cameras



ECCO - 102234 Digital Dental Camera 3

US\$ 715.00

- Camera 12mpixel
- Imaging Software
- Direct Connect Downloading Cord
- Hard Carrying Case
- Cheek Refractors
- Intraoral Mirror
- Lithium Rechargeable Battery
- Charger
- 2 Gig Memory Card
- Memory Card Reader
- DVD Tutorial and Written Instructions

Simplest intra-oral camera ever! The camera weighs under a pound and fits into the palm of almost any hand. The camera is capable of accomplishing the full range of dental photography with ease from the 1:1 extreme close up of a posterior quadrant to anterior, lateral and all mirror photos including occlusals and linguals. This camera works great for portraits and personal photography as well



ECCO - 101951 Int Oral Cam II Video output

US\$ 175.00

ECCO-CAM-2 Intraoral Camera II

- direct video output
 - Size: 8.5 x 1.2x x0.3in
 - Weight 2.3oz - [batteries included]
 - 0.02in Pinhole Lens
 - Digital Image Processing
 - 2 White LED lamp x 2 -405nm LED x 2
 - Output Video Cable [7ft]- Simply plug the video input jack into any display device and you can view the intraoral condition immediately on the video display screen
 - 2 Blue leds for plaque detection : the remarkable plaque detector mode helps your patients to visualize dental plaque and tartar - which is illuminated in pink
 - 1 Magnifier Lens for enhanced micro view
 - 30 disposable bags - and 4 AAA bateries included
- Made by MIHARU - JAPAN - for Us Dental Depot -



ECCO - 105800 Int.Oral Cam 1.3Mp Comps + USB

US\$ 249.00

ECCO-Cam- 1.3mega pixels Comps + Usb

Intraoral Camera - Includes software
This camera is design with dentist in mind. You can put this mini camera into the camera's focus and form the picture automatically.

Product accessories

- Handset
- Cable
- Handset Braket
- Camera barrier

Features:

- Cd and image freezing feature
- Lamp: 6pcs white LED [5600k]
- Freeze button
- Light on/off button
- Handpiece holder to add directly to the dental unit
- Automatic power-off function
- Superwide automatic adusting lens
- Quickk connect
- Weight: 34g - Power input: DC 5V 300ma [provided by computer directly]





Specifications:

Power: DC 5V+0.5 (1.5A)
Picture: 3800000 Coms
Image Resolution: 1280 x 960 pixels - 640 x 480 pixels - 320 x 240 pixels
Video output: USB 2.0 [All digitize]
Signal System NTSC. PAL
Type of Focus: 5mm ~ 50mm
Length of handset: 200mm
Cable Length: 2m connect port: 5 pinhole

ECCO-106944 Intraoral Camera

US\$ 334.00

- 2 MEGA PIXELS WIRELESS USB
- To use wireless intraoral camera, we can capture teeth images of patient and then transmit into computer on wireless LAN.
- Without any USB cable, it let dentist easily operate it, save images more quickly and conveniently.
- The quality of image and picture will not interrupted and will not lose truth.
- There are 4 channel selections to switch.
- Fast recharge and easily to operate.



SPECIFICATION

Image Sensor - 1/4 sony CCD
Image Resolution - 2.0 mega pixels
Lamp - 6pcs white (5600K)
Photo visual angle - 70
Image Resolution - 640 x 480 pixels
Frequency - 5725~5875 MHz
Transmit Distance - About 10M
Receiving Sensitivity - 95dB



Intraoral Cameras **Catalog Index / Índice del Catálogo** 1