



## Camera Intraoral Intraoral Cameras



### ECCO- 113495 Ecco-Cam 5 Wireless Mem + LCD

US\$ 495.00



Ecco-Cam-5 Wireless Intraoral Camera - 2 megapixel - + TF Memory card + 2.5in LCD Screen --- -- Resolution 2 mega pixels  
Accumulation point 5mm 50mm Angle 70-- Image Resolution 640 x 480 pixels  
**Ecco-Cam-5 Wireless Intraoral Camera - 2 megapixel - + TF Memory card + 2.5in LCD Screen**  
Camera with TF ( + SD ADAPTER ) card, Video output directly ON 2.5in Screen , or connect with computer through USB card reader.

You can also stream video simultaneously with the wireless receiver ( sold separately )

Reliable light source with 6 White LED lights

Easy to operate and convenient to carry, the images display on 2.5 inch LCD and there is no cable attached or maximum distance allowed from the screen.

Adjust the Brightness and contrast in the LCD directly.

5 Buttons to control:

Freez, Save, Previous picture, next picture and delete.

Different kinds of capture angle to display clear pictures, easy operation and automatically focus, camera can be used without docking station, pictures saved permanent.

High power capacity lithium battery, do not need to connect with other adapter.

The box includes :

Main handpiece ( Camera )

2.5 Inches LCD

Wireless transmitter ( and power source for their camera )

Ac Adapter

TF memory card

USB card reader ( your connection to the PC )

Advantage of the unit:

1- The machine has 2.5 Inches LCD already incorporated .

2 -TF memory card for memory 500 PCS picture.

\* Wireless Receiver is available upon request - 90,00 usd- (# 114067 )

### ECCO-105800 Int.Oral Cam 1.3Mp Comps + USB

US\$ 261.42

ECCO-Cam- 1.3mega pixels Comps + Usb

Intraoral Camera - Includes software

This camera is design with dentist in mind. You can put this mini camera into the camera's focus and form the picture automatically.

Product accessories

Handset

Cable

Handset Braket

Camera barrier

Features:

Cd and image freezing feature

Lamp: 6pcs white LED [5600k]

Freeze button

Light on/off button

Handpiece holder to add directly to the dental unit

Automatic power-off function

Superwide automatic adusting lens

Quickk connect

Weight: 34g - Power input: DC 5V 300ma [provided by computer directly]

Specifications:

Power: DC 5V+0.5 (1.5A)





**Picture: 3800000 Coms**

Image Resolution: 1280 x 960 pixels - 640 x 480 pixels - 320 x 240 pixels

Video output: USB 2.0 [ All digitize]

Signal System NTSC. PAL

Type of Focus: 5mm ~ 50mm

Length of handset: 200mm


Cable Length: 2m connect port: 5 pinhole

**ECCO-101951 Intra Oral Camera II Video out**

US\$ 175.00

**ECCO-CAM-2 Intraoral Camera II**

- Direct video output
- Size: Length approx. 200x25x25mm - Weight 35gr - [without battery and cable]
- CMOS/Pinhole Lens
- Digital Image Processing
- Output Video Cable [7ft] - Simply plug the video input jack into any display device and you can view the intraoral condition immediately on the video display screen
- User Manual
- JAPAN - for Us Dental Depot -



**DINE - 113528 Led background Photography**

US\$ 899.00

Dine ILED background Photography Dines new LED Light Panel is designed to completely eliminate shadows in portrait and profile photography. The high intensity light panel provides a constant light making it functional with virtually any digital SLR camera flash -eliminating back flash delay and synchronization issues. **Say Goodbye to Background Shadows!**

Dine's new LED Light Panel is designed to completely eliminate shadows in portrait and profile photography. The high intensity light panel provides a constant light making it functional with virtually any digital SLR camera flash, eliminating back flash delay and synchronization issues.



**LED Light Panel- Shadowless Portraits and Profiles**

Dines new LED light panel is lightweight and extremely thin allowing it to be placed anywhere in the office. The Dine LED Light Panel weighs 8 lbs with lighted dimensions of 33? x 21 and a thickness of 1 . Unlike slave light boxes that need to be built into walls or occupy office space by extending outward, the flat light panel requires no elaborate installation and requires nothing more than wall space. The panel runs on ac power. Also included is a wireless remote control that allows users to control light intensity for brighter or a more subtly lit background. The panel is designed to create esthetically pleasing portraits by eliminating all background shadows. Dines NEW LED panel is available with white and blue backgrounds.

**DINE - 113527 PENTAX Intra-Oral Photography**

US\$ 795.00

Dine Intra-Oral Photography modifier without a ring Flash - Full package Dine Digital Solution Camera Kit

The Dine Digital Solution Camera/Accessory Kit -- Full Package Camera/Accessory Kit including - Full Package

Our simplest intra-oral camera ever!! featuring a Dine modified small format digital camera. The camera weighs under a pound and fits into the palm of almost any hand. The camera is capable of accomplishing the full range of dental photography with ease from the 1:1 extreme close-up of a posterior quadrant to anterior, lateral, and all mirror photos including occlusals and linguals. The camera also works great for portraits and personal photography.

The latest model is waterproof and shock proof making it the most durable dental camera on the market.

Camera features include :





14-megapixels of resolution for large prints  
5X optical zoom (28-140mm equivalent)  
2.7 inch LCD screen with 16:9 wide-screen aspect ratio and anti-reflective coating

Enhanced Digital Microscope mode  
Detachable Macro stand

5 LED macro lights mounted around the lens bezel

Waterproof to depths of up to 33 feet  
Shockproof from drops as high as 5 feet from the ground  
Crushproof ti withstand up to 220 pounds of force

Coldproof/Freezeproof down to 14F/-10C  
Dustproof

720p HD video at 30 frames per second

Pixel Track and Digital Shake Reduction  
Face Detection, Smile Capture, Blink Detection  
Autofocus assist lamp  
Infrared remote control capabilities  
Rechargeable Li-Ion battery  
HDMI port for playing back HD videos and still images on you HDTV  
SD/SDHC/SDXC card slot; Eye-fi compliant

**TE FULL PACKAGE INCLUDEDES**

- ?Dine modified 14mp Camera
- ?Direct Connect downloading cord
- ?Hard Carrying Case
- ?Cheek Retractors
- ?Intra-oral Mirror
- ?Lithium Rechargeable Battery and charger
- ?2 gig Memory Card (SD)
- ?Memory Card Reader
- ?Lifetime Telephone Support and full written instructions
- ?US Warranties

**DINE - 113529 Pixelease Imaging Software**

US\$ 100.00



PixelEase Imaging Software The days of overpriced dental imaging software have ended!!!Dines PixelEase is powerful dental software designed to make using digital images on your practices computer easy The system consists of three components:Photo Base- Photo Studio and Photo Printer. Each component is linked together at the touch of button. This powerful imaging software is designed to rival the highend dental imaging softwares but at a fraction \$The days of overpriced dental imaging software have ended!!!

Dines PixelEase is powerful dental software designed to make using digital images on your practices computer easy. The system consists of three components:Photo Base, Photo Studio and Photo Printer. Each component is linked together at the touch of button. This powerful imaging software is designed to rival the highend dental imaging softwares but at a fraction of the price. Plus lifetime telephone support is included for free.

**ECCO-106944 Wireless Intraoral Cam 2 Mp**

US\$ 338.00



Ecco-Cam-3 Wireless Intraoral Camera #106944

- 2 MEGA PIXELS WIRELESS USB
- Capture teeth images of patient and then transmit into computer on wireless LAN.
- Without any USB cable, lets the dentist easily capture, save images .
- Freeze image and save them for posterior control
- The Image quality and picture will not be interrupted during operation .
- 4 Transmission channels selection .
- Fast Charger
- Transmitter connects to regular USB
- CD with software included - ( save customers name, and re call the pictures later )
- Antifog lenses



**SPECIFICATION**

Image Sensor - 1/4 Sony CCD  
Image Resolution - 2.0 mega pixels  
Led Lamp ( illumination ) - 6pcs white (5600K)  
Angle of Vision - 70 degrees  
Image Resolution - 640 x 480 pixels  
Frequency - 5725-5875 MHz  
Transmit Distance - About 10 meters  
Receiving Sensitivity - 95dB

**1 year warranty**

**ECCO- 113494 Wireless Intraoral Cam 2 Mp+SD**

**US\$ 450.00**



Ecco-Cam-4 Wireless Intraoral Camera - 2 megapixel - + 2gb SD card --- -- Resolution 2 mega pixels Accumulation point 5mm 50mm  
Angle 70-- Image Resolution 640 x 480 pixels **Ecco-Cam-4 Wireless Intraoral Camera - 2 megapixel - + 2gb SD card #113494**

**High Quality images, simple operation and installation. Light and easy mobility, non-spherical lens, superwide automatically adjustin-focus lens.**

**Hand piece holder with automatic power off function. To use wireless intraoral camera, you can caPture teeth images of patient and then transmit into computer through wireless LAN.**

**Without USB cable, it let dentist easily operate it, save images more quickly and conveniently.**

- High Resolution
- High Sensivity
- High Performance
- Compact Size, easy to use and focus.

**US DENTAL DEPOT Inc is an Authorized ECCO Dealer in the US.**



## Catalog Index / Índice del Catálogo

Intraoral Cameras ..... 1