



## Core Build-Up Material



### ALPHA CORE BUILD UPS

The materials exhibits ideal mechanical properties with high compressive and tensile strength yet maintaining a dentin like feel when cutting.



### ALPHA CORE DC DUAL CURE

- This highly filled radiopaque, dual cure composite is designed specifically for core build-ups and cementations of posts and pins.
- The material exhibits ideal mechanical properties with high compressive and tensile strength yet maintaining a dentin like feel when cutting.
- Its ideal handling characteristics allow the materials to be placed with a syringe and depending on the shade of the base paste it can be used to match or contrast the tooth color.
- Available in Natural or Contrast Blue.

Code	Name	Price	+
106664	ALPHA - Alpha-Core Dual Core Build Up Kit - Natural color - 50gr cartridge - automix tip...	US\$ 66.50	+Info
106665	ALPHA - Alpha-Core Dual Core Build Up Kit - Blue color - 50gr cartridge - automix tips 1...	US\$ 66.50	+Info

### ALPHA CORE LC LIGHT CURE

- Alpha-Core LC is a highly radiopaque light cure composite having 5 mm average depth of cure designed specifically for core build-ups and cementations of posts and pins.
- The material exhibits ideal mechanical properties with high compressive and tensile strength yet maintaining a dentin like feel when cutting.
- Available in a unique color changing shade that is blue initially and turns pink when cured then back to blue in 3 to 7 days.



Code	Name	Price	+
106663	ALPHA - Alpha-Core Light Cure Core Build Up - 4.5gr syringe	US\$ 12.25	+Info

### ALPHA CORE SC SELF CURE

- Alpha-Core SC is a self cure composite specifically designed for core build ups and cementations of posts and pins.
- The material exhibits ideal mechanical properties with high compressive and tensile strength yet maintaining a dentin like feel when cutting.
- Its ideal handling characteristics allows the materials to be placed with a syringe and depending on the shade of the base paste it can be used to match or contrast the color.



### Alpha Core SC Self Cure Kit

Code	Name	Price	+
106666	ALPHA - Alpha-Core Self Cure Core Build Up Kit - White color - Catalyst 14gr - Base 14gr...	US\$ 28.00	+Info
106667	ALPHA - Alpha-Core Self Cure Core Build Up Kit - Natural color - Catalyst 14gr - Base 14...	US\$ 28.00	+Info
106668	ALPHA - Alpha-Core Self Cure Core Build Up Kit - Blue color - Catalyst 14gr - Base 14gr ...	US\$ 28.00	+Info

### Alpha Core SC Self Cure Refill

Code	Name	Price	+
106669	ALPHA - Alpha-Core Self Cure Core Build Up - Catalyst 14gr - Refill - Dual cure resin ba...	US\$ 14.00	+Info
106670	ALPHA - Alpha-Core Self Cure Core Build Up - Base White 14gr - Refill - Dual cure resin ...	US\$ 14.00	+Info
106671	ALPHA - Alpha-Core Self Cure Core Build Up - Base Natural 14gr - Refill - Dual cure res...	US\$ 14.00	+Info



106672	ALPHA - Alpha-Core Self Cure Core Build Up - Base Blue 14gr - Refill - Dual cure resin b...	US\$ 14.00	<a href="#">+Info</a>
--------	---	------------	-----------------------

## ECCO-557793 Core Build Up 5gm Syringe 1pk

US\$ 10.87

Ecco-Light Cure Core Build Up Individual 5gm Syringe - [J# 557793 ] - place and shape material to desired form - light cure and then finish with handpiece -



## CORE SHADE CEMENT ( Shofu )

Formulated to provide a dependable, easily detectable, metal-free core build-up.

- Supportive and protective base under stress-bearing fillings and restorations.
- High compressive strength (147 MPa).
- Excellent bonding to tooth structure, metal pins and posts.
- Gray in color for easy identification during placement and finishing
- High fluoride release
- Radiopaque for quick detection on radiographs
- Biocompatible
- Economical
- Superior handling



US DENTAL DEPOT Inc is an Authorized SHOFU Dealer in the US.

### Kit

- Contains:**
- 25g gray powder
  - 10ml liquid
  - Mixing pad
  - Dispenser
  - Instructions

Code	Name	Price	+
111905	SHOFU - Coreshade Base Cement Kit 25gr	US\$ 102.47	<a href="#">+Info</a>

### Refills

- Powder: 25 gr
- Liquid: 10 ml

Code	Name	Price	+
111911	SHOFU - Base Cement Liquid 10ml	US\$ 44.46	<a href="#">+Info</a>
111910	SHOFU - Base Cement Powder 25gr	US\$ 62.31	<a href="#">+Info</a>

## PULPDENT- 039-HC Dual Build Up Kit 20&20gr

US\$ 107.66

- Hardcore - Core Build Up Material Kit - 20gm Base + 20Gr Dual Catalyst - 20 Core Forms
- Mixing Pad - Spatulas - [P# HC ]
- Hard core is a 2 paste dual cure composite restorative and core build up material with fluoride release the off-white color is easily distinguished from tooth structure and does not create shadows incrowns - - Handling is ideal - mixes easily - resists slumping - sets exceptionally hard and is remarkably strong -cures both chemically and by light
- Chemical cure is complete in 8 minutes from beginning of mix
- can be accelerated by using a curing light



**DenTASTIC UNO™**  
5<sup>TH</sup> Generation  
Single Component Dental Adhesive

*DenTASTIC UNO combines the patented PMGDM adhesive primer with hydrophilic resins and acetone solvent into an exceptionally strong single component adhesive.*





## ECCO-557795 LC Core BuildUp 25carp 0.35gr

US\$ 65.21

Ecco-Light Cured Core Build Up 25 carpules 0.35gr each - [ J# - 557795 ] - place and shape material to desired form  
- light cure and then finish with handpiece - - Light Cured Core Build Up 25 carpules 0.35gr each  
- Place and shape material to desired form  
- Light cure and then finish with handpiece



## Core Build-Up Material Core Build-up Forms



## CORE BUILD UP FORMS ( Pulpdent )



### CORE BUILD UP FORMS

Code	Name	Price	+
039-HCF-1_A	PULPDENT - Forms Small 64pk	US\$ 38.06	<a href="#">+Info</a>
039-HCF-2	PULPDENT - Forms Medium 64pk	US\$ 38.06	<a href="#">+Info</a>
039-HCF-3	PULPDENT - Forms Large 64pk	US\$ 38.06	<a href="#">+Info</a>
039-HCF-4	PULPDENT - Forms X-Large 64pk	US\$ 38.06	<a href="#">+Info</a>
039-HCF-AS_A	PULPDENT - Forms Assorted 64pk	US\$ 38.06	<a href="#">+Info</a>



## Catalog Index / Índice del Catálogo

Core Build-Up Material .....	1
Core Build-up Forms .....	3