

Curing Units

Light curing units -



ECCO-102161 Bleaching Tip for Led lamp

Ecco-Bleaching Tip for Led lamp ECCO-Light-S-2 and S-1 - 8mm diameter -

US\$ 47.24



ART- 101253 LC Unit ART-L-3 LED 110v

US\$ 232.96

Bonart Light Curing Unit ART L-3 LED 110volts [B# L3FARTBA1PS0] Plastic Handle - Cordless - Comes with Rechargeable Base - 18 month warranty



ART- 101231 LC Unit ART-L2 110v

US\$ 232.96

Bonart Light Curing Unit ART- L2 [B#] 110vol Handle - Power Cord - Eye Protector - 9.5mm Light Guide Probe - High output light intensity 550W/cm2 - 4 optional digital timing setting- timer memorial function- 10sec Beeper Reminder - Selective Colors of Blue - Purple- Green - Phillip Light Bulb 12v/75w - 1 yr on the PC Board - 1 yr on the accessories -



ECCO- 400106 LIGHT-S-0 Ortho Led 2000nw US\$ 210.00

[ECCO-LIGHT-S-0 LED Curing 2000](#)

Portable Curing Light Wireless 2000
For Orthodontics cementation
Excellent design - Stainless steal body
CE stamp certificate
Lithium Battery

2. Standard accessories:

Curing light
Connect line (each one set)
Special power supply
Operation manual and maintenance card

3. Technical parameters

Power voltage AC100-240V 50-60Hz
Input power 5W
Light intensity 2200mw/cm2
Wavelength 400-500nm
1 Year warranty



ECCO-102285 Light-S-IV Cur/Bleach Lite LED

US\$ 265.00

[ECCO-LIGHT-S-4](#)

Excellent For Bleaching and Curing Orthodontics

- 600-1400-2200mW/cms
- Curing Light
- 5 Watt LED
- Handpiece
- Charging Base
- Lithium Battery Charger
- Eye Protector Screen
- Glass Fiberoptic
- User Manual
- Beep At Each 5 Second to Indicate Curing Time





- Seven Pre-Set **CURING TIME** - 5/10/15/20/25/30/40 Seconds
- Cord/Cordless (Both)
- Three **PRE SET PROGRAMS** (Flat - Intermittent - Incremental)
- Fast-Soft Mode 1- Soft Mode 2
- Light meter for checking wave length anytime



ECCO-101381 LIGHT-S-I Led 1200mw/cm2

US\$ 183.73

LIGHT-S-I Led 1200mw/cm2

1. Product introduction - This product is used in fast solidification by lighting dental resin. It is a new generation of standard developed LED curing light with humanity intelligent chip, the world advanced high intensity LED light and Li-ion rechargeable battery.

It can cure all the dental resin in the market nowadays. Its performance is further surpass the halogen curing light's, which can **solidify 4-5mm deepness in 20 seconds.**

2. Standard accessories:

Curing light

Connect line (each one set)

Special power supply

Operation manual and maintenance card

3. Technical parameters

Power voltage AC100-240V 50-60Hz

Input power 5W

Light intensity 600mw/cm2 and 1200mw/cm2

Wavelength 400-500nm

1 Year warranty



ECCO-101921 LIGHT-S-II Led 1400mw -Display

US\$ 165.00

ECCO LIGHT-S-II 1400 mW Item # 101921

NOW IN SPECIAL OFFER !!!

Carton LED curing light S-II is a perfect combination of the highly efficient LED tube, the micro-control electric circuit and the rechargeable lithium battery. So it is featured with high efficiency, outstanding function, convenient operation, and longevity, making it become the latest top product in the global dentistry market. 1.

It can produce a strong blue light directly.

The depth of the polymerization of resin material can reach 2mm.

2. It is a kind of cableless machine, equipped with rechargeable battery. When charged once, the machine can be used for 40 minutes continuously.

3. The machine can hang directly on the toothchair or the rack for special use for safety and convenience.

4. Convenient to operate:

5. 3 preset curing light time from 10- 20- 30 seconds.

6. The light intensity can reach 500mW/cm2 - 900 mW and 1400 mW

SPECIFICATIONS AND EQUIPMENT CLASS

Power supply: Input voltage: 110 & 220VAC automatic Adaptor

Output voltage: 4.2VDC+/-0.1 Output current: 600+/-50 mA

Internal power supply lithium battery: 3.7 v saturated electricity measures 700mAh charging time: about 2.5 hours

Equipment class: Class A (Medical Appliance Classification Catalog)

Protection from ingress of liquids: None

Lamp: 3.6V AC, 3W

Output wavelength range: 420-480nm (nanometers)

Feature size: 20*198mm





Net weight: about 105g

Warranty: one year Standard accessory include: adapter 1 pc , (European or American type) , sterilize sheet 1pack (50pcs) user's manual



Catalog Index / Índice del Catálogo

Light curing units - 1