



Disposables



SURGICAL ASPIRATOR TIPS (Defend)

Made with smooth edges for maximum patient comfort. The surgical aspirator is molded at a 30 degree angle and fits 11mm suction holders. This item is also autoclavable up to 135 C or may be disposed after single use. Available in 1/8" or 1/4" diameters.



Size: 25/bag
Weight: 1 lbs

Code	Name	Price	+
113653	DEFEND - Surg Aspirator Tips 1/4 25pk	US\$ 4.10	+Info
113652	DEFEND - Surg Aspirator Tips 1/8 25pk	US\$ 4.10	+Info
113651	DEFEND - Surg Aspirator Tips 1/16 25pk	US\$ 4.10	+Info

VENTED ORAL EVACUATION TIPS (Defend)

These evacuators fit all standard evacuator hose ends and contain smooth edges for maximum patient comfort. The straight handle allows for increased control and are convenient single-use disposable tips. One end is vented and the other is non-vented (dual use). Non sterile, 1/8" diameter. Available in standard and short (3") lengths.



Code	Name	Price	+
113649	DEFEND - Vented Oral Evac Tips 6 100pk	US\$ 4.80	+Info
113650	DEFEND - Vented Oral Evac Tips 3 100pk	US\$ 4.80	+Info

X-RAY SENSOR SLEEVES (Defend)

All the popular sensor sleeves are now available. Protect your investment in digital imaging with Defend sleeves which boast leak proof construction with extra soft, tear resistant plastic.



Size: 500/box

Code	Name	Price	+
113634	DEFEND - X-Ray Sensor Sleeve 1 3/8x 8	US\$ 4.41	+Info
111353	DEFEND - X-Ray Sensor Sleeve 1 5/8x 8	US\$ 4.41	+Info
113633	DEFEND - X-Ray Sensor Sleeve 2 1/2x10	US\$ 5.88	+Info
113639	DEFEND - X-Ray Head Sleeves 15x26 250pk	US\$ 21.10	+Info

Disposables

Bibs for Patients



DenTASTIC UNO™ 5TH Generation

Single Component Dental Adhesive

DenTASTIC UNO combines the patented PMGDM adhesive primer with hydrophilic resins and acetone solvent into an exceptionally strong single component adhesive.



EMBRACE™ WetBond™ Class V

Cervical Restorative Resin

Light Cure
Wet-Bondig
Radiopaque

- Bonds to the slightly moist tooth.
- Self-etching to dentin.
- Cures with all lights.
- 4 shades: A2, A3.5, B2, D2





PATIENT TOWELS (Defend)

The unique technique used to bond the layers of tissue and poly securely together eliminates separation of the layers and the reinforced edges provide added strength and durability. 2-ply tissue and 1-ply poly for absorption and patient protection. 13" x 18", available in 5 popular colors

Unit/Case: 500/cs
Weight: 12 lbs



Code	Name	Price	+
113644	DEFEND - Patient Bibs Lavender 500 bx	US\$ 22.95	+Info
113643	DEFEND - Patient Bibs Green 500 bx	US\$ 22.95	+Info
113642	DEFEND - Patient Bibs Pink 500 bx	US\$ 22.95	+Info
113641	DEFEND - Patient Bibs White 500 bx	US\$ 22.95	+Info
111361	DEFEND - Patient Bibs Blue 500 bx	US\$ 22.95	+Info

Disposables Cotton Rolls



DEFEND- 111358 Cotton Rolls Bx 2000

US\$ 13.95

Cotton Rolls # CS-0200

These quality cotton rolls are made from 100% cotton, are non-linting, non-sterile, soft texture, pliable and 20% more absorbent than the leading brand. The virgin cotton holds together better, maintains its shape when wet and will not adhere to mucous membranes.

Shipping:

Size: 1.5" x 3.75" - 2000/bxUnit

Weight: 6 lbs



Disposables Disposable Cups



DEFEND- 111360 Drinking Cups 5oz Blue 1000Bg

US\$ 25.05

Disposable Cups # DC-7001

This drinking cups have a special embossed grip that provides greater strength and makes dispensing easier. These recyclable cups have a lustrous, high-gloss finish and are available in five designer colors.

Size: 5oz

Blue

1000/Bags



Disposables Disposable Mask



→ **FREE SHIPPING IN USA** for orders over 500 us\$ (placed online only)
→ **15 us\$ FLAT RATE** for orders over 100 us\$ → **FREE SHIPPPING in Backorders**

Does not include equipment, heavy powders and hazzardous materials



FACE MASKS (Defend)

All masks meets all criteria for ASTM standards for Moderate Barrier protection. Perfect for dust, airborne particles, moisture, spray and splatter. Great feature, increased protection. One of our most popular infection control products.



MOLDED

The combination of a fiberglass-free with a latex-free elastic band help to ease the mind of people with allergies.

Contains : 50/bx Unit / Weight: 0.5 lbs
(also available by the case : 20 box / 14lb)



Breathe-EZ



Code	Name	Price	+
113605	DEFEND - Defend- Breath-EZ Valueline Molded Mask Blue 50 Box - (No Ear-loop) MK-1106 ...	US\$ 4.70	+Info

Economy



Code	Name	Price	+
113658	DEFEND - Defend- Economy Ear-loop Molded Mask Pink 50 Box - MK-1008 ...	US\$ 8.30	+Info
111370	DEFEND - Defend- Economy Ear-loop Molded Mask Blue 50 Box - MK-1007 ...	US\$ 8.30	+Info

Premium



Code	Name	Price	+
113606	DEFEND - Defend- Premium Ear-Loop Molded Mask Blue 50 Box - MK-1036 ...	US\$ 7.70	+Info
113606-P	DEFEND - Defend- Premium Ear-Loop Molded Mask Pink 50 Box - 10% Smaller - MK-1066 ...	US\$ 7.70	+Info

PLEATED

These fluid-resistant masks contain a soft, 3-ply, pleated material for ultimate "cool-breathing" comfort. Shipping Information -

Contains : 50/bx Unit / Weight: 0.5 lbs
(also available by the case : 20 box / 14lb)





Breathe-EZ



Code	Name	Price	+
111369	DEFEND - Defend- Breath-EZ Ear-loop Pleated Mask Blue 50 Box - MK-1246 ...	US\$ 4.10	+Info

Diffuser (Anti-Fog)

Our most popular mask contains an antifog foam nosepiece that absorbs moisture to prevent eyewear from fogging. Constructed of a soft, 3-ply, pleated material for ultimate comfort. The fiberglass-free material combined with a latex-free elastic ear-loop is ideal for those affected by allergies. The fluid-resistant masks have a BFE and PFE of greater than 99%. Provides primary barrier filtration.



Contains : 50/bx Unit / Weight: 0.5 lbs
(Also available by case: 20 box / 14 lbs)

Code	Name	Price	+
113692	DEFEND - Defend- Diffuser (Anti-Fog) Ear-loop Pleated Mask Pink 50 box - MK-1056 ...	US\$ 7.70	+Info
113691	DEFEND - Defend- Diffuser (Anti-Fog) Ear-loop Pleated Mask Blue 50 box - MK-1046 ...	US\$ 7.70	+Info

ECCO-101734 Molded Face Mask - Blue 50pk

US\$ 175.50



Ecco-Molded Face Mask - Blue - 50 masks per box- sold only by the case [20boxes per case] - [D# 871B]

ECCO-101735 Molded Face Mask - Pink 50pk

US\$ 175.50



Ecco-Molded Face Mask - Pink - 50 masks per box- sold only by the case [20boxes per case] - [D# 871P]

IBRITE- 101579 Non-Woven Face Cover x 100

US\$ 20.62

IBRITE- Non-Woven Face Cover with opening for mouth x 100 - [P# IB100]

ECCO-101736 Pleated Face Mask - Blue 50pk

US\$ 175.50



Ecco-Polymas Pleated Face Mask - Blue - Ear Loop - 50 masks per box- sold by the case [20boxes per case]- [D# 873A]

ECCO-112428 Pleated Face Mask - Blue 50pk

US\$ 4.50



Ecco-Polymas Pleated Face Mask - Blue - Ear Loop - 50 masks per box- - [# 10-18010]

ECCO-101737 Pleated Face Mask - Pink 50pk

US\$ 175.50



Ecco-Polymas Pleated Face Mask - Pink - Ear Loop - 50 masks per box- sold by the case [20boxes per case]- [D# 873P]



**All equipments are also available
in 220 volt. in our website...**

WWW.USDENTALDEPOT.COM





ECCO-112429 Pleated Face Mask - Pink 50pk

US\$ 4.50



Ecco-Polymas Pleated Face Mask - Pink - Ear Loop - 50 masks per box- -(# 10-18110)

Disposables

Disposable needles



DISPOSABLE NEEDLES (C-K Meta)



Code	Name	Price	+
107600	META - META - C-K Ject- Disposable Dental Needles 27g Long 100pk	US\$ 6.21	+Info
107601	META - META - C-K Ject- Disposable Dental Needles 27G Short 100pk	US\$ 6.21	+Info
107602	META - META - C-K Ject- Disposable Dental Needles30G Short 100pk	US\$ 6.21	+Info

WAND- 255-N-W-P-30G Needle Wand Plus (B-D) Box 30G

US\$ 135.00



WAND - STA Handpiece with 30Ga.- 1 Needles- 50 box - Mfg # 5050-301 [[Anestesia Computarizada Aguja]]

WAND - 111599 Wand HP w/o Needle 50bx

US\$ 135.00

Wand HP- / Needle 50 bx (Mfg # WA-2050A)

Disposables

Face Shield



FACE SHIELD (Defend)

This Face Shield is made to fit over glasses and face masks. The adjustable velcro headband allows for a comfortable fit. The shield protects against splashing. It is anti-fog treated and latex and fiberglass free. Available in full and 1/2 shield sizes.



Code	Name	Price	+
113607	DEFEND - Half Face Shield 25pk	US\$ 50.80	+Info
113615	DEFEND - Full Face Shield 25pk	US\$ 53.10	+Info



→ **FREE SHIPPING IN USA** for orders over 500 us\$ (placed online only)
→ **15 us\$ FLAT RATE** for orders over 100 us\$ → **FREE SHIPPPING in Backorders**

Does not include equipment, heavy powders and hazzardous materials



Disposables
Gauze & Sponges



DEFEND- 111357 Cotton Sponges 2x2 S 5000Bg

US\$ 57.00

Cotton Filled Gauze Sponges # CS- 0101

DEFEND has improved the performance of its cotton filled gauze sponges by adding 20% more cotton. Now more absorbent than before, these pure white, 8 ply sponges are shred-resistant and made of 100% cotton filled USP gauze. The 2" x 2" pads are available either sterile or non-sterile.

Shipping:

Size: 5000/csUnit - 2" x 2"

Weight: 6 lbs



NON WOVEN SPONGES (Defend)

A great alternative to cotton filled sponges, these are made from a 4-ply, 32 gram rayon/polyester blend which provides for maximum absorbency and little adhesion to wounds. These are soft, smooth, strong, non-linting and non-sterile. Available in 2" x 2", 3" x 3" and 4" x 4"

Size: 5000/cs / 4000cs / 2000cs

Weight: 2 lbs



Code	Name	Price	+
111359	DEFEND - Non Woven Sponges 2x2 5000Bg	US\$ 27.00	+Info
113646	DEFEND - Non Woven Sponges 4x4 2000Bg	US\$ 39.00	+Info
113645	DEFEND - Non Woven Sponges 3x3 5000Bg	US\$ 38.00	+Info

ETHICON- 113814 Surgifoam Oral Sponge 24pk

US\$ 140.00

Ethicon- Surgifoam Oral Sponge 24pk - Absorbable Gelatin Sponge 1cm x 1cm x 1cm #1969 **Ethicon- Surgifoam Oral Sponge 24pk #1969**

Absorbable Gelatin Sponge 1cm x 1cm x 1cm

A sterile, water-insoluble, malleable, porcine gelatin absorbable sponge intended for hemostatic use by applying to bleeding surface.



Disposables
Gloves Latex



LATEX GLOVE (Defend)

Our Powder-Free Textured Latex Gloves are polymer coated for an easy on / off. The black color and textured finish for enhanced grip make the Blackjack gloves great for the tattoo professional.



Code	Name	Price	+
111363	DEFEND - Latex Glove Pwdr Free S 100Box	US\$ 8.50	+Info
111364	DEFEND - Latex Glove Pwdr Free M 100Box	US\$ 8.50	+Info
111365	DEFEND - Latex Glove Pwdr Free L 100Box	US\$ 8.50	+Info



Black

Box with 100 gloves



Code	Name	Price	+
112893	DEFEND - Latex Black P-Free X-S 100 Box	US\$ 8.50	+Info
112884	DEFEND - Latex Black P-Free M 100 Box	US\$ 8.50	+Info
113807	DEFEND - Latex Black P-Free L 100 Box	US\$ 8.50	+Info

Disposables

Gloves Nitril



NITRILE GLOVES POWDER FREE (Defend)

Nitrile Protection Meets Superior Comfort

Our Powder Free Nitrile Textured Exam Gloves are lightweight, yet highly puncture resistant, have great tactile strength and dexterity and are a great alternative for individuals sensitive to latex materials.

Wearing Defend Nitrile Gloves reduces the risk of developing certain health conditions such as contact dermatitis and latex allergy from latex gloves.

Defend Powder Free Nitrile Textured Exam Gloves Offer:

- Significantly more tactile sensitivity
- Less hand fatigue
- Easier finger and hand movement
- Soft, comfortable feel
- Keeps skin cooler



Code	Name	Price	+
111366	DEFEND - Nitril Glove Pwdr Free S 100bx	US\$ 7.90	+Info
111367	DEFEND - Nitril Glove Pwdr Free M 100bx	US\$ 7.90	+Info
111368	DEFEND - Nitril Glove Pwdr Free L 100bx	US\$ 7.90	+Info

Disposables

Infection Control



DEFEND- 111371 Defend Disinfectant Wipes

US\$ 8.99

DEFEND + PLUS Disinfecting/Deodorizing/Cleaning Wipes # SO-9000

These disinfecting wipes contain a stable, low pH formulated disinfectant and deodorant for use on hard, non-porous surfaces, in hospitals, emergency medical settings, laboratories, medical offices, dental offices where control of cross contamination is required. Containing a dual chain quaternary ammonium/alcohol formula. Nontoxic, Nonallergenic.

Shipping:

Size: 160/canister
Weight: 2 lbs





Disposables
Napkin Holder



ECCO-101738 Bib Chains 5pk

US\$ 11.50



Ecco-Bib chains - 5pk - assorted colors - [Metal Chain - Napkin Holder - Plastic Sleeves Patient Chains clips for Bibs]

Disposables
Saliva Ejectors



ECCO- 113708 Evacuator Tips 100pk

US\$ 7.00



Ecco- High Volume Evacuator Tips 100pk - Product Description: Softer end for patient comfort- opaque plastic Disposable. Vented one side only. Box of 100. **Ecco- High Volume Evacuator Tips 100pk**

Disposable
Color: White
Softer and for patient comfort
Opaque plastic
One side vented

Box of 100.

DEFEND- 111362 Saliva Ejec Clear Reinfor100pk

US\$ 5.00



Defend- Saliva Ejector Clear w/ Blue # SE-7001

Fit all standard saliva ejector hose ends. Pliable, holds their shape and the tip minimizes tissue impingement. Wire reinforced

Shipping:
Size: 100 pk
Weight: 1 lbs

Disposables
Sleeves & Covers



CHAIR SLEEVES

DEFEND Chair Sleeves are made of a clear plastic and are designed to easily fit over most chairs. These disposable chair sleeves are highly durable and easy to use.

#BF-4000 (1/2 Chair) 27.5" X 24 " - 225/box

#BF-4500 (Full Chair) 29" x 79" - 125/box



Code	Name	Price	+
111355	DEFEND - Defend- Half Chair Sleeve 27.5 x 24 Clear 225Box (69.96cm x 61cm) # BF-4000 ...	US\$ 24.45	+Info
113638	DEFEND - Defend- Half Chair Sleeve 27.5 x 24 Clear 125Box (69.96cm x 61cm) - # BF-4500...	US\$ 30.30	+Info



All equipments are also available
in 220 volt. in our website...

WWW.USDENTALDEPOT.COM





DEFEND-113635 Handpiece Sleeve 1 1/2x8 500p

US\$ 8.82



Defend- Handpiece Sleeve #7006

These Hand Piece Barrier Sleeves are made to fit over most handpiece instruments. Made of high quality clear plastic, the Handpiece Barrier Sleeves are simple to use and easy to dispose. 1.5" x 8"

Weight: 0.5 lbs
Contains: 500 Box

HEADREST COVERS (Defend)

These Headrest Covers are designed to fit over the top of headrests. They are easy to apply and dispose of.

Made with soft and comfortable textured tissue outer covering and a poly backing to help prevent moisture soak through. The headrest covers are fluid-resistant and disposable. Available in 10" x 10", 10" x 13", 11" x 9.5" and 14" x 9.5" .



Headrest Cover Clear Plastic



Code	Name	Price	+
113636	DEFEND - Headrest Covers 11x9 1/2 225p	US\$ 6.50	+Info
113637	DEFEND - Headrest Covers 14x9 1/2 225p	US\$ 7.65	+Info

Headrest Cover White



Code	Name	Price	+
113647	DEFEND - Headrest Covers 10x10 500pk	US\$ 22.40	+Info
113648	DEFEND - Headrest Covers 10x13 500pk	US\$ 27.10	+Info

DEFEND- 113640 Light Handle Sleeves 4x5 3/4

US\$ 8.35



Defend- Light Handle Sleeves #BF5000

These Light Handle Sleeves are made of clear plastic and are easy to apply. Made to fit most light handles.

Size: 4" x 5.75"
500/box
Weight: 1 lb

KODAK - 113780 Sensor Sheath Kodak 6100 Size

US\$ 21.00

KODAK- Sensor Sheath Kodak 6100 Size 2 100/Bx (TIDI) HS #8954624 - [Sensor Sleeves]

DEFEND- 113779 Tray Covers White 1000bx

US\$ 28.00

Defend- Tray 8.5 x 12.25in White Bx 1000 # Mfg 1000934

EMBRACE™ WetBond™ Pit & Fissure Sealant

Light Cure • Fluoride
Releasing • Radiopaque

- Contains No Bisphenol A and No Bis-GMA
- Wet-Bonding - Moisture tolerant resin chemistry
- Tooth Integrating - Seals against microleakage
- Margin-free - No chipping, pitting, or staining.





DEFEND- 111354 Syringe Sleeves 2.5x10in 500bx

US\$ 9.00

Syringe Sleeve Covers # BF-3000

Fits most 3-way air/water syringes and have a pre-punched opening making application easy. They can also be used to cover handpieces, saliva ejectors, ultrasonic scalers and more. 2.5" x 10"

Shipping:

Size: 500/boxUnit

Weight: 1 lb



DEFEND- 111356 Tray Sleeves Covers 500bx

US\$ 16.10

Defend- Tray Sleeves # BF-7000

These Tray Sleeves cover the entire instrument tray and are made of clear plastic which is easy to apply and easier to dispose of.

Size:10.5 x 14 (26.67cm x 35.56cm)

Color: Clear

Weight: 1lb -Bx 500



Disposables

Sterilization Pouches



STERILIZATION POUCHES (Defend)

Superior design and function are the trademarks of these self-sealing pouches. Triple sealed seams and strong materials aid in preventing tears.

The blue tinted, transparent film assists in detecting tears should they occur. Indicators for steam, chemical vapor and ETO ensure proper sterilization temperature has been met. Constructed of medical grade paper.



Sterization Self Sealing

Code	Name	Price	+
113616	DEFEND - Self Seal Pouch 2.2x2.7 200pk	US\$ 5.20	+Info
113617	DEFEND - Self Seal Pouch 2.0x2.75 200pk	US\$ 5.35	+Info
113614	DEFEND - Self Seal Pouch 2.25x 4in 200p	US\$ 6.05	+Info
113619	DEFEND - Self Seal Pouch 3.5x5.25 200p	US\$ 6.10	+Info
113618	DEFEND - Self Seal Pouch 2.75x10 200pk	US\$ 6.40	+Info
111375	DEFEND - Self Seal Pouch 3.5x 10in 200p	US\$ 6.55	+Info
113620	DEFEND - Self Seal Pouch 5.25x10 200pk	US\$ 11.50	+Info
113621	DEFEND - Self Seal Pouch 7.5x 13 200pk	US\$ 23.10	+Info
113622	DEFEND - Self Seal Pouch 9 x 15 200pk	US\$ 29.40	+Info

Sterization PLUS Pouches

With all the superior design characteristics of our standard DEFEND pouches, DEFEND+PLUS pouches are **constructed with a built-in, lead-free Internal Indicator**, giving you peace of mind when it comes to the sterilization of your instruments. **The Internal Indicator ensures the correct sterilization temperature is reached where you really need it.**

This eliminates the need for an internal indicator strip inside the pouch with your instruments. In addition, the external indicator will still change color when the optimal temperature is reached.



Code	Name	Price	+
113610	DEFEND - Steril Pouch 2.75x 10in 200p	US\$ 8.20	+Info
113611	DEFEND - Steril Pouch 3.5x 5.25in 200p	US\$ 9.20	+Info
111374	DEFEND - Steril Pouch 3.5x 10in 200p	US\$ 9.20	+Info
113612	DEFEND - Steril Pouch 5.25 x 10in 200p	US\$ 14.50	+Info
113613	DEFEND - Steril Pouch 7.5x 13in 200p	US\$ 28.00	+Info

STERILIZATION POUCHES (Professional Choice)

Each box contains 200 pouches.
used for sterilizing instruments in autoclave.

They are self sealing

See Through

Easy to Use with instructions:

1. Place item in pouch.
2. Peel off strip from adhesive band.
3. Fold over paper flap along dotted line and press down firmly from the center outward.
4. When Vapor or Steam processed, arrow indicator turns brown or when Ethylene Oxide processed, arrow indicator turn yellow.



PROFESSIONAL CHOICE

Code	Name	Price	+
101537	PROF CHOICE - Professional Choice Pouches - Self Sealing - 2-1/4x 4inches 200pkg- 57mm x102mm ...	US\$ 6.30	+Info
101538	PROF CHOICE - Professional Choice Pouches - Self Sealing - 2-3/4x 10inches -200pkg 70mm x 254m...	US\$ 6.30	+Info
101539	PROF CHOICE - Professional Choice Pouches - Self Sealing - 3-1/2x 6-1/2inches -200pkg 89mm x 1...	US\$ 11.20	+Info
101540	PROF CHOICE - Professional Choice Pouches - Self Sealing - 3-1/2x 9inches -200pkg 89mm x 254mm...	US\$ 6.30	+Info



Catalog Index / Índice del Catálogo

Disposables	1
Bibs for Patients	1
Cotton Rolls	2
Disposable Cups	2
Disposable Mask	2
Disposable needles	5
Face Shield	5
Gauze & Sponges	6
Gloves Latex	6
Gloves Nitril	7
Infection Control	7
Napkin Holder	8
Saliva Ejectors	8
Sleeves & Covers	8
Sterilization Pouches	10