

## Etching Materials

### Acid Neutralizer



#### ECCO-567320 Acid Neutralizer 4x1.2cc Syrin

US\$ 20.53



Ecco-Acid Neutralizer Mini Dental kit of 4 x 1.2cc Syringes + 20 tips - [ J# - 567320 ] - In additin to rubber dam- use acid neutralizer to surround and protect tissue areas where hydrofluoric acid is used - use acid neutralizer remaining acid before removal -

### Etching Materials

### Endo-Prep



#### ECCO-567270 Endo Prep 4x1.2cc Syringes

US\$ 27.77



Ecco-Endo Prep 4 x 1.2cc syringes + 20 tips - etch - [ J# - 567270 ] - endo prep is a 20% edta viscous solution used to chelate- lubricate and debrite root canal preparations - apply endo prep directly into the canal orifice or onto endodontic files

### Etching Materials

### Eugenol



#### ECCO-585895 Eugenol 30 cc bottle

US\$ 0.00



Ecco-Eugenol 30cc bottle - etch - [ J# - 585895 ] - oil of eugenol is used as a pulp sedative and as a temporary cement -

### Etching Materials

### Fluorhydric Acid for Porcelain



#### ECCO-567310 Porcelain Etch 4x1.2cc Syringe

US\$ 44.68



Ecco-Porcelain Etching 4 x 1.2cc syringes + 20 tips - etch - [ J# - 567310 ] - 9.6% hydrofluoric acid is used to etch porcelain and composites - neutralize acid- wash- dry- apply silane- and then bond restoration - - caution: both dentist and patient must be properly protected from hydrofluoric acid - neutralize completely with acid neutralizer before removing -

### Etching Materials

### Glycolic Acid



#### DenTASTIC UNO™ 5<sup>TH</sup> Generation

Single Component Dental Adhesive

*DenTASTIC UNO combines the patented PMGDM adhesive primer with hydrophilic resins and acetone solvent into an exceptionally strong single component adhesive.*



→ **FREE SHIPPING IN USA** for orders over 500 us\$ (placed online only)  
→ **15 us\$ FLAT RATE** for orders over 100 us\$ → **FREE SHIPPPING in Backorders**

Does not include equipment, heavy powders and hazzardous materials



### ECCO-567330 Glycolic Acid 2% 4x1.2cc Syrin

US\$ 12.68

Ecco-Glycolic Acid 2% 4 x1.2cc Syringes + 20 tips - etch - [ J# - 567330 ] - 2% glycolic acid is used to remove the outer dentin smear layer from crown preparations before taking an impression - apply to preparation and then rinse with water -



### Etching Materials Phosphoric Acid



### PHOSPHORIC ACID ( Ecco )



Code	Name	Price	+
557618	ECCO - Ecco-Phosphoric Acid 35% 12gr Mint Syringe - etch - ...	US\$ 5.44	+Info
557418D	ECCO - Ecco-Phosphoric Acid 35% 12gr Syringe + 25 tips - etch - ...	US\$ 9.66	+Info
557415	ECCO - Ecco-Phosphoric Acid 35% 2oz Mint Viscous - etch - ...	US\$ 15.70	+Info
557619	ECCO - Ecco-Phosphoric Acid 35% 8oz Mint Viscous - etch - ...	US\$ 42.26	+Info
557419	ECCO - Ecco-Phosphoric Acid 35% 8oz Viscous - etch - ...	US\$ 38.64	+Info
567040	ECCO - Ecco-Phosphoric Acid 35% 4 x 1.2cc syringes + 20 tips - etch - ...	US\$ 12.68	+Info
567055	ECCO - Ecco-Phosphoric Acid 35% 20 x 1.2cc Syringes +40 tips - etch - ...	US\$ 45.89	+Info
557418	ECCO - Ecco-Phosphoric Acid 35% 12gm syringe - etch - ...	US\$ 5.44	+Info
557413	ECCO - Ecco-Red Etching Agent - 38% Phosphoric Acid - 30cc - ...	US\$ 5.44	+Info



## Catalog Index / Índice del Catálogo

Acid Neutralizer .....	1
Endo-Prep .....	1
Eugenol .....	1
Fluorhydric Acid for Porcelain .....	1
Glycolic Acid .....	1
Phosphoric Acid .....	2