



Polishing and Finishing Composite Polishing



GC-FUJI - 108-000252 Refill A X-Coarse (2.6mm) (12)

US\$ 59.95

GC-Fuji-Refill A Coarse (2.6mm) 12 -Ask For Colours- Metal strips electro-plated with aluminium-oxide particles for approximal contouring- finishing and polishing of all kind of restoratives. [p G 05-09-10]



SUPER SNAP (SHOFU)

SuperSnap is designed for easy, fast and safe contouring in addition to finishing and polishing of all microfilled and hybrid composites.

Super-Snap is designed for easy, fast and safe contouring, finishing and polishing of all microfilled and hybrid composites.

- Elastic silicone shank mount for easy placement of the disk onto the mandrel
- Manufactured without a metal center
- Very flexible, ultra-thin disks allow easy access to interproximal areas
- Two sizes ? four grits (silicone carbide and aluminum oxide)

*This photo belongs to item #111398 - Super Snap Rainbow Tech Kit



US DENTAL DEPOT Inc is an Authorized SHOFU Dealer in the US.

Kits

Code	Name	Price	+
111399	SHOFU - Shofu - Super Snap Mini Kit 48 Assorted Disks - Kit Contains: 48 Assorted Disks ...	US\$ 25.45	+Info
111398	SHOFU - Shofu - Super Snap Rainbow Tech Kit - 180pcs- () Contains: 180 Assorted Disk (...)	US\$ 86.39	+Info

48 Single Packs



Code	Name	Price	+
111720	SHOFU - Shofu- Super Snap Singles Mini 48pk - Kit (Mfg # 0503)	US\$ 99.66	+Info
111719	SHOFU - Shofu - Super Snap Singles Standard 48pk Kit (Mfg # 0502)	US\$ 96.30	+Info



Super-Snap Disk Refills

SHAPE CATEGORY	Contouring		Finishing				Polishing		Super-Polishing		
SHAPE											
COLOR/GRIT	Black (coarse)		Violet (medium)				Green (fine)		Red (super fine)		
REFILLS (50 pk)	PN	L506	L507	L508	L528	L509	L522	L501	L519	L502	L521

Code	Name	Price	+
110704	SHOFU - Shofu- Super Snap Dark Violet (Medium) Countouring Disk 50 box - Refill - () ...	US\$ 20.55	+Info
110699	SHOFU - Shofu Super Snap Dark Violet Double Sided (Medium) Disk Bx/50- Refill - (m...	US\$ 22.25	+Info
110687	SHOFU - Shofu Super Snap Black (Coarse) Contour Disk Bx/50 - Refill - (mfg # L506)	US\$ 22.25	+Info
110685	SHOFU - Shofu Super Snap Red (Super Fine) Double Sided Disk Bx/50 - Refill - (mfg ...	US\$ 22.25	+Info
110683	SHOFU - Shofu Super Snap Dark Green (Fine) Double Sided Disk Bx/50 - Refill - (mf...	US\$ 22.25	+Info
107-L508	SHOFU - Shofu- Super Snap Dark Violet (Medium) Disk 50 box - Refill - ...	US\$ 20.86	+Info

Super-Snap Mini Disk Refill

Code	Name	Price	+
110697	SHOFU - Shofu Super Snap Dark Red (Super Fine) Mini Disk Bx/50- Refill - (mfg # L5...	US\$ 22.25	+Info
110689	SHOFU - Shofu Super Snap Black (Course) Mini Disk Bx/50 - Refill - (mfg # L507)	US\$ 22.25	+Info
110693	SHOFU - Shofu Super Snap Light Violet (Medium) Mini Disk Bx/50- Refill - (mfg # L5...	US\$ 22.25	+Info
110695	SHOFU - Shofu Super Snap Dark Green (Fine) Mini Disk Bx/50- Refill - (mfg # L519)	US\$ 22.25	+Info

Super-Snap Buff Disk

Use SuperBuff disks for the final polish on all composites. They are paste-impregnated disks that, when exposed to water, release a high-quality polishing paste.

- Paste-impregnated with aluminum oxide
 - Disposable
 - Eliminates the potential mess associated with polishing paste syringe systems
- 25 Disks



Kits

25 Disks plus 2 Mandrels

Code	Name	Price	+
111400	SHOFU - Shofu- Super Snap Super Buffs Kit 25pk Disks + 2 Ca mandrels (Set) Contains: 2...	US\$ 27.00	+Info

Refills

Code	Name	Price	+
110700	SHOFU - Shofu- Super Snap Buffs Disks 25pk - Refill () (mfg # L523) [discs]	US\$ 22.25	+Info
110701	SHOFU - Shofu- Super Snap Buffs Mini Disks 25pk - Refill () (mfg # L524) [discs]	US\$ 22.25	+Info
111723	SHOFU - Shofu- Super Snap Super Buffs Discks 50pk - Refill - (Mfg # 0536) [discs]	US\$ 30.04	+Info

Super-Snap Mandrels

The improved designed metal mandrel for Super-Snap Disks and Buff Disks creates a secure seating of the disk onto the mandrel.

Available in HP, CA and FG shanks and packaged 6 to a box. Disposable Super-Snap Plastic CA Mandrels are packaged 100 per box.



Code	Name	Price	+
111701	SHOFU - Shofu - Super Snap Hp Mandrels (6 per bx) - Refill (Mfg # 0437)	US\$ 15.00	+Info
111702	SHOFU - Shofu - Super Snap FG Mandrels - Refill (Mfg# 0438)	US\$ 15.52	+Info
111703	SHOFU - Shofu - Super Snap CA Mandrels - Refill (Mfg# 0439)	US\$ 15.52	+Info
111704	SHOFU - Shofu - Super Snap Plastic Mandrels (100 per bx) - Refill (Mfg # 0440)	US\$ 24.02	+Info



Super-Snap Polystrips

Polystrips are color-coded with matching grit size for easy identification and offer the same results as the Super-Snap disks for interproximal areas.

SHAPE CATEGORY	Contouring	Finishing	Polishing	Super-Polishing
SHAPE				
COLOR/GRIT	Black/Violet (coarse/medium)		Green/Red (fine/super fine)	
(100 PK) PN	L525		L526	

Code	Name	Price	+
111438	SHOFU - Shofu - Super Snap Polystrips Green/ Red 100pk (Interproximal Strips) () ...	US\$ 23.02	+Info
113738	SHOFU - Shofu - Super Snap Polystrips Black/ Violet 100pk (Interproximal Strips) () ...	US\$ 23.02	+Info

Polishing and Finishing Diamond Pastes



DIAMOND PASTES (Pulpdent)

For those seeking the absolutely sharpest edge, Diamond Paste is the way to go when "good enough" isn't good enough for you. With grits as fine as 1 micron, it is possible to polish an edge so fine that you'll need a microscope to see anything but perfection.



Code	Name	Price	+
039-SPARK-3	PULPDENT - Spark Diam Pol Paste 3ml syr	US\$ 17.89	+Info
039-SPARK	PULPDENT - Spark Diam Pol Paste 1.2ml 4pk	US\$ 31.98	+Info

ADS- 102403 Diamond Polishing Paste 5gr

US\$ 31.50

Diamond Polishing Paste 5gr Syringe

- For Porcelain and Chrome - Heavy concentration of diamond crystals in a non-greasy.
- Water soluble base.
- Glaze adjusted areas and contacts on porcelain without refiring.
- Polish hard partial allos to a final super luster and rinse clean.



SHOFU- 111725 Directdia Pste-3G Syr 2 tips

US\$ 16.98

Shofu - Direct Dia Paste (3 Grams) # 0558

Recommended use with Super-Snap Buff Disk and Buff Mini Disk.
 Superfine diamond particles in a smooth paste
 Replaces Ultra II porcelain polishing paste
 Use Direct Dia Polishing Paste for: enamel or porcelains

US DENTAL DEPOT Inc is an Authorized SHOFU Dealer in the US.



ECCO-527001 ECCO Diamond Paste Syr 3gr

US\$ 22.94

Ecco-Diamond Polishing Paste Syringe - 3gr -



Polishing and Finishing Ox Alum Polishing Paste



ECCO-527002 Restor. Polishing Paste 4gm

US\$ 9.66

Restorative Polishing Paste

- Individual 4gm syringe Restorative Composites Polishing Paste contains fine aluminum oxide in a water soluble base
- for polishing composite - amalgam - and castings -





SHOFU- 112168 Dura-Polish Ox Alum 20gr

US\$ 91.66

Dura-Polish 20g #0552

Description:

Dura-Polish (20g) is impregnated with up to 73% by weight Al₂O₃ particles, and when pre-polishing dental restoration efficiently achieves an extremely homogenous surface.

For use on Metals, Hard Resin, Hybrid Ceramic Lab Composit and Denture Base Resin.



ECCO-567290 Restorative Polishing 2x1.2cc

US\$ 20.53

Restorative Polishing Paste Mini Kit J# 567290

- Restorative polishing paste contains fine aluminum oxide in a water soluble base - for polishing composite amalgam and castings.

Contains: 2x 1.2cc syringes + 20 tips



Polishing and Finishing Polishing Wheels



RENFERT- 023-90-0000 Polisoft A Wheel 22X3mm-50

US\$ 32.00

Renfert - Polisoft A Wheel 22X3mm

For gentle polishing and powerful smoothing of hard alloys. Extremely effective optimum alternative to rubber and silicone polishers.

Grinding ridges are smoothed away before the final surface treatment. great durability. Results in a matt shine and sealed surface.

Contains: 50 pcs



RENFERT- 023-91-0000 Polisoft- Wheel 22X3Mm-50 Pcs

US\$ 36.00

Renfert - Polisoft Wheel 22X3mm

Effective polishing of precious and non-precious metal alloys.

-Optimum pre-polish alternative to rubber or silicone polishers. Quickly smoothes all metal surfaces.

-Virtually no dust or odour produced. Surface is cooled due to porous nature.

-Achieves a silky matt shine, long life-span.

Contains: 50 pcs



DenTASTIC UNO™ 5TH Generation Single Component Dental Adhesive

DenTASTIC UNO combines the patented PMGDM adhesive primer with hydrophilic resins and acetone solvent into an exceptionally strong single component adhesive.





UNMOUNTED FIN. DISK(Shofu)

Use Dura-Green Wheels (silicon carbide) and Chipless Wheels (silicon carbide) for: porcelains and precious alloys
Use Cut-Off Disks (resin-bonded) and SpruCut Wheels (silicon carbide) for: non-precious alloys and model castings

US DENTAL DEPOT Inc is an Authorized SHOFU Dealer in the US.

SHAPE CATEGORY	No. 7	No. 7a	No. 8	Cut-Off Disk	SpruCut No. 8
SHAPE					
ISO #	190	190	190	250	195
THICKNESS	1.7mm	2.4mm	3mm	-	3mm
DURA-GREEN (12 pk)	PN 0150	0151	0152	-	-
CHIPLESS (12 pk)	PN 0455	0456	0457	-	-
UNM. WHEELS (12 pk)	PN -	-	-	0441 (100 pk)	0160

Code	Name	Price	+
111522	SHOFU - Shofu - Spru Cut Unmounted Wheels 12pk (Mfg # 0160)	US\$ 12.62	+Info
111519	SHOFU - Shofu - Dura-Green Unmounted Wheels #7 12pk (Mfg # 0150)	US\$ 15.54	+Info
111521	SHOFU - Shofu - Dura-Green Unmounted Wheels #8 12pk (Mfg # 0152)	US\$ 15.54	+Info
111710	SHOFU - Shofu - Unmounted Chipless Wheel #7 12pk (Mfg # 0455)	US\$ 16.00	+Info
111711	SHOFU - Shofu - Unmounted Chipless Wheel #8 12pk (Mfg # 0457)	US\$ 16.00	+Info
111705	SHOFU - Shofu - Cutoff Disks 100pk Unmounted Finishing wheels (Mfg # 0441)	US\$ 19.57	+Info



Catalog Index / Índice del Catálogo

Composite Polishing	1
Diamond Pastes	3
Ox Alum Polishing Paste	3
Polishing Wheels	4