



Preventive Products

Caries Detector



ECCO-567080 Caries Detector 4 x 1.2cc syr

US\$ 19.50



Ecco-Caries Detector Mini Dental Kit 4 x 1.2cc syringes + 20 tips - [J# 567080] A red dye which stains soft outer bacteria infected dentine to aid identification and removal - Caries should be verified and removed with hand instruments before using slow speed bur - Ecco-Caries Detector Mini Dental Kit

4 x 1.2cc syringes + 20 tips - [J# 567080]

A red dye which stains soft outer bacteria infected dentine to aid identification and removal of caries (Caries should be verified and removed with hand instruments before using slow speed bur)

ECCO-541081 Caries Detector bottle 15cc

US\$ 15.70



Ecco-Caries Detector Individual 15cc bottle - [J# - 541081] A red dye which stains soft outer bacteria infected dentine to aid identification and removal - Caries should be verified and removed with hand instruments before using slow speed bur -

Preventive Products

Denture Cleaners



REGENT-110545 Kleenite Dental Cleaser

US\$ 7.50

- Fast, Powerfull Cleaning Action
- Brightens and Whitens
- Powerful, long-lasting cleaning action
- 90 Cleanings Net Wt 9oz



Preventive Products

Waterpik Irrigators



EMBRACE™ WetBond™ Seal-n-Shine

*Light Cure • Finishes and
Polishes Restorations • Seals
Composites and Margins*

- *Eliminates final finishing and polishing steps*
- *Does not alter occlusion or anatomy*
- *Bonds to the moist tooth*
- *Cures with all lights*



→ **FREE SHIPPING IN USA** for orders over 500 us\$ (placed online only)
→ **15 us\$ FLAT RATE** for orders over 100 us\$ → **FREE SHIPPPING in Backorders**

Does not include equipment, heavy powders and hazzardous materials



WATERPIK- 102432 WP-360- Waterpik Cordele Sys

US\$ 40.00

Waterpik Cordless Dental System

- 2 jet tips for low pressure
- 2 jet tips for high pressure
- Tips rotates 360 grades for easy access
- Accepts most antimicrobial agents
- 110V
- Pressure Range: 45-75 psi
- Reservoir Capacity: 150 ml
- Pulse Rate (Optimal efficacy range 1200-1600): 1450
- Shipping Weight - 1.4 lbs
- 1 years warranty





Catalog Index / Índice del Catálogo

Caries Detector	1
Denture Cleaners	1
Waterpik Irrigators	1