



Surgical Products

Surgical Dressing



PULPDENT- 101441 Periocare Periodontal Dress 2p

US\$ 41.88

Periocare Periodontal Dressing - 2 Paste System

- [P# PC] Highly elastic periodontal dressing which sets resiliently hard and does not chip or fall apart in the mouth - It assists in tissue placement after periodontal surgery and provides durable protection of tissue.
- Does not support the growth of bacteria - When removing from surgical site the site is clean and the Periocare comes off neatly in large pieces.
- After mixing it is ready to pick up with wet fingers 75-90 seconds - 7 minute working time - sets in 15 minutes.





Catalog Index / Índice del Catálogo

Surgical Dressing 1