



Waxes All types

Light Cure wax



ECCO-557794 Light Cured Wax Kit

US\$ 42.26

Light Cured Wax Kit Contais 2 x 4.5gr Syringes & Separator liquid 5gr - [J# - 557794]

- Light cured wax is a unique resin for fabricating patterns for inlay/onlays and endodontic posts and cores - the material can be handled similarly to inlay wax and then light cured to produce a pattern which can not be deformed
- Lc wax will can be burnt out without any residue -

Waxes All types

Wax Accesories



RENFERT- 1132-0000 Wax Knife Large

US\$ 27.30

Wax Knife Large Pair 132-0000

RENFERT- 1130-0000 Wax Knife Small

US\$ 21.61

Wax Knife Small Pair 130-0000



Catalog Index / Índice del Catálogo

Light Cure wax	1
Wax Accesories	1