



Curing Units



AQUAPRES HYDRAULIC CURING UNIT - Parts (Lang)

Improve the quality of the fabrication process. When combining Jet? Acrylics and Lang Dental's Aquapres?, you can decrease undesirable results, which can be time consuming and costly. The complete fabrication system consistently generates dependable, high-quality results by eliminating porosity in methacrylate acrylics and achieving long-lasting color.

Features:

- Hydraulic, not pneumatic
- 3 Easy-Steps to generate pressure
- Compact design

Benefits:

- Provides a safer work environment over air pressure without the need for an air compressor, heat or electricity
- Creates a dense acrylic for color stable results and increases efficiency
- Requires no extra attachments, space saving

Aquapres Unit Code: #114112



Code	Name	Price	+
114114	LANG - Aquapres 5in O Ring	US\$ 6.90	+Info
114117	LANG - Aquapres Pressure Release Valv	US\$ 28.00	+Info
114115	LANG - Aquapres Bowl	US\$ 32.50	+Info
114119	LANG - Aquapres Handle	US\$ 32.65	+Info
114113	LANG - Aquapres Gauge	US\$ 35.60	+Info
114118	LANG - Aquapres Frame	US\$ 87.10	+Info
114116	LANG - Aquapres Lid	US\$ 98.80	+Info

Curing Units

Curing Unit Parts



KEYSTONE- 114180 Air line Adapter

US\$ 16.40

Keystone- Air line Adapter #9622780

ECCO-103876 Base for ECCO Light

US\$ 50.39

Ecco-Base for ECCO Light - replacement part

ECCO-102238 Light Meter For LED Lamps

US\$ 125.00

Ecco-Light Meter For light curing LED Lamps replacement part



→ **FREE SHIPPING IN USA** for orders over 500 us\$ (placed online only)
 → **15 us\$ FLAT RATE** for orders over 100 us\$ → **FREE SHIPPPING in Backorders**

Does not include equipment, heavy powders and hazzardous materials



LIGHT PARTS (Ecco)

All ECCO branded products, regardless of point of manufacture, are backed by the US DENTAL DEPOT certification of quality control, as well as our direct warranty service and our free technical support in 2 languages, on the phone or by email.



Code	Name	Price	+
101943	ECCO - Charger ECCO Light 110/220V	US\$ 18.90	+Info
101938	ECCO - Fiber Optic For ECCO Lights	US\$ 39.90	+Info
101941	ECCO - Green Button for ECCO Light	US\$ 3.15	+Info
101942	ECCO - Replac Batt for ECCO Light	US\$ 51.44	+Info



Catalog Index / Índice del Catálogo

Curing Units	1
Curing Unit Parts	1