



Model Trimmers Replacement Parts



HANDLER- 103346 #3207GK Gasket Kit

US\$ 29.82

#3207GK Gasket Kit - Includes gasket - cement - use on 10in and 12in models -

RAY FOSTER- 110533 Automatic Solenoid Water Valve

US\$ 111.28

- Automatically turns the water on and off simultaneously with the motor.
A great convenience and time saver.
- Accessorie for Model Trimmer.



RAY FOSTER- 101709 Drain Hose Replacement 2ft

US\$ 8.23

Drain Hose Replacement for model trimmer - 2 ft - [R# M130]

ECCO-103074 ECCO- Front Platform screw

US\$ 5.25

ECCO -Trim plataform screw - [for model trimmers ECCO].
- Tornillo de plataforma delantera.



RAY FOSTER- 110534 Foot Switch - Press On/Off

US\$ 43.87

- Maintained Operation. Press on, Press off.
- Accessorie for Model Trimmer.



RAY FOSTER- 110528 Lock Button w/ Screws & Gasket

US\$ 18.19

- REPLACEMENT PARTS FOR RAY FOSTER MODEL TRIMMER

RAY FOSTER- 110527 Motor Switch

US\$ 11.24

REPLACEMENT PARTS FOR RAY FOSTER MODEL TRIMMER

RAY FOSTER- 110536 Nu-Life Abrasive Cleaning Bar

US\$ 10.70

Accessorie for MODEL TRIMMER

RAY FOSTER- 110535 Orthodontic Work Table - 12in

US\$ 222.56

- Convert your existing 12" Foster Model Trimmer to an Orthodontic Model Trimmer capable of trimming accurate and attractive study models.
- Accessorie for Model Trimmer.





RAY FOSTER- 110529 Rubber Door Gasket w/ Adhesive

US\$ 18.19

- REPLACEMENT PARTS FOR RAY FOSTER MODEL TRIMMER

RAY FOSTER- 110531 Sediment Tray

US\$ 132.68

- Helps Keep the counter free from splatter and drippings.
- Handy tray in front holds models.
- Accessorie for Model Trimmer.



RAY FOSTER- 110530 Splash Shield

US\$ 55.64

- Protect your face and clothing from splatter with an easily installed.
- Fully adjustable splash shield.
- Accessories for Model Trimmer



RAY FOSTER- 110532 Water Spray Attachement

US\$ 39.59

- Great for washing casts or flushing the model trimmer after use.
- Accessorie for Model Trimmer

RAY FOSTER- 110526 Water Spray Tube Assem 10in

US\$ 48.15

- Water Spray Tube Assembly for 10" Model Trimmer
REPLACEMENT PARTS FOR RAY FOSTER MODEL TRIMMER

RAY FOSTER- 110523 Water Spray Tube Assem 12in

US\$ 48.15

Water Spray Tube Assembly for 12" Model Trimmer
REPLACEMENT PARTS FOR RAY FOSTER MODEL TRIMMER



Catalog Index / Índice del Catálogo

Replacement Parts 1