




DENT. UNIT

PORTABLE

 <p>ASEPTICO - AEU-425 Electronic Portable Systems Tr</p> <p><i>Electronic Portable Systems Transport II Portable Electric Dental System - With ASC-10 Ultrasonic piezo Scaler <Price Does Not Include Handpiece >- Features Included - 3Way Air/Water Syringe - Autoclavable 30K Motor that accepts All E-Type Handpie- (High and low speed) - High And Low Volume evacuation - self contained water supply - self containde compressor - on-off foot switch - easy transportable with pull out handle and wheels - 60lbs weight - 110v/220v compatable - - W/O ASC-10 piezo electric Scaler Included [equipo portatil unidad] (0)</i></p> <p>u\$s 3.820,00</p>	 <p>ASEPTICO - 100881 ASEPTICO ADC-01</p> <p><i>ASEPTICO ADC-01 - Portable chair- Comfortable and safe for up to 500 pounds - High strength steel construction - Back adjusts with patient in chair from full upright to "supine position - Adjustable seat height from 8" to 22" - Can be used as standard sit-up chair" "Folded dimensions: 23" x 34" x 8" - Shipping Weight: 70 lb" (0)</i></p> <p>u\$s 1.150,00</p>
 <p>ASEPTICO - 100879 ASEPTICO ADC-08</p> <p><i>ASEPTICO ADC-08 - Portable stool to go with ADC-01 AseptiChair - It folds neatly into a carrying case Shipping Weight: 19 lb (0)</i></p> <p>u\$s 495,00</p>	

Catalog Index / Índice del Catálogo

DENT. UNIT	2
PORTABLE	2