

PINSETTERS



USDD - 100931

USDD Laser Pinsetter Power

"USDD Laser Pinsetter - Laser Light beam indicator for accurate pin location - High Speed motor - Uses the regular standard Pindex drills - Counter sunk "P" type drill - Standard Pindex drill bur - All Metal construction - Sealed Housing - [120Volts Adj. Depth Drilling- 1/8 Id Chuck] - D-pin (Idem Pindex kit) -Pindex Maquina) [P# D24-0605] [R# 1835-4000] Similar RENFERT Top Spin - (0)

u\$s 430,00



BUFFALO - 008-39000

Buffalo Laser Pinsetter Power

"Buffalo Laser Pinsetter Powerrite - Laser Light beam indicator for accurate pin location - Ergonomic design for both standing and sitting operators - High Speed motor - Automatic On/Off Switch - Uses all types of Drills - Counter sunk "P" type drill" includes - All Metal construction - Sealed Housing 120vol (Adj. Depth Drilling- Auto Start- 1/8 Id Chuck) (Idem Pindex kit) -Pindex Maquina)((0)

u\$s 457,00

Catalog Index / Índice del Catálogo

PINSETTERS

2