

STEAMER



WORLIDENT - 100202

Hot Shot Worldent

“Hot Shot Worldent Steamer Vapor Jet jelrus Extra Large Capacity Over one gallon- 4.5 litre “Inox” stainless steel tank.” Cool Dry Steam -Wet Steam - Maximum Steam - See MORE INFO FILE (0)

u\$s 1.028,00

USDD - 100880

USDD STEAMER ST350

USDD Steamer ST350 - No tank - Six Minute Start up - Unlimited liquid capacity - Removes polishing compounds - oil - wax - State of the art electronic controls 1 Year Warranty - 190 Psi with Wand Assembly - 110 and 220 Volts - Made in USA (0)

u\$s 810,00



RELIABLE - 100497

RELIABLE I500b Dental Steam

RELIABLE High preassure steamer cleaner - Stainless steel chamber - Safety double chamber cover - Highly flexible hose - Safety features with redundant systems - 2.5 liters - 50psi pressure - Made in Italy - - US Warantee 1 year - Professional Dental Steam Cleaner - The i500B has been designed from the ground up to deliver superior quality dry steam. Provide the best possible user experience. The i500B offers a quick/efficient way of cleaning away left over wax plaster-waste a FEATURES: i37 steam gun | Automatic shut-off - with low water light indicator | Regular tap water | Pressure gauge manometer | 50 p.s.i. operating pressure = DRY STEAM | 12-gauge heavy-duty wiring | Ceramic terminal block | Stainless steel tank and o (0)

u\$s 595,00



EDG - 100464

VAPORJET Steam cleaner

EDG- VAPORJET High preassure steam cleaner - Aluiminum Steam Chamber - Safety double chamber cover - Highly flexible hose - Safety features with redundant systems - Digital temperature control - Stainless Steel body - US Warantee - MADE IN BRASIL - WWW.EDG.COM.BR - (EDG is the biggest furnace manufacturer in Brasil) Furnaces : also with ONE year warranty on the muffles - SINTERIZATION FURNACES - 2 years warranty on the muffles no questions asked - (0)

u\$s 800,00

Catalog Index / Índice del Catálogo

STEAMER

2
