



**US
DENTAL
DEPOT**


USA

US DENTAL DEPOT
ceo@usdentaldepot.com

Ph (305) 300 7249
Fax (786) 513 2868
Mob (786) 210 3284
Skype: usdentaldepot

MAIL
247 SW 8 STREET #301
Miami - FL 33130

WAREHOUSE
6555 NW 9th Av - Suite 211
Fort Lauderdale, FL 33309



VAC FORMER

 <p>VH - 100883 VH Vacuum Forming (no motor) <i>VH- Vacuum Forming Venturi [neumatic] - Just conect to a compressor - No motor needed - Used for copies of models in plaster in diverse works- Aluminum Structure - Uses round and cquare plaques - 5x 5 - [termoformadora thermoformer vacuumforming] Dimensions 360 x 200x 250 mm - 14 x 8 x 10 in - 110 or 220 volts (0)</i></p> <p>u\$S 228,00</p>	 <p>VH - 100882 VH Vacuum Forming with motor <i>VH- Vacuum Forming with Motor - No need compressor - Uses round and cquare plaques - 5 x 5 inches - Used for copies of models in plaster in diverse works- Aluminum Structure [termoformadora thermoformer vacuumforming] Dimensions 360 x 200x 350 mm - 14 x 8 x 14 in - Motor 900watts - 110 or 220 volts (0)</i></p> <p>u\$S 333,00</p>
 <p>KEYSTONE - 034-7000470 Proline Air Vac 2000- 110V <i>Proline Air Vac 2000 - 110volts - vacuum former - made in usa - 5 year warranty on parts and labor [K# 7000470] Termo Formado Termoformado Maquina --- vacuum forming machine -no motor means unit is very compact and very quiet - works on all types of plastics and appliances- very economical - can be used to layer appliances - unit utilizes the labs internal air compressor to provide maximum amount of air pressure - (0)</i></p> <p>u\$S 360,85</p>	 <p>KEYSTONE - 034-1560010 Vacuum Forming Machine III <i>Vacuum Forming Machine III (Ex - Kevac] 110V - [K# 7000330] 5x5inches Sheets - Original Keystone - Made in USA - 5 Years warranty in parts and labor - all aluminum construction - heavy duty spring and latch mechanism to hold all types of thermoplastics- symmetrically engineered heating element for uniform heating distribution- Vacuum former - [Termo Formado Termoformado Maquina vacuum former] --- [Keystone Abrasives - Keystone Dental Tools] [Us Dental Depot is specialized business in Dental lab equipment & Dental Porcelain Furnaces] (0)</i></p> <p>u\$S 282,60</p>
 <p>KEYSTONE - 101209 Pro-Form Vac. Form Dual 110v <i>Pro-form Dual Chambered Vacuum Former - 110volt - made in USA - 5 year parts and labor [K# 9594130] -vacuum forming dual-chambered deavy duty vacuum especially designed for dental industry all aluminum - two post sliding frame to assure ease in operation- vacuum motor strategically located below a perforated plate to create direct downward suction- (0)</i></p> <p>u\$S 312,25</p>	 <p>KEYSTONE - 101208 Pro-Form Vac. Form Single 110v <i>Pro-form Single Chambered Vacuum Former - 110 voltz - made in USA - 5 year parts and labor - [K# 9594470] - vacuum forming is economical unit that is excellent for fabricating most thermoplastic appliances - hexagonal post that keeps sliding frame from rotating- all aluminum construction- uniform heating distribution- (0)</i></p> <p>u\$S 251,10</p>
 <p>KEYSTONE - 100872 KEYSTONE Pressure Dome <i>KEYSTONE Pressure Dome - Accesory for vaccum forming machine - Adds new extra pressure to the vaccum. Fits in any vacuum unit that makes 5x5 plastics - weight 1lb - sizes package 7 x7 x7 inches - [K#9594050] [Keystone Abrasives - Keystone Dental Tools] [Us Dental Depot is specialized business in Dental lab equipment & Dental Porcelain Furnaces] (0)</i></p> <p>u\$S 113,76</p>	 <p>BUFFALO - 008-80186 Econo-Vac Heavy Duty Vacuum Fo <i>Econo-Vac Heavy Duty Vacuum Forming Machine Econo-Vac- 120volts Vac- High Powered Termo Formado Maquina (0)</i></p> <p>u\$S 384,75</p>

Catalog Index / Índice del Catálogo

VAC FORMER

2