




ZIRCON.ABR

ABRASIVES

 <p>USDD - 100754 USDD Zirconia Grinder - Wheel <i>ECCO Zirc. Grinder - Sintered Alumina Wheel - medium grit - 1pk [J# CX4150] [abrasives- disk-tool - diamond -] [simil dupont - jota- meisinger- shofu] (0)</i></p> <p>u\$ 25,20</p>	 <p>USDD - 100755 USDD Zirconia Grinder Point <i>ECCO Zirc. Grinder - Sintered Alumina Point - medium grit - 1pk ---(J# CX4160) [abrasives- disk-tool - diamond -] [simil dupont - jota- meisinger- shofu] (0)</i></p> <p>u\$ 14,40</p>
 <p>USDD - 100756 USDD Zirconia Grinder Cylinder <i>ECCO Zirc. Grinder - Sintered Alumina Cylinder - medium grit - 1pk--(J# CX4170) [abrasives- disk-tool - diamond -] [simil dupont - jota- meisinger- shofu] [cilinder barrel] (0)</i></p> <p>u\$ 16,20</p>	

POLISHERS

 <p>USDD - 100722 USDD Dura Disc - Coarse - Blue <i>USDD Dura Disc Coarse Disc - BLUE - Silicone Diamond disc - for fast shaping -scratch removal - leaves satin finish on all porcelains and final stage gloss finish on cercon - Alumina and Zirconia copins [Rubber disk] [zirconium zirconia sinterizado diamante] [K#1900522 - J#CX4100] [zirconia polishers] (0)</i></p> <p>u\$ 22,00</p>	 <p>USDD - 100723 USDD Dura Disc - X fine- Gray <i>USDD Dura Disc - X fine- GRAY -RSilicone Diamond disc - for fast shaping -scratch removal - leaves satin finish on all porcelains and final stage gloss finish on cercon - Alumina and Zirconia copins [Rubber disk] [zirconium zirconia sinterizado diamante] [K#19005 - J#CX4120] [zirconia polishers] (0)</i></p> <p>u\$ 22,00</p>
 <p>USDD - 100763 USDD Dura Disc - Fine - Pink <i>USDD Dura Disc Fine Disc - PINK -Silicone Diamond disc - for fast shaping -scratch removal - leaves satin finish on all porcelains and final stage gloss finish on cercon - Alumina and Zirconia copins [Rubber disk] [zirconium zirconia sinterizado diamante] [K#1900522 - J#CX4110] [zirconia polishers] (0)</i></p> <p>u\$ 22,00</p>	 <p>USDD - 100930 USDD Dura disc Ex Coarse- GREE <i>USDD Dura Disc Extra Coarse Disc - GREEN - Silicone Diamond disc - for fast shaping -scratch removal - leaves satin finish on all porcelains and final stage gloss finish on cercon - Alumina and Zirconia copins [Rubber disk] [zirconium zirconia sinterizado diamante] [K#1900522 - J#CX4130] [zirconia polishers] (0)</i></p> <p>u\$ 22,00</p>
 <p>USDD - 100932 USDD Dura Point - Fine - Pink <i>USDD Dura Point Fine Disc - PINK -Silicone Diamond point - for fast shaping -scratch removal - leaves satin finish on all porcelains and final stage gloss finish on cercon - Alumina and Zirconia copins [Rubber disk] [zirconium zirconia sinterizado diamante] J#CX4113] [zirconia polishers] (0)</i></p> <p>u\$ 17,00</p>	 <p>USDD - 100933 USDD Dura Point - X Fine -GRAY <i>USDD Dura Point Extra Fine Disc - GRAY -Silicone Diamond point - for fast shaping -scratch removal - leaves satin finish on all porcelains and final stage gloss finish on cercon - Alumina and Zirconia copins [Rubber disk] [zirconium zirconia sinterizado diamante] J#CX4123] [zirconia polishers] (0)</i></p> <p>u\$ 17,00</p>



USDD - 100934

USDD Dura Point - Coarse-BLUE

USDD Dura Point Coarse - BLUE- Silicone Diamond point - for fast shaping -scratch removal - leaves satin finish on all porcelains and final stage gloss finish on cercon - Alumina and Zirconia copins [Rubber disk] [zirconium zirconia sinterizado diamante] J#CX4103] [zirconia polishers] (0)

u\$s 17,00

Catalog Index / Índice del Catálogo

ZIRCON.ABR	2
ABRASIVES	2
POLISHERS	2