



# US DENTAL DEPOT

## US DENTAL DEPOT

6555 NW 9<sup>th</sup> Avenue suite 212  
Fort Lauderdale 33309  
Florida USA

USA - 1 1866 625 5558 – 1 954 874 6325

Colombia 57 1 489 9051

Peru 51 1 707 3444

Venezuela 58 212 729 1639

Argentina 54 11 5258 7148 Cba 54 351 569 1409

[www.USDENTALDEPOT.com](http://www.USDENTALDEPOT.com)

### Description

Premium Microscope Inspection system with bright 150 watt fiber optic illuminator, long working distance and wide zoom range.

### Zoom Optics

3.25x to 22.5x with included 1/2x reducing lens Long 7" working distance with 1/2x reducing lens 6.5x to 45x with reducing lens removed

### Illumination

The Alpha 1501 provides bright 150 watt light through a fiber optic ring light guide. This projects a 360 degree ring of cool light to virtually eliminate shadows.



## Customize Your Inspection System

### Z2 Zoom Scope



- Magnification from 6.5X to 45X
- Range Allows Full Field of View and Close-Up Inspection

### Vision Dual-Mag Scope




- Choose 10-20x, 10-30x, or 20-40x
- Long Working Distance and Wide Field of View

### Articulating Boom Stand



- 3 1/2 Foot Reach
- Seven Pivot / Rotation Points
- Choose C-clamp or Heavy Base

### Stainless Steel Boom Stand



- Sturdy with Focus Rack
- Stainless Steel Stand
- Choose C-clamp or Heavy Base

### Supplementary Lens (0.5x or 2x)



- 0.5x or 2x lenses available
- Change the Working Distance or Field of View

### Eyepieces



- 10x, 16x, 20x, or 25x
- HighPoint 10x Allows Users to Wear Glasses

### Alpha 1501 Fiber Optic Light



- 150 watt EKE bulb
- Choose Ring Light or Dual Gooseneck

### LED or Fluorescent Ring Light



- Cool White illumination
- Both Have Variable Intensity
- LED is Brighter, Focused Light