

Dappens



SILICONE MIXING CUPS FOR ACRYLICS - KEYSTONE

SILICONE MIXING CUPS FOR ACRYLICS - KEYSTONE

These Silicone Mixing Bowls are made of a special silicone rubber surface which is adaptable for dental acrylics, resins and cements. The special shape with a wide base prevents spillage and aids in dispensing.



Code	Name	Price	+
5920490	KEYSTONE - Silicone Mixing Cups - Small 2 per pack - silicon dappen (K#5920490)...	US\$ 8.25	+Info
5920500	KEYSTONE - Silicone Mixing Cups - Medium -silicon dappen - (K#5920500) tazones d...	US\$ 9.81	+Info
5920510	KEYSTONE - Silicone Mixing Cups - Large - silicon dappen - (K#5920510) tazones d...	US\$ 12.66	+Info

ADS - 118213 Mix Cups-30cc pkg 25

US\$ 6.50

ADS DENTAL - Mixing Cups - 30 Cc pkg 25 - # 8000-140

ADS - Disposable Mixing Cups - 1 OZ

Dappens

#118213 - Mfg # B895-1

Features

- * Economical enough for one time use.
- * Durable enough for continued reuse. Flexible, graduated, autoclavable.
- * Withstands temperatures from 32F to 250F.

Quantity: 25 x 1 oz. or 30cc cups



Dappens

Glass Dappens



ADS - 118214 Mix Jars Acrylic & CAP large

US\$ 12.57

ADS DENTAL - Mixing Jars Acrylic & CAP large - # J959-2ADS - Mixing Jar

CAP - Large 8 Oz

Dappens

118214 - Mfg # J595-2

Features

- * Highly glazed, white porcelain jars used by denture technicians for mixing and holding acrylic monomer dough.



* Large, 8 oz. jar with cap

GLASS DAPPENS

GLASS DAPPENS -ECCO

Disposable dishes are designed for quick and easy mixing of all materials such as: prophylactic paste, fluoride products, composite materials, acrylics, cements and amalgams.

Made of special thermoplastic material which is impervious to all dental materials.

Disposable to prevent cross contamination and reduce clean.



Code	Name	Price	+
0600510	Glass Dappen Dish - Blue - ...	US\$ 2.50	+Info
0600610	Glass Dappen Dish - Clear - ...	US\$ 2.50	+Info
0600710	Glass Dappen Dish - Green ...	US\$ 2.50	+Info

Dappens

Silicone Dappens



PLASDENT - 001-901DD-A Dappen Dishes Assorted Colors 200pk

US\$ 12.82

PLASDENT - Dapdish Dental Dappen Dishes Assorted Colors

Silicon

001-901DD-A - Mfg # 901DD-A



Description:

- * Dapdish.
- * Dental Dappen Dishes Assorted.
- * Manufacturer Code: 901DD-A.
- * Brand: Plasdent.

Color: Colors In White - Blue - Yellow & Red.

Quantity: 200Pcs/Box.

SILICONE ACRYLIC MIXING BOWL - ECCO
SILICONE ACRYLIC MIXING BOWL - NON STICK CUPS - ECCO

Resistant to high temperatures up to 134°C. Flexible silicone rubber material.
Five different sizes. Washable for repeated use. Sucker-cup easily secures bowls
to smooth surfaces.



Code	Name	Price	+
115049	ECCO - Mixing Cup X-Small - 25mm x 25mm Sold by the unit (also called Acrylic N...	US\$ 2.72	+Info
115051	ECCO - Mixing Cup Small - 35mm x 38mm (M# 506md21) - Sold by the unit (also ca...	US\$ 2.72	+Info
101784	ECCO - Mixing Cup Medium - 44mm x 31mm (M# 506md21) - Sold by the unit (also ...	US\$ 3.77	+Info
115050	ECCO - Mixing Cup Large 65mm x 50mm (M# 506md21) -(also called Acrylic Non St...	US\$ 3.77	+Info

Catalog Index / Indice

Dappens	1
Glass Dappens	1
Silicone Dappens	2