

## Die Spacers

### Light Cure die spacers



#### LC DIE SPACERS - YETI

##### YETI - LIGHT CURING - LIGHT SPACER

Advantages of the light curing die spacer LIGHT SPACER:

- \* homogeneous curing
- \* creates a very hard, scratch-resistant and durable surface
- \* extremely smooth surface - no problem when lifting the coping
- \* excellent adhesion on plaster base
- \* unlimited application time cures only when put into the light curing unit



Curing in any available standard unit

Curing time: approx. 2 - 5 min

Delivery: in black UVA protection bottles of 10 ml

***Die LIGHT SPACER can only be used with each other! Not to be mixed with other Spacers!***

Code	Name	Price	+
34-1860213	YETI - Light Curing liquid Die Spacer - 3 micron Clear - 10ml bottle - K# 186021...	US\$ 19.72	<a href="#">+Info</a>
34-1860214	YETI - Light Curing liquid - Red Die Spacer - 10 micron Red - 10ml bottle - # 52...	US\$ 19.72	<a href="#">+Info</a>
34-1860215	YETI - Light Curing liquid - Blue Die Spacer - 10 micron Blue - 10ml bottle - # ...	US\$ 19.72	<a href="#">+Info</a>

#### Die Spacers

### Paint on die spacers



#### DIE SPACER - AMANN GIRRBACH

##### AMANN GIRRBACH - GIROFORM DIE LINK

##### Paint on Die Spacers

##### Description:

Giroform Die Link is an assortment of air-drying die spacers containing innovative chemical primers and nano-materials.

Giroform Die Link adheres well to dental plasters and gives the die a distinctive, perfectly smooth and scratch-resistant surface.

#UsddAmann

**AMANN GIRRBACH is a trademark of Amann Girschbach AG**



Code	Name	Price	+
114581	Giroform Die Link Starter Kit	US\$ 118.46	+Info
114578	Giroform Die Link Silver 15ml	US\$ 19.51	+Info
114577	Giroform Die Link Red 15ml/Bt	US\$ 19.51	+Info
114576	Giroform Die Link Grey 15ml/Bt	US\$ 19.51	+Info
114575	Giroform Die Link Gold15ml/Bt	US\$ 19.51	+Info
114574	Giroform Die Link Blue 15ml/Bt	US\$ 19.51	+Info
114582	Giroform Die Link Thinner 15ml	US\$ 18.12	+Info

## DIE SPACER - KEYSTONE

### DIE SPACER - KEYSTONE

#### Features:

- \* Fast drying
- \* Heat resistant
- \* 1oz bottles
- \* Provides excellent uniform seating of the restoration



Code	Name	Price	+
034-1028525	KEYSTONE - Silver Die Spacer - 1oz bottle - (K# 1028525 ) Fast drying- Heat resi...	US\$ 19.37	+Info
034-1028523	KEYSTONE - Blue Die Spacer - 1oz bottle K# 1028523 [Espaciador Separador Troquel]	US\$ 19.37	+Info
034-1028524	KEYSTONE - Yellow Die Spacer - 1 oz bottles - (K# 1028524 ) Fast drying- Heat re...	US\$ 19.37	+Info
034-1028527	KEYSTONE - Thinner Die Spacer - 1oz bottles - (K# 1028527 ) #NOMBRE? ...	US\$ 19.37	+Info

## DIE SPACER - YETI

### DIE SPACERS - YETI KEYSTONE

#### Description:

Yeti Die Spacers seals plaster and serves to stimulate the space needed for fixing cement between the crown and the die. This item is able to fulfill even the highest demands.

#### Features:

- \* Available in Clear and Silver.
- \* Comes in 20ml bottles.



Code	Name	Price	+
186022	Yeti Silver Die Spacer 20ml	US\$ 19.52	+Info
100718	Clear Transp Spacer 20ml	US\$ 19.52	+Info
34-1860007	Yeti die spacer Thinner 20ml	US\$ 15.38	+Info

## DIE SPACER REFILL - ADS

### DIE SPACER REFILL - ADS

**Features:**

- \* Creating the exact spacing before wax-up is essential.
- \* Produce controlled relief of dies & consistent film thickness.
- \* Die Spacer provides excellent uniform seating of the restorations.
- \* Heat resistant - lies flat.
- \* Easy on your eyes.
- \* One even coat provides approximately 12 microns relief.
- \* Dries in minutes

**Specify color**

- \* Blue, Clear, Cream, Gray, Green, Orange, Red, White, Yellow, Silver and Thinner



## DIE SPACER AND THINNER KIT

### 1 OZ KIT

Code	Name	Price	+
117862	ADS DENTAL - Die Spacer 1 oz. BLUE & 1 oz. thinner kit - # D750-1-BE	US\$ 16.37	+Info
117863	ADS DENTAL - Die Spacer 1 oz. CLEAR & 1 oz. thinner kit - # D750-1-CR	US\$ 16.37	+Info
117864	ADS DENTAL - Die Spacer 1 oz. CREAM & 1 oz. thinner kit - # D750-1-CM	US\$ 16.37	+Info
117865	ADS DENTAL - Die Spacer 1 oz. GREEN & 1 oz. thinner kit - # D750-1-GN	US\$ 16.37	+Info
117866	ADS DENTAL - Die Spacer 1 oz. ORANGE & 1 oz. thinner kit - # 0750-1-OE	US\$ 16.37	+Info
117867	ADS DENTAL - Die Spacer 1 oz. RED & 1 oz. thinner kit - # 0750-1-RD	US\$ 16.37	+Info
117868	ADS DENTAL - Die Spacer 1 oz. WHITE & 1 oz. thinner kit - # D750-1-WE	US\$ 16.37	+Info
117869	ADS DENTAL - Die Spacer 1 oz. YELLOW & 1 oz. thinner kit - # D750-1-YW	US\$ 16.37	+Info
117870	ADS DENTAL - Die Spacer 1 oz. GRAY & 1 oz. thinner kit - # 0750-1-GY	US\$ 16.37	+Info

### ECONOMY 4 BOTTLE KIT

Code	Name	Price	+
117871	ADS DENTAL - Die Spacer economy - BLUE kit ( 4 btls. ) - # D750-6-BE	US\$ 37.85	+Info
117872	ADS DENTAL - Die Spacer economy - CLEAR kit ( 4 btls. ) - # D750-6-CR	US\$ 37.85	+Info
117873	ADS DENTAL - Die Spacer economy - CREAM kit ( 4 btls. ) - # D750-6-CM	US\$ 37.85	+Info
117874	ADS DENTAL - Die Spacer economy - GREEN kit ( 4 btls. ) - # D750-6-GN	US\$ 37.85	+Info
117875	ADS DENTAL - Die Spacer economy - ORANGE kit ( 4 btls. ) - # D750-6-OE	US\$ 37.85	+Info
117876	ADS DENTAL - Die Spacer economy - RED kit ( 4 btls. ) - # D750-6-RD	US\$ 37.85	+Info
117877	ADS DENTAL - Die Spacer economy - WHITE kit ( 4 btls. ) - # D750-6-WE	US\$ 37.85	+Info
117878	ADS DENTAL - Die Spacer economy - YELLOW kit ( 4 btls. ) - # D750-6-YW	US\$ 37.85	+Info

## DIE SPACER KITS STANDARD

Code	Name	Price	+
103077	ADS DENTAL - Die Spacer Kit - Standard - 1oz ADS DENTAL - Die Spacer + 1oz Thinn...	US\$ 13.30	+Info
103078	ADS DENTAL - Die Spacer Kit - Standard - 1oz ADS DENTAL - Die Spacer + 1oz Thinn...	US\$ 14.15	+Info
102404	ADS DENTAL - Die Spacer Kit - Standard - 1oz ADS DENTAL - Die Spacer + 1oz Thinn...	US\$ 14.15	+Info
103079	ADS DENTAL - Die Spacer Kit - Standard - 1oz ADS DENTAL - Die Spacer + 1oz Thinn...	US\$ 14.15	+Info
103076	ADS DENTAL - Die Spacer Kit - Standard - 1oz ADS DENTAL - Die Spacer + 1oz Thinn...	US\$ 13.30	+Info
103081	ADS DENTAL - Die Spacer Kit - Standard - 1oz ADS DENTAL - Die Spacer + 1oz Thinn...	US\$ 13.30	+Info
103080	ADS DENTAL - Die Spacer Kit - Standard - 1oz ADS DENTAL - Die Spacer + 1oz Thinn...	US\$ 13.30	+Info

## DIE SPACER REFILL

### WITHOUT BRUSH CAP

Code	Name	Price	+
102443	ADS DENTAL - Die Spacer Refill Without Brush Cap - Orange 1qt # D750-06 Creating...	US\$ 117.32	+Info
103091	ADS DENTAL - Die Spacer Refill Without Brush Cap - Yellow 1qt Creating the exact...	US\$ 117.32	+Info
103092	ADS DENTAL - Die Spacer Refill Without Brush Cap - Red 1qt Creating the exact sp...	US\$ 117.32	+Info
103093	ADS DENTAL - Die Spacer Refill Without Brush Cap - Clear 1qt Creating the exact ...	US\$ 117.32	+Info
103094	ADS DENTAL - Die Spacer Refill Without Brush Cap - White 1qt Creating the exact ...	US\$ 117.32	+Info
103095	ADS DENTAL - Die Spacer Refill Without Brush Cap - Creamy ...	US\$ 117.32	+Info
103090	ADS DENTAL - Die Spacer Refill Without Brush Cap - Blue 1qt Creating the exact s...	US\$ 117.32	+Info

## DIE SPACER REFILL STANDARD

### 1 OZ

Code	Name	Price	+
117879	ADS DENTAL - Die Spacer refill - BLUE 1 oz. - # D750-2-BE	US\$ 14.54	+Info
117883	ADS DENTAL - Die Spacer refill - CLEAR 1 oz. - # D750-2-CR	US\$ 14.54	+Info
117887	ADS DENTAL - Die Spacer refill - CREAM 1 oz. - # D750-2-CM	US\$ 14.54	+Info
117891	ADS DENTAL - Die Spacer refill - GRAY 1 oz. - # D750-2-GY	US\$ 14.54	+Info
117893	ADS DENTAL - Die Spacer refill - GREEN 1 oz. - # D750-2-GN	US\$ 14.54	+Info
117897	ADS DENTAL - Die Spacer refill - ORANGE 1 oz. - # D750-2-OE	US\$ 14.54	+Info
117901	ADS DENTAL - Die Spacer refill - RED 1 oz. - # D750-2-RD	US\$ 14.54	+Info
117905	ADS DENTAL - Die Spacer refill - WHITE 1 oz. - # D750-2-WE	US\$ 14.54	+Info
117909	ADS DENTAL - Die Spacer refill - YELLOW 1 oz. - # D750-2-YW	US\$ 14.54	+Info

### 4 OZ

Code	Name	Price	+
117880	ADS DENTAL - Die Spacer refill - BLUE 4 oz. - # D750-4-BE	US\$ 32.16	+Info
117884	ADS DENTAL - Die Spacer refill - CLEAR 4 oz. - # D750-4-CR	US\$ 32.16	+Info
117888	ADS DENTAL - Die Spacer refill - CREAM 4 oz. - # D750-4-CM	US\$ 32.16	+Info
117892	ADS DENTAL - Die Spacer refill - GRAY 4 oz. - # D750-4-GY	US\$ 32.16	+Info
117894	ADS DENTAL - Die Spacer refill - GREEN 4 oz. - # D750-4-GN	US\$ 32.16	+Info

117898	ADS DENTAL - Die Spacer refill - ORANGE 4 oz. - # D750-4-OE	US\$ 32.16	<a href="#">+Info</a>
117902	ADS DENTAL - Die Spacer refill - RED 4 oz. - # D750-4-RD	US\$ 32.16	<a href="#">+Info</a>
117906	ADS DENTAL - Die Spacer refill - WHITE 4 oz. - # ID750-4-WE	US\$ 32.16	<a href="#">+Info</a>
117910	ADS DENTAL - Die Spacer refill - YELLOW 4 oz. - # D750-4-YW	US\$ 32.16	<a href="#">+Info</a>
117913	ADS DENTAL - Die Spacer refill - SILVER 4 oz. - # D750-4-SR	US\$ 32.16	<a href="#">+Info</a>

## 8 OZ

Code	Name	Price	+
117881	ADS DENTAL - Die Spacer refill - BLUE 8 oz. - # D750-5-BE	US\$ 50.49	<a href="#">+Info</a>
117885	ADS DENTAL - Die Spacer refill - CLEAR 8 oz. - # D750-6-CR	US\$ 50.49	<a href="#">+Info</a>
117889	ADS DENTAL - Die Spacer refill - CREAM 8 oz. - # D750-5-CM	US\$ 50.49	<a href="#">+Info</a>
117895	ADS DENTAL - Die Spacer refill - GREEN 8 oz. - # D750-5-GN	US\$ 50.49	<a href="#">+Info</a>
117899	ADS DENTAL - Die Spacer refill - ORANGE 8 oz. - # D750-5-OE	US\$ 50.49	<a href="#">+Info</a>
117903	ADS DENTAL - Die Spacer refill - RED 8 oz. - # D750-5-RD	US\$ 50.49	<a href="#">+Info</a>
117907	ADS DENTAL - Die Spacer refill - WHITE 8 oz. - # 0750-5-WE	US\$ 50.49	<a href="#">+Info</a>
117911	ADS DENTAL - Die Spacer refill - YELLOW 8 oz. - # D750-5-YW	US\$ 50.49	<a href="#">+Info</a>

## QUART

Code	Name	Price	+
117882	ADS DENTAL - Die Spacer refill - BLUE qt - # D750-06-BE	US\$ 123.79	<a href="#">+Info</a>
117886	ADS DENTAL - Die Spacer refill - CLEAR qt - # D750-06-CR	US\$ 123.79	<a href="#">+Info</a>
117890	ADS DENTAL - Die Spacer refill - CREAM qt - # D750-06-CM	US\$ 123.79	<a href="#">+Info</a>
117896	ADS DENTAL - Die Spacer refill - GREEN qt - # D750-06-GN	US\$ 123.79	<a href="#">+Info</a>
117900	ADS DENTAL - Die Spacer refill - ORANGE qt - # D750-06-OE	US\$ 123.79	<a href="#">+Info</a>
117904	ADS DENTAL - Die Spacer refill - RED qt - # D750-06-RD	US\$ 123.79	<a href="#">+Info</a>
117908	ADS DENTAL - Die Spacer refill - WHITE qt - # D750-06-WE	US\$ 123.79	<a href="#">+Info</a>
117912	ADS DENTAL - Die Spacer refill - YELLOW qt - # D750-06-YW	US\$ 123.79	<a href="#">+Info</a>

## THINNER

Code	Name	Price	+
117915	ADS DENTAL - Die Spacer standard - THINNER 1 oz. - # D750-3	US\$ 5.06	<a href="#">+Info</a>
117916	ADS DENTAL - Die Spacer standard - THINNER 4 oz. - # D750-7	US\$ 12.57	<a href="#">+Info</a>
117917	ADS DENTAL - Die Spacer standard - THINNER 8 oz. - # D750-10	US\$ 19.08	<a href="#">+Info</a>
117918	ADS DENTAL - Die Spacer standard - THINNER 16 oz. pt - # D750-11	US\$ 29.96	<a href="#">+Info</a>
117919	ADS DENTAL - Die Spacer standard - THINNER qt - # D750-12	US\$ 43.17	<a href="#">+Info</a>

### RENFERT - 023-1954-0500 Pico-Fit Gold 15ml

US\$ 16.66

RENFERT - Pico-Fit Gold 15ml-# 1954-0500 # 19540500 - Espaciador de Modelos-Dowel Separador Die Spacer)RENFERT - Pico Fit Gold 15 ML

Paint on die spacers

# 023-1954-0500 - Mfg # 1950500



**Function & Performance**

- \* High quality die spacer varnish with good contrast on the most popular plasters.

**Advantages**

- \* Short drying period.
- \* Flow property individually controlled using the thinner.
- \* Builds an even layer thickness and a smooth surface.

**Pico-Fit, silver**

- \* Colors should be used in combination with each other.
- \* Purpose: trouble spots will be visible by the wear-off of the gold layer without damaging the model.

**Key functions**

- \* Very resistant to physical impact
- \* Creates an even layer thickness and smooth surfaces
- \* Short drying times

**Technical data**

- \* Description: Pico-Fit gold
- \* Quantity: 15 ml (0.51 fl.oz)
- \* Layer thickness: approx. 14-20 m

US DENTAL DEPOT is an Authorized RENFERT dealer in US

**DIE:MASTER - SYSTEMATIC DIE PREPARATION - RENFERT**

**DIE:MASTER - SYSTEMATIC DIE PREPARATION - RENFERT**

**Function & Performance**

With the die:master system, Renfert offers air-drying die varnish in four layer thicknesses (10, 13, 15 and 20 m) and six opaque colors.

**Die:Master Set**

- \* Faster and clearer access to all components of the system thanks to the stand.
- \* Practical, firm base for drying all types of die.
- \* Quick and clean storage with supplied high quality separating brush with magnetic fixation.

**Die:Master Die Spacers**

- \* Extremely steam-proof and scratch-resistant thanks to high density, glass-like die spacer layer.
- \* Generates maximum bond strength between hardener and die spacer due to primer technology.





\* Die spacer will not puddle on the brush or die, yet it sticks well to line angles due to its extremely thixotropic flowability.

### Die:Master Duo

- \* Cross links the die stone structure for in-depth hardness and high edge-strength.
- \* Maximum bond strength thanks to chemical bond to die spacer layer.
- \* 0 m thickness for a 100% dimensionally-stable die.

### Die:Master Aqua Thinner

- \* Long-term, sustainable quality,- even after frequent use.
- \* Controlled regulation of the viscosity because it can be dispensed accurately drop-by-drop.

### Picosep

Silicone based separator completes the die:master system components in an optimal way.

### #UsddRenfert

## STAND KITS

Code	Name	Price	+
1956-1000	RENFERT - DieMaster Stand KiT - Gold/Silver Color-# 1956-1000 # 19561000 #NAME?	US\$ 86.18	+Info
1956-0900	RENFERT - DieMaster Stand KiT - Blue Color-# 1956-0900 # 19560900 #NAME?	US\$ 86.18	+Info
1956-0800	RENFERT - DieMaster Stand KiT - Red Color-# 1956-0800 # 19560800 #NAME?	US\$ 86.18	+Info
120705	RENFERT - Die:master Set (ivory) (new) - # 19571000	US\$ 86.18	+Info

## EACH

Code	Name	Price	+
120704	RENFERT - Die:master Ivory (new) - # 19570200	US\$ 21.78	+Info
1956-0500	RENFERT - DieMaster Gold Color Die Spacer-# 1956-0500 # 19560500 #NAME?	US\$ 21.77	+Info
1956-0300	RENFERT - DieMaster Blue Color Die Spacer-# 1956-0300 # 19560300 #NAME?	US\$ 21.77	+Info
1956-0200	RENFERT - DieMaster Red Color Die Spacer-# 1956-0200 # 19560200 #NAME?	US\$ 21.77	+Info
1956-0100	RENFERT - DieMaster Gray Color Die Spacer-# 1956-0100 # 19560100 #NAME?	US\$ 21.77	+Info
1956-0400	RENFERT - DieMaster Silver Color Die Spacer-# 1956-0400 # 19560400 #NAME?	US\$ 21.77	+Info
1956-0600	RENFERT - DieMaster Duo Hardener and Primer in One-15ml-# 1956-0600 # 19560600 #...	US\$ 17.07	+Info
1956-0700	RENFERT - DieMaster Thinner 15ml-# 1956-0700 # 19560700	US\$ 17.07	+Info

## FLASH DIE SPACER - ADS



## GOLD

Code	Name	Price	+
117924	ADS DENTAL - Die Spacer Flash Dry refill - gold 8 oz. - # D751-5-GD	US\$ 50.49	+Info
117923	ADS DENTAL - Die Spacer Flash Dry refill - gold 4 oz. - # D751-4-GD	US\$ 32.16	+Info
117922	ADS DENTAL - Die Spacer Flash Dry refill - gold 1 oz. - # D751-2-GD	US\$ 14.54	+Info
117921	ADS DENTAL - Die Spacer Flash Dry economy - GOLD 4 [As. kit - # D751-6-GD	US\$ 37.86	+Info
117920	ADS DENTAL - Die Spacer Flash Dry GOLD ( spacer & thinner ) kit - # D751-1-GD	US\$ 16.37	+Info

## KIT

Code	Name	Price	+
111135	ADS Flash Die Spacer Refill without cap -8oz # D751-5	US\$ 49.05	+Info
102491	ADS DENTAL - Flash Die Spacer Kit - Sliver 1oz + thiner # D751-1 For controlled ...	US\$ 17.30	+Info

## SILVER

Code	Name	Price	+
117929	ADS DENTAL - Die Spacer Flash Dry SILVER 8 oz. - # D751-5-SR	US\$ 50.49	+Info
117928	ADS DENTAL - Die Spacer Flash Dry SILVER 4 oz. - # D751-4-SR	US\$ 32.16	+Info
117927	ADS DENTAL - Die Spacer Flash Dry SILVER 1 oz. - # D751-2-SR	US\$ 14.54	+Info
117926	ADS DENTAL - Die Spacer Flash Dry economy - SILVER 4 btls kit - # D751-6-SR	US\$ 37.86	+Info
117925	ADS DENTAL - Die Spacer Flash Dry SILVER ( spacer & thinner ) kit - # D751-1-SR	US\$ 16.37	+Info

## THINNER

Code	Name	Price	+
117933	ADS DENTAL - Die Spacer Flash Dry THINNER qt - # D751-12	US\$ 41.96	+Info
117932	ADS DENTAL - Die Spacer Flash Dry THINNER 8 oz. - # D751-10	US\$ 17.20	+Info
117931	ADS DENTAL - Die Spacer Flash Dry THINNER 4 oz. - # D751-7	US\$ 12.57	+Info
117930	ADS DENTAL - Die Spacer Flash Dry THINNER 1 oz. - # D751-3	US\$ 5.06	+Info

## PICO-FIT - RENFERT

### PICO-FIT - RENFERT

#### Function & Performance

High quality die spacer varnish with good contrast on the most popular plasters.

#### Advantages

- \* Short drying period.
- \* Flow property individually controlled using the thinner.
- \* Builds an even layer thickness and a smooth surface.

**US DENTAL DEPOT is an Authorized RENFERT dealer in US.**



Code

Name

Price

+



023-1954-0400	RENFERT - Pico-Fit Gold/Silver/Thinner-SeT - #19540400 # 19540400 -2 x 15ml pain...	US\$ 44.34	<a href="#">+Info</a>
023-1954-0600	RENFERT - Pico-Fit Silver 15ml-# 1954-0600 # 19540600 -(Espaciador de Modelos-Do...	US\$ 16.66	<a href="#">+Info</a>

## REFRACTORY DIE SPACER - ADS

### REFRACTORY DIE SPACER - ADS

Compatible with most common best selling porcelains and matched refractory investments. Provides essential die relief for cement during firing. Construct precise fitting porcelain jackets, crowns, and veneers. Use as a block-out material to repair refractory dies and fill in undercuts. One component liquid - no mixing or measuring of powders and liquids. Ready to use. Porcelain restorations release easily from dies. Prevents breakage of veneers, onlays or inlays. Maintains the integrity of the original refractory die. No more grinding of porcelain after firing to adjust space for bonding materials. Easy to see - paints on blue then fades to white when degassed.



### BLUE

Code	Name	Price	+
118321	ADS DENTAL - Refractory Die Spacer Blue 1 oz. - # H750-4	US\$ 22.07	<a href="#">+Info</a>
118322	ADS DENTAL - Refractory Die Spacer Blue 4 oz. - # H750-5	US\$ 65.08	<a href="#">+Info</a>

### WHITE

Code	Name	Price	+
118323	ADS DENTAL - Refractory Die Spacer White 1 oz. - # H750-2	US\$ 32.47	<a href="#">+Info</a>
118324	ADS DENTAL - Refractory Die Spacer White 4 oz. - # H750-3	US\$ 60.83	<a href="#">+Info</a>

## TRU-FIT DIE SPACER - TAUB

### TRU-FIT DIE SPACER - TAUB

A controlled technique for relieving dies before wax-up, by alternating coatings of gold and silver colored liquids over the die for uniform coverage of 25 microns space. Each coat is approximately 6 microns. Liquids are fast setting. Use on stone, amalgam, plated or epoxy dies. It is a high spot indicator. Tru-Fit's gold and silver colors can be used as a visual aid on models for patient education. Gold solution can be used to simulate crowns, inlays and bridgework. Silver solution can be used to simulate amalgam restorations, partial denture framework.

According to studies at a leading university, precision castings will not seat, be as much as 200 microns, unless relieved. To provide adequate space for cement and common internal discrepancies for ALL crowns and inlays, four coats of fast-drying Tru-Fit will provide much better marginal fit by allowing castings to seat more completely. Tru-Fit also greatly reduces the need for occlusal adjusting.



Code	Name	Price	+
138-139-4OZ	Taub-Tru-Fit Kit Cement Spacer 4Oz Gold/Silver & 8Oz Tinner K# 139-4OZ Espacia...	<b>US\$ 95.00</b>	<a href="#">+Info</a>
138-139-2OZ	Taub-Tru-Fit Kit Cement Spacer 2Oz Gold/Silver & 4Oz Tinner K# 139-2OZ Espacia...	<b>US\$ 55.42</b>	<a href="#">+Info</a>
138-139	Taub-Tru-Fit Kit Cement Die Spacer 1oz Gold/Silver & 2oz Thinner k# 03-1001 Esp...	<b>US\$ 33.02</b>	<a href="#">+Info</a>
138-138	Taub-Tru-Fit Kit Cement Spacer 1/2 Oz Gold/Silver & 1Oz Tinner K# 03-1000 Espa...	<b>US\$ 31.95</b>	<a href="#">+Info</a>

## Die Spacers

### Tooth Color Die Spacers



<p><b>YETI - 417079 Color Spacer 7my Red 20ml</b></p> <p>YETI - COLOR SPACER red 7my 1x20ml # 500-0001</p>	<p>US\$ 19.52</p>
<p><b>YETI - 417083 Die Spacer 10my Blue 20ml</b></p> <p>YETI - DIE SPACER blue 10my 1x20ml # 505-00018</p>	<p>US\$ 19.52</p>
<p><b>YETI - 417084 Die Spacer 10my Dentin 20ml</b></p> <p>YETI - DIE SPACER dentin 10my 1x20ml # 506-00018</p>	<p>US\$ 19.52</p>
<p><b>YETI - 417082 Die Spacer 13my Gold 20ml</b></p> <p>YETI - DIE SPACER gold 13my 1x20ml # 504-00018</p>	<p>US\$ 19.52</p>
<p><b>ADS - 118120 Lab Die-Aid 1/2oz</b></p> <p>ADS - Spacer Die-Aid - 0.5 Oz Paint on die Spacers #118120 - Mfg: #PM101-12</p> <p><b>Feature</b></p> <ul style="list-style-type: none"> <li>* Sensation SL Perfect Tooth Color Spacer Die Aid Refill</li> <li>* Specifically designed to be applied over Sensation SL Perfect Tooth Color Spacers.</li> <li>* Non-acetone, non-alcohol hardener seals Tooth Color Spacer and prevents from peeling when waxing the case.</li> <li>* DO NOT USE a normal die hardener or the Sensation SL Perfect Tooth Color Spacer will dissolve.</li> </ul> <p><b>Quantity:</b> 0.5 oz. bottle</p>	<p>US\$ 13.37</p>

### SENSATION SL TOOTH COLOR REFILLS - ADS

#### AMERICAN DENTAL SUPPLY - SENSATION SL TOOTH

- \* Sensation SL Perfect Match Divex Liquid reduces the surface reaction during the pressing procedure on Lithium Disilicate (e.max)
- \* Sensation SL Perfect Match Glaze Paste is available in 2 formulas: "MPC" for porcelain-to-metal and leucite-reinforced glass pressable ceramics and "EZ" for zirconia and lithium disilicate (e.max)
- \* Perfect Match MPC Glaze Paste is a ready-to-use premixed glaze specifically designed for porcelain-to-metal and leucite-reinforced glass pressable ceramics (Cerpress SL)
- \* Is a smooth glaze that allows the restoration to provide a beautiful glaze in one firing
- \* The pre-mix provides the exact consistency for every application
- \* Perfect Match EZ Glaze Paste is a ready-to-use premixed glaze specifically designed for zirconia and lithium disilicate all-ceramic restorations
- \* Is a smooth glaze that allows the restoration to provide a beautiful glaze in one firing
- \* The pre-mix provides the exact consistency for every application, is wear resistant, creating a protective coat and may be applied to polished zirconia
- \* Perfect Match MPC and EZ Glaze Paste will not clump, get bumpy or leave streaks and provides a smooth glaze every time, allowing the ceramist to achieve a perfect glaze without losing the surface anatomy on the restoration



### PERFECT DIE SPACER

Code	Name	Price	+
115025	ADS DENTAL - SENSATION SL Starter Kit - Perfect Match Includes: Tooth COLOR SPAC...	US\$ 165.00	+Info
118244	ADS DENTAL - Perfect Die Spacer kit - # PM101-1	US\$ 165.97	+Info
118245	ADS DENTAL - Perfect Die Spacer refill shade #1 1/2 oz. - # PM101-2	US\$ 17.59	+Info
118246	ADS DENTAL - Perfect Die Spacer #2 0.5oz. - # PM101-3	US\$ 17.59	+Info
118247	ADS DENTAL - Perfect Die Spacer #3 0.5oz. - # PM101-4	US\$ 17.59	+Info
118248	ADS DENTAL - Perfect Die Spacer #4 0.5oz. - # PM101-5	US\$ 17.59	+Info
118249	ADS DENTAL - Perfect Die Spacer #5 0.5oz. - # PM101-6	US\$ 17.59	+Info
118250	ADS DENTAL - Perfect Die Spacer #6 0.5oz. - # PM101-7	US\$ 17.59	+Info
118251	ADS DENTAL - Perfect Die Spacer #7 0.5oz. - # PM101-8	US\$ 17.59	+Info
118252	ADS DENTAL - Perfect Die Spacer #8 0.5oz. - # PM101-9	US\$ 17.59	+Info
118253	ADS DENTAL - Perfect Die Spacer #9 0.5oz. - # PM101-10	US\$ 17.59	+Info
118254	ADS DENTAL - Perfect Die Spacer #10 0.5oz. - # PM101-11	US\$ 17.59	+Info
118255	ADS DENTAL - Perfect Die Spacer #11 0.5oz. - # PM101-13	US\$ 17.59	+Info
118256	ADS DENTAL - Perfect Die Aid refill 0.5oz. - # PM101-12	US\$ 13.37	+Info

### REFILL

Code	Name	Price	+
115020	ADS DENTAL - Tooth Color Spacer Kit 11 natural shades- Sensation SL Perfect Matc...	US\$ 166.59	+Info
115021	ADS DENTAL - Tooth Color Spacer Refill natural shades- Sensation SL Perfect Matc...	US\$ 17.75	+Info
115027	ADS DENTAL - Tooth Color Spacer Refill natural shades- Sensation SL Perfect Matc...	US\$ 17.75	+Info

115028	ADS DENTAL - Tooth Color Spacer Refill natural shades- Sensation SL Perfect Matc...	<b>US\$ 17.75</b>	<a href="#">+Info</a>
115029	ADS DENTAL - Tooth Color Spacer Refill natural shades- Sensation SL Perfect Matc...	<b>US\$ 17.75</b>	<a href="#">+Info</a>
115030	ADS DENTAL - Tooth Color Spacer Refill natural shades- Sensation SL Perfect Matc...	<b>US\$ 17.75</b>	<a href="#">+Info</a>
115031	ADS DENTAL - Tooth Color Spacer Refill natural shades- Sensation SL Perfect Matc...	<b>US\$ 17.75</b>	<a href="#">+Info</a>
115032	ADS DENTAL - Tooth Color Spacer Refill natural shades- Sensation SL Perfect Matc...	<b>US\$ 17.75</b>	<a href="#">+Info</a>
115022	ADS DENTAL - Die Aid Refill 0.5oz Natural shades- Sensation SL Perfect Match # P...	<b>US\$ 13.12</b>	<a href="#">+Info</a>

## Catalog Index / Indice

Light Cure die spacers .....	1
Paint on die spacers .....	1
Tooth Color Die Spacers .....	10