

Wires Clasps & Arch wires

Ball Clasps



CLASPS - KEYSTONE

CLASPS - KEYSTONE

* Anchorage of removable appliances, frequently applied on incisors and small posterior teeth.

* The ball clasp can be bent in different manners under the condition that the ball is always inserted gingivally in the interdental area.

* Its application is normally with short clinical crowns or with teeth that don't have retentive undercuts (for example lower incisors).

* It should lie in the interdental area, gingival of the contact point. For this, it is necessary to slightly carve the papila.



Our Stainless Steel Ball Clasps come in two different sizes:

* Small - .028" dia. (.7mm)

* Large - .032" dia. (.8mm)

Comes in either: 10/Pkg or 100/Pkg.

10PK

Code	Name	Price	+
034-1270004-10	KEYSTONE - 0.024in / 0.7 mm Ball Clasps Pk/10 - (MFG # 1270004-10)	US\$ 8.51	+Info
034-1270005	KEYSTONE - 0.028in / 0.7mm Ball Clasps Pk/10 - (MFG # 1270005)	US\$ 8.51	+Info
034-1270010	KEYSTONE - 0.032in / 0.8mm Ball Clasps Pk/10 - (MFG #1270010)	US\$ 8.51	+Info
034-1270021-10	KEYSTONE - 0.036in / 0.9mm Ball Clasps Pk/10 - (MFG # 1270021-10)	US\$ 8.51	+Info
034-1270022-10	KEYSTONE - 0.040in / 1.0 mm Ball Clasps Pk/10 - [MFG #1270022-10]	US\$ 9.10	+Info

100PK

Code	Name	Price	+
034-1270004	KEYSTONE - 0.024in / 0.6 mm Ball Clasps Pk/100 - (Mfg # OD1620)	US\$ 32.39	+Info
034-1270006	KEYSTONE - 0.028in / 0.7mm Ball Clasps Pk/100- (MFG # OD1600) OB28/100	US\$ 32.39	+Info
034-1270020	KEYSTONE - 0.032in (0.8mm) Ball Clasps Pk/100 - (MFG # OD1621)	US\$ 32.39	+Info
034-1270021	KEYSTONE - 0.036in / 0.9mm Ball Clasps Pk/100 - (MFG # OD1622) OB36/100	US\$ 32.39	+Info
034-1270022	KEYSTONE - 0.040in / 1.0 mm Ball Clasps Pk/100 - (MFG #OD1623)	US\$ 32.39	+Info

CLASP WIRE - KEYSTONE

STAINLESS STEEL CLASP WIRE - KEYSTONE

These Stainless Steel Clasp Wires are available in half-round and round configurations.



Features:

- * Stainless Steel Clasp Wire - Gauge Half Round
- * Stainless Steel Clasp Wire - Gauge Round

HALF-ROUND CLASP WIRE - KEYSTONE

Code	Name	Price	+
034-1900720	KEYSTONE - 12 gauge Clasp Wire - Stainless HALF-ROUND 5pk x (30cm - 12in) (...	US\$ 7.58	+Info
034-1900730	KEYSTONE - 14gauge Clasp Wire - Stainless HALF-ROUND 5pk x (30cm - 12in) # O...	US\$ 7.58	+Info
034-1900750	KEYSTONE - 16gauge Clasp Wire- Stainless HALF-ROUND 5pk x (1.2mm (120) 30c...	US\$ 7.58	+Info

ROUND CLASP WIRE - KEYSTONE

Code	Name	Price	+
034-1900760	Keystone - 17gauge Clasp Wire - Stainless Round 10pk - (1.10mm (110) 300mm -...	US\$ 7.58	+Info
034-1900770	KEYSTONE - 18gauge Clasp Wire - Stainless Round 10pk - (1.00mm (100) 300mm ...	US\$ 7.58	+Info
034-1900780	Keystone - 19gauge Clasp Wire - Stainless Round 10pk - (0.90mm (090) 300mm -...	US\$ 11.50	+Info
034-1900790	KEYSTONE - 20gauge Clasp Wire - Stainless Round 10pk x (0.80mm (080) 30cm...	US\$ 11.50	+Info

Wires Clasps & Arch wires

Retainer Wire



RETAINER WIRE - ECCO

RETAINER WIRE - ECCO

- * Diamond drawn medium-hard temper for retainers and functional appliances
- * Easy to bend all loops and will not fracture
- * Fifty or ten 14 lights per tube
- * Also 50gm and 500gm spools



10pk

Code	Name	Price	+
106977	ECCO - 30gauge Ligature Retainer Wire .25mm x14 inch log - 10pk (Alambre de or...	US\$ 5.90	+Info
106979	ECCO - 28gauge Ligature Retainer Wire .30mm x14 inch log - 10pk (Alambre de or...	US\$ 6.13	+Info
106981	ECCO - 27gauge Ligature Retainer Wire .36mm x14 inch log - 10pk (Alambre de or...	US\$ 6.13	+Info

106983	ECCO - 25gauge Ligature Retainer Wire .45mm x14 inch log - 10pk (Alambre de or...	US\$ 6.13	+Info
--------	--	-----------	-----------------------

50pk

Code	Name	Price	+
106970	ECCO - 30gauge Ligature Retainer Wire .25mm size x14 inch log - 50pk (Alambre ...	US\$ 17.98	+Info
106972	ECCO - 28gauge Ligature Retainer Wire .30mm size x14 inch log - 50pk (Alambre ...	US\$ 17.98	+Info
106974	ECCO - 27gauge Ligature Retainer Wire .36mm size x14 inch log - 50pk (Alambre ...	US\$ 17.98	+Info
106976	ECCO - 25gauge Ligature Retainer Wire .45mm size x14 inch log - 50pk (Alambre...	US\$ 22.47	+Info


Wires Clasps & Arch wires

ROACH CLASPS



KEYSTONE - 034-1271160 Roach Clasp Stra Bicuspid sMALL

US\$ 18.09



KEYSTONE - Roach Clasp Straight Angle Bicuspid SMALL- Stainless - 10pk # 1271160
 #1271160KEYSTONE - Roach Clasp Med Bic Small Med Straight
 Roach Clasps
 # 034-1271160 - Mfg # 1271160

Features:

- * This Stainless Steel Wire Clasp is bicuspid medium. Items come in straight, left, and right.
- * Package 10/Pk


Us Dental Depot Is An authorized Keystone Leader In Us

ROACH CLASPS - KEYSTONE

STAINLESS STEEL WIRE CLASPS

This Stainless Steel Wire Clasp is bicuspid medium. Items come in straight, left, and right.

Available in package of 10
 Size: Medium



Code	Name	Price	+
034-1271210	KEYSTONE - Roach Clasp Right Angle Bicuspid Medium - Stainless - 10pk #1271210 #...	US\$ 18.09	+Info
034-1271190	KEYSTONE - Roach Clasp Straight Angle Bicuspid Medium - Stainless - 10pk # 12711...	US\$ 18.09	+Info
034-1271200	KEYSTONE - Roach Clasp Left Angle Bicuspid Medium - Stainless - 10pk #1271200 #1...	US\$ 18.09	+Info

Wires Clasps & Arch wires

Spooled Ligat & Retainer Wire



SPOOLED LIGATURE WIRE - BESQUAL

SPOOLED LIGATURE WIRE - BESQUAL

High quality stainless steel ligature wire is super smooth for easy tying. Soft stainless steel with a bright, diamond die high quality finish. Spools for the office that prefers to form their own ligatures.

Available in 0.8in; 0.9in and 1.0in

Weight: 50gr



Code	Name	Price	+
123011	BESQUAL - Spooled Ligature Wire .008in Diameter - 20GA 0.8mm (0.032in) x 11m - (...	US\$ 9.98	+Info
123012	BESQUAL - Spooled Ligature Wire .008in Diameter - 19GA 0.9mm (0.036in) x 9m - (...	US\$ 9.98	+Info
123013	BESQUAL - Spooled Ligature Wire .008in Diameter - 18GA 1.0mm (0.040in) x 7m - (...	US\$ 9.98	+Info

SPOOLED LIGATURE WIRE - ECCO

SPOOLED LIGATURE WIRE - ECCO

Bright Finish with Soft Temper for Easy Use

Weight: 1pound/455grams



Code	Name	Price	+
113556	ECCO - 32gauge Spooled Ligature Wire .008in Diameter T304/304L #6 - (0.20mm Dia...	US\$ 37.08	+Info
113557	ECCO - 31gauge Spooled Ligature Wire .009in Diameter T304/304L #7 - (0.23mm Dia...	US\$ 37.08	+Info
113558	ECCO - 30gauge Spooled Ligature Wire .010in Diameter T304/304L #8- (0.25mm Diam...	US\$ 37.08	+Info
113559	ECCO - 29gauge Spooled Ligature Wire .011in Diameter - T304/304L #9 (0.28 mm Di...	US\$ 37.08	+Info
113560	ECCO - 28gauge Spooled Ligature Wire .012in Diameter T304/304L #10 - (0.31mm Di...	US\$ 37.08	+Info
113561	ECCO - 26gauge Spooled Ligature Wire .014in Diameter T304/304L #11 - (0.36mm Di...	US\$ 37.08	+Info

Wires Clasps & Arch wires

Spring Tensile Wire



SPRING TENSILE ROUND WIRE - ECCO

SPRING TENSILE ROUND WIRE - ECCO



Code	Name	Price	+
113562	ECCO - Spring Tensile Round Wire .012Diam T304 SPR 14 - Round Orthodontic Wire ...	US\$ 8.99	+Info
113563	ECCO - Spring Tensile Round Wire .014 Diam - Round Orthodontic Wire [# OD114]	US\$ 8.99	+Info
113564	ECCO - Spring Tensile Round Wire .016Diam T304 SPR 14 - Round Orthodontic Wire ...	US\$ 8.99	+Info
113565	ECCO - Spring Tensile Round Wire .018Diam T304 SPR 14 - Round Orthodontic Wire ...	US\$ 8.99	+Info
113566	ECCO - Spring Tensile Round Wire .020Diam T304 SPR 14 - Round Orthodontic Wire ...	US\$ 8.99	+Info
113567	ECCO - Spring Tensile Round Wire .022 DiamT304 SPR 14 - Round Orthodontic Wire ...	US\$ 8.99	+Info
113568	ECCO - Spring Tensile Round Wire .024Diam T304 SPR 14 - Round Orthodontic Wire ...	US\$ 8.99	+Info
113569	ECCO - Spring Tensile Round Wire .025Diam T304 SPR 14 - Round Orthodontic Wire ...	US\$ 8.99	+Info

Catalog Index / Indice

Ball Clasps	1
Clasp Wires	1
Retainer Wire	2
ROACH CLASPS	3
Spoiled Ligat & Retainer Wire	4
Spring Tensile Wire	4