

## Zirconia abrasives

### CAD/CAM Burs



### AMANN GIRRBACH DIAMOND CAD/CA

#### M MILLING BURS - ATTRITOR

AMANN GIRRBACH DIAMOND CAD/CAM MILLING BURS - ATTRITOR



Available in 0.4mm - 1mm - 1.4mm - 1.8mm

Code	Name	Price	+
124124	ATTRITOR - Amann Girschbach Diamond Cad Cam Bur 0.6mm for disilicate and feldspar	US\$ 55.00	+Info
124125	ATTRITOR - Amann Girschbach Diamond Cad Cam Bur 1mm for disilicate and feldspar	US\$ 55.00	+Info
124126	ATTRITOR - Amann Girschbach Diamond Cad Cam Bur 1.4mm for disilicate and feldspar	US\$ 55.00	+Info
124127	ATTRITOR - Amann Girschbach Diamond Cad Cam Bur 1.8mm for disilicate and feldspar	US\$ 55.00	+Info

### CERAMILL MOTION ORIGINAL - AMANN GIRRBACH

High Precision, long life sensing and milling tools for cost effective machining of ceramill blanks.



#### Available sizes:

- \* 1.2mm
- \* 0.2mm
- \* 0.3mm and additional special patterns.

Material removal and dust removal are optimally controlled due to the tool geometry tailored to the material.

Extensive assortment of milling cutters for the management of a wide range of indications. Custom selection for Zirconia Abutments selected from Amann Girschbach Burs and Readers

Sold by unit or as a set

#UsddAmann

AMANN GIRRBACH is a trademark of Amann Girschbach AG

Code	Name	Price	+
760751	AMANN GIRRBACH- Zirconia Bur - F0.6 - #760751 ( for item T0.6 ) Used for Implant...	US\$ 73.43	+Info

760709	AMANN GIRRBACH- Zirconia Bur - F0.9mm #10 - #760709 ( for item T0.9mm ) Used for...	US\$ 73.43	+Info
760711	AMANN GIRRBACH- Zirconia Bur - F1.2 - #760711 ( for item T1.2 ) Used for Implant...	US\$ 63.68	+Info
760721	AMANN GIRRBACH- Zirconia Bur - F2 - #760721 ( for item T2 ) Used for Implants & ...	US\$ 67.68	+Info
760720	AMANN GIRRBACH- Zirconia Bur - UF2mm #28 - #760720 ( for item UT2mm ) Used for I...	US\$ 73.43	+Info
760730	AMANN GIRRBACH- Zirconia Bur UF3mm - #760730( for item UT3mm ) Conical tools wi...	US\$ 73.43	+Info
760701	AMANN GIRRBACH- Zirconia Bur - F4 - #760701 ( for item T4 ) Used for Implants & ...	US\$ 73.43	+Info
760740	AMANN GIRRBACH- Zirconia Bur - CF2mm #14 - #760740 ( for item CT2mm ) Used for I...	US\$ 73.43	+Info
760761	AMANN GIRRBACH- Zirconia Bur - SF2mm - #760761( for item ST2mm ) Whool without ...	US\$ 73.43	+Info
760771	AMANN GIRRBACH- Zirconia Reader SF1.2mm - #760771( for item ST1.2mm ) Shank-ty...	US\$ 73.43	+Info
760781	AMANN GIRRBACH- Zirconia Bur SF0.6mm - #760781 ( for item ST0.6mm ) Thin shank-...	US\$ 73.43	+Info
760791	AMANN GIRRBACH- Zirconia Bur KF2.5mm #760791- ( for item KT2.5mm ) Ball tool f...	US\$ 73.43	+Info
760892	AMANN GIRRBACH- Zirconia Reader - KT2mm #17 ( KT1.8mm ) - #760892 ( for item KF2...	US\$ 63.68	+Info
760881	AMANN GIRRBACH- Zirconia Reader ST0.6mm - #760881 ( for item SF0.6mm ) Thin sha...	US\$ 63.68	+Info
760891	AMANN GIRRBACH- Zirconia Reader KT2.5mm - #760891( for item KF2.5mm ) Ball tool...	US\$ 63.68	+Info
760871	AMANN GIRRBACH- Zirconia Reader ST1.2mm - #760871 ( for item SF1.2mm ) Shank-ty...	US\$ 63.68	+Info
760861	AMANN GIRRBACH- Zirconia Reader ST2mm - #760861 ( for item SF2mm ) Whool withou...	US\$ 63.68	+Info
760851	AMANN GIRRBACH- Zirconia Reader - T0.6 - #760851 ( for item F0.6 ) Used for Impl...	US\$ 63.68	+Info
760840	AMANN GIRRBACH- Zirconia Reader - CT2mm #13 - #760840 ( for item CF2mm ) Used fo...	US\$ 63.68	+Info
760831	AMANN GIRRBACH- Zirconia Reader - T3 - #760831 ( for item F3 ) Used for Implants...	US\$ 63.68	+Info
760830	AMANN GIRRBACH - Zirconia Reader UT3mm - #760830 ( for item UF3mm ) Conical too...	US\$ 63.68	+Info
760821	AMANN GIRRBACH- Zirconia Reader - T2 - #760821 ( for item F2 ) Used for Implants...	US\$ 62.23	+Info
760820	AMANN GIRRBACH- Zirconia Reader - UT2mm #27 - # 760820 ( for item UF2mm ) Used f...	US\$ 63.68	+Info
760811	AMANN GIRRBACH- Zirconia Reader - T1.2 - #760811 ( for item F1.2 ) Used for Impl...	US\$ 62.23	+Info
760809	AMANN GIRRBACH- Zirconia Reader - T0.9mm #9 - #760809 ( for item F0.9mm ) Used f...	US\$ 63.68	+Info
760801	AMANN GIRRBACH- Zirconia Reader - T4 - #760801 ( for item F4 ) Used for Implants...	US\$ 63.68	+Info
760792	AMANN GIRRBACH- Zirconia Bur - KF1.8mm ( KF2mm ) #18 - #760792 ( for item KT1.8...	US\$ 73.43	+Info

## SHARK - 416754 AMANN GIRRBACH Cad Bur 0.3mm

US\$ 21.00



AMANN GIRRBACH Cad/Cam Coated Milling Bur 0.3mm - each - for Ceramill RotoMotion -  
Made in Germany **SHARK - Amann Girrbach CAD/CAM - Ceramill Motion Milling Bur -  
0.3mm**

**Burs**

**# 416754**

### Description:

- \* Ultra-precision.
- \* Hard, wear-resistant material.
- \* Last longer.
- \* Can be applied to milling the Zirconia disc, PMMA disc, Acrylic disc for CAD/CAM

**Quantity: 1 pc**

## AMANN GIRRBACH DLC CAD/CAM MILLING BUR S - AMANN GIRRBACH

### DLC CAD/CAM MILLING BURS - AMANN GIRRBACH

Ceramill Roto Motion Burs with Magnet Holder

**Features:**

- \* Ultra-precision.
- \* Hard, wear-resistant material.
- \* Last longer.
- \* Can be applied to milling the Zirconia disc, PMMA disc, Acrylic disc for CAD/CAM

**Available in Size**

- \* 0.9mm
- \* 1.0mm
- \* 1.2mm
- \* 2.0mm
- \* 2.5mm

#UsddAmann

AMANN GIRRBACH is a trademark of Amann Girrbach AG



Code	Name	Price	+
114391	AMANN GIRRBACH- Ceramill Roto Motion 1.0mm bur with magnetic holder 1pk #76060...	US\$ 110.52	+Info
114394	AMANN GIRRBACH- Ceramill Roto Motion 1.2mm bur with magnetic holder 1pk #7606...	US\$ 110.52	+Info
114395	AMANN GIRRBACH- Ceramill Roto Motion 2.0mm bur with magnetic holder 1pk #7606...	US\$ 110.52	+Info
114392	AMANN GIRRBACH- Ceramill Roto Motion 2.5mm bur with magnetic holder 1pk #760605...	US\$ 110.52	+Info

## AMANN GIRRBACH DLC CAD/CAM MILLING BURS - S HARK

### AMANN GIRRBACH DLC CAD/CAM MILLING BURS - SHARK

**Features**

- \* Ultra-precisi
  - \* Hard, wear-resistant material
  - \* Last longer
  - \* Milling Bur 2.5mm - 1.0mm - 0.6mm
- Can be applied to milling the Zirconia disc, PMMA disc, Acrylic disc for CAD/CAM

#UsddAmann

AMANN GIRRBACH is a trademark of Amann Girrbach AG



Code	Name	Price	+
------	------	-------	---

115577	AMANN GIRRBACH Cad/Cam Coated Milling Bur 0.6mm - each - for Ceramill RotoMotio...	US\$ 21.00	<a href="#">+Info</a>
115571	AMANN GIRRBACH Cad/Cam Coated Milling Bur 1mm - each - for Ceramill RotoMotion-...	US\$ 21.00	<a href="#">+Info</a>
115565	AMANN GIRRBACH Cad/Cam Coated Milling Bur 2.5mm - each - for Ceramill RotoMotio...	US\$ 21.00	<a href="#">+Info</a>

## CEREC COMPATIBLE COMPACT MILLING BURS - SHARK

### CEREC COMPACT COMPATIBLE MILLING BURS- SHARK

The Cerec Inlab Sirona Milling Instruments Are The Machine Parts, Known As Burs, That Do The Actual Milling Of The Composite Blocks In Their Respective Unit.

For Cerec Inlab Compact Cad/Cam Milling Diamond



Code	Name	Price	+
113474	CEREC compatble Compact Cylinder Pointed Bur #20 - 6pk -5855734 for Cerec In...	US\$ 153.00	<a href="#">+Info</a>
113475	CEREC compatble Compact Step Bur #12mm - 6 pk -6052265 for Cerec Inlab Compac...	US\$ 153.00	<a href="#">+Info</a>

## CEREC MCXL COMPATIBLE MILLING BURS - SHARK

### CEREC MCXL COMPATIBLE MILING BURS- SHARK

Details:

- \* Full Compatible Sinter Diamond Burs For Cerec Sirona
- \* Available For Ceramic And Zirconia

CEREC is a registered trademark of Sirona Dental Systems.

#UsddCerec



### Cylinder Burs

Code	Name	Price	+
112250	Cylinder Ptd Bur #12S for MCXL	US\$ 19.55	<a href="#">+Info</a>
116130	Cylinder Ptd Bur #20 for MCXL	US\$ 19.55	<a href="#">+Info</a>


### Step Burs

Code	Name	Price	+
116222	Step Bur #12 for MCXL ea	US\$ 19.55	<a href="#">+Info</a>
112251	Step Bur #12S for MCXL ea	US\$ 19.55	<a href="#">+Info</a>
116131	Step Bur #20 for MCXL ea	US\$ 19.55	<a href="#">+Info</a>

## CEREC X5 COMPATIBLE MILLING BURS - SHARK

### CEREC X5 COMPATIBLE MILLING BURS - SHARK

Dental milling bur compatible with Sirona CEREC inLab MC X5 for milling zirconia.  
Special coating to ensure long and accurate milling.



Code	Name	Price	+
123456	CEREC - MC X5 Bur 0.5mm Blue Ring Composite # SIC05MM Compatible# 6478171 ...	US\$ 42.00	+Info
123455	CEREC - MC X5 Bur 1.0mm Blue Ring Composite # SIC10MM Compatible# 6478163 ...	US\$ 42.00	+Info
123454	CEREC - MC X5 Bur 2.5mm Blue Ring Composite # SIC25MM Compatible# 6478155 ...	US\$ 42.00	+Info
123453	CEREC X5 Bur 0.5mm Red Ring P.MMA # SIP05MM Compatible# 6478114 [ fresas sirona ]	US\$ 42.00	+Info
123452	CEREC X5 Bur 1.0mm Red Zr # SIP10MM Compatible# 6478106 [ fresas sirona ]	US\$ 42.00	+Info
123451	CEREC - MC X5 Bur 2.5mm Red Ring P.MMA (Zirconio) # SIP25MM Compatible# 6478098 ...	US\$ 42.00	+Info
123450	CEREC - MC X5 Bur 0.5mm Yellow Ring ZrO2 (Zirconio) # SI05MM Compatible# 6478056...	US\$ 42.00	+Info
123449	CEREC - MC X5 Bur 1.0mm Yellow Ring ZrO2 (Zirconio) # SI10MM Compatible# 647804...	US\$ 42.00	+Info
123448	CEREC - MC X5 Bur 2.5mm Yellow Ring ZrO2 (Zirconio) # SI25MM Compatible# 647803...	US\$ 42.00	+Info

## DIAMOND COATED (DC) CAD/CAM MILLING BURS - ATTRITOR

### DIAMOND COATED (DC) CAD/CAM MILLING BURS - ATTRITOR

#### Features:

Compatible with AMANN GIRRBACH - VHF - WIELAND - ROLAND and ZIRKONZAHN



#### AMANN GIRRBACH

Code	Name	Price	+
124070	ATTRITOR - AMANN GIRRBACH Cad/Cam Bur DC round 0.6mm head - 3mm Shank - 47 mm le...	US\$ 68.00	+Info
124071	ATTRITOR - AMANN GIRRBACH Cad/Cam Bur DC round 1mm head - 3mm Shank - 47 mm leng...	US\$ 68.00	+Info
124072	ATTRITOR - AMANN GIRRBACH Cad/Cam Bur DC round 2-5mm head - 3mm Shank - 47 mm le...	US\$ 68.00	+Info

#### ROLAND

Code	Name	Price	+
124076	ATTRITOR - ROLAND DC Cad/Cam Bur Diamond round 0.6mm head - 4mm Shank - 50 mm le...	US\$ 58.00	+Info
124077	ATTRITOR - ROLAND DC Cad/Cam Bur Diamond round 1mm head - 4mm Shank - 50 mm leng...	US\$ 58.00	+Info
124078	ATTRITOR - ROLAND DC Cad/Cam Bur Diamond round 2mm head - 4mm Shank - 50 mm leng...	US\$ 58.00	+Info

#### VHF

Code	Name	Price	+
124063	ATTRITOR - VHF & YETI Cad/Cam Bur DC Diamond coated 0.6mm head - 3mm Shank - 35 ...	US\$ 68.00	+Info
124064	ATTRITOR - VHF & YETI Cad/Cam Bur DC Diamond coated 1mm head - 3mm Shank - 35 mm...	US\$ 68.00	+Info
124065	ATTRITOR - VHF & YETI Cad/Cam Bur DC Diamond coated 2mm head - 3mm Shank - 35 mm...	US\$ 68.00	+Info

## WIELAND

Code	Name	Price	+
124067	ATTRITOR - WIELAND Cad/Cam Bur Diamond coated 0.7mm head - 3mm Shank - 35 mm len...	US\$ 78.00	+Info
124068	ATTRITOR - WIELAND Cad/Cam Bur Diamond coated 1mm head - 3mm Shank - 35 mm lengh...	US\$ 78.00	+Info
124069	ATTRITOR - WIELAND Cad/Cam Bur Diamond coated 2-5mm head - 3mm Shank - 35 mm len...	US\$ 78.00	+Info

## ZIRKONZAHN

Code	Name	Price	+
124073	ATTRITOR - ZIRKONZAHN Cad/Cam Bur Diamond coated 0.6mm head - 3mm Shank - 47 mm ...	US\$ 68.00	+Info
124074	ATTRITOR - ZIRKONZAHN Cad/Cam Bur Diamond coated 1mm head - 3mm Shank - 47 mm le...	US\$ 68.00	+Info
124075	ATTRITOR - ZIRKONZAHN Cad/Cam Bur Diamond coated 2mm head - 3mm Shank - 47 mm le...	US\$ 68.00	+Info

## YETI - 417061 Digi MILL Cal. blank 10mm

US\$ 48.90

YETI - YETI digi MILL Calibration blank 10mm # 481-1001

## MANUAL ZIRCONIA BURS - ECCO

### MANUAL ZIRCONIA BURS - ECCO

- \* Compatible with Zirconcshan
- \* Comatible with Amann Girrbach Ceramill Manual
- \* Compatiblewith Ecco-Mill System
- \* 03 mm shank
- \* 47 mm Lengh
- \* Made in Germany



If you change bur brands we strongly recomend to buy the correspondent Reader the first time.

## KIT

Code	Name	Price	+
102214	ECCO - Zircon Drill bur/read Kit - 10 pieces - 5 readers / 5 burs. 3mm shank for...	US\$ 225.00	+Info

## ECCO

Code	Name	Price	+
102352	ECCO - Zircon Bur 1mm diam x 8 mm length - each - 3mm shank Total length 47mm ( ...	US\$ 35.70	+Info



102351	ECCO - Zircon Bur 1mm diam x 15 mm length - each - 3mm shank - Total length 47mm...	US\$ 35.70	<a href="#">+Info</a>
102350	ECCO - Zircon Bur 2mm diam - each - 3mm shank - Total length 47mm ( compatible w...	US\$ 35.70	<a href="#">+Info</a>
103404	ECCO - Zircon Bur 3mm diam - each - 3mm shank Total length 47mm ( compatible wit...	US\$ 35.70	<a href="#">+Info</a>
102349	ECCO - Zircon Bur 4mm diam - each - 3mm shank Total length 47mm ( compatible wit...	US\$ 35.70	<a href="#">+Info</a>

## READERS

Code	Name	Price	+
102355	ECCO - Zircon Reader for 1mm diam x 12 mm length ( 0.8mm diam ) each - 3mm shan...	US\$ 29.40	<a href="#">+Info</a>
103405	ECCO - Zircon Reader for 3mm diam -( 2.4mm diam ) each - 3mm shank Total length...	US\$ 29.40	<a href="#">+Info</a>
102353	ECCO - Zircon Reader for 4mm diam -( 3.2mm diam ) each - 3mm shank Total length...	US\$ 29.40	<a href="#">+Info</a>

## ROLAND DIAMOND CAD/CAM MILLING BURS - ATTRITOR

### ROLAND DIAMOND CAD/CAM MILLING BURS - ATTRITOR

#### Features:

A standard milling tool for finishing and adding detail to dental prosthetics, is designed for use on the Roland DWX-4W wet mill. It is approved for grinding glass ceramics and composite resins.

Available in 1mm - 2.5mm - 0.6mm



Code	Name	Price	+
124121	ATTRITOR - Roland Disilicate Diamond Cad Cam Bur 1mm for disilicate and feldspar	US\$ 48.00	<a href="#">+Info</a>
124122	ATTRITOR - Roland Disilicate Diamond Cad Cam Bur 0.6mm for disilicate and feld...	US\$ 48.00	<a href="#">+Info</a>
124123	ATTRITOR - Roland Disilicate Diamond Cad Cam Bur 2.5mm for disilicate and feld...	US\$ 48.00	<a href="#">+Info</a>

## ROLAND DLC D50 - D51 CAD/CAM MILLING BURS - SHARK

### ROLAND DLC D50-D51 CAD/CAM MILLING BURS - SHARK

#### Dental Milling Burs for its advanced D50 - D51 zirconia milling machine.

Each tool is specially engineered for precise performance, maximum durability and longer life.

#### Features:

- \* Available in 4 Different Sizes
- \* Carbide Milling Burs
- \* Compatible with roland existing Carbide Tools

Rolands milling burs are available in four different sizes allowing dental professionals to choose an ideal tool for every milling strategy.

Because they are fully compatible with Rolands existing carbide tools, users can switch to the new tools without any changes in the CAM software, making the upgrade a true "plug-n-play" solution.



Code	Name	Price	+
115583	ROLAND D50 - D51 Cad/Cam Coated Milling Burs 0.6mm - each - for ROLAND DWX-50 ...	US\$ 21.50	+Info
115578	ROLAND D50 - D51 Cad/Cam Coated Milling Burs 0.8mm - each - for ROLAND DWX-50 ...	US\$ 21.50	+Info
115572	ROLAND D50 - D51 Cad/Cam Coated Milling Burs 1mm - each - for ROLAND DWX-50 Ma...	US\$ 21.50	+Info
115566	ROLAND D50 D51 Cad/Cam DLC Milling Burs 2mm - each - for ROLAND DWX-50 Made...	US\$ 21.50	+Info

## VHF & YETI DLC CAD/CAM MILLING BURS - SHARK

### VHF & YETI DLC CAD/CAM COATED MILLING BUR - SHARK

#### Features:

- \* Ultra-precision.
- \* Hard, wear-resistant material.
- \* Last longer.
- \* Can be applied to milling the Zirconia disc, PMMA disc, Acrylic disc for CAD/CAM

**Quantity:** 1 pc

**Sizes:** 0.6 mm - 1mm -2mm

**Made in Germany**



### Set

Code	Name	Price	+
416999	VHF Cad/Cam DLC Coated Milling Bur - 40mm long 15 pcs set = 7 x 2mm - 5 1mm & 3...	US\$ 272.69	+Info

### 35mm 4 Axes

Code	Name	Price	+
115582	VHF & YETI Cad/Cam DLC Coated Milling Bur 0.6mm - 35mm long Each - Made in Ger...	US\$ 20.00	+Info
115575	VHF Cad/Cam DLC Coated Milling Bur 1mm - 35mm long Each - Made in Germany Z10...	US\$ 20.00	+Info
115569	VHF Cad/Cam DLC Coated Milling Bur 2mm - 35mm long Each - Made in Germany . Z20...	US\$ 20.00	+Info

### 40mm 5 Axes

Code	Name	Price	+
124307	VHF & YETI Cad/Cam DLC Coated Milling Bur 0.6mm - 40mm long Each - Made in Ge...	US\$ 20.00	+Info
124306	VHF & YETI Cad/Cam DLC Coated Milling Bur 1mm - 40mm long Each - Made in Germ...	US\$ 20.00	+Info
124305	VHF & YETI Cad/Cam DLC Coated Milling Bur 2mm - 40mm long Each - Made in German...	US\$ 20.00	+Info



## WIELAND DLC CAD/CAM MILLING BURS - SHARK

### WIELAND DLC CAD/CAM MILLING BURS - SHARK

WIELAND Cad Cam Coated Bur.

The wieland Cad/Cam Burs are sold individually and for use with the Ceramill Motion

**Features:**

- \* Milling Bur 2.0mm - 1.0mm - 0.6mm
- \* Super Accuracy.
- \* Hard ware-resistant material.
- \* Longer life.
- \* Can be applied to milling zirconia disc, wax block disc, acrylic disc CAD / CAM use.
- \* Quantity: 1 pc

**Made In Germany**



Code	Name	Price	+
115580	WIELAND Cad/Cam Coated Milling Bur 0.6mm - each -for Wieland & - Made in Germany	US\$ 20.00	+Info
115574	WIELAND Cad/Cam Coated Milling Bur 1mm - each - for Wieland & - Made in Germany	US\$ 20.00	+Info
115568	WIELAND Cad/Cam Coated Milling Bur 2.0mm - each - for Wieland & - Made in Germany	US\$ 20.00	+Info

## WIELAND DLC MINI CAD/CAM MILLING BURS - SHARK

### WIELAND DLC MINI CAD/CAM MILLING BURS - SHARK

WIELAND MINI Cad Cam Coated Bur. The Wieland Mini Cad / Cam Burs are sold individually and for use with the Ceramill Motion - Made in Germany

**Variable In:**

Milling Bur 2.5mm - 1.0mm - 0.6mm

**Features:**

- \* Super Accuracy.
- \* Hard Ware-Resistant Material.
- \* Longer Life.
- \* Can Be Applied To Milling Zirconia Disc, Wax Block Disc, Acrylic Disc CAD / CAM Use.

**Quantity: 1 Pc**

**Made In Germany**



Code	Name	Price	+
115579	WIELAND MINI Cad/Cam Coated Milling Bur 0.6mm - each - Made in Germany	US\$ 20.00	+Info
115573	WIELAND MINI Cad/Cam Coated Milling Bur 1mm - each - for Ceramill Motion - Made...	US\$ 20.00	+Info
115567	WIELAND MINI Cad/Cam Coated Milling Bur 2.5mm - each - for Ceramill Motion - M...	US\$ 20.00	+Info

**YETI & VHF GLASS CERAMIC MILLING BURs 4 & 5 Axes - ATTRITOR**



Code	Name	Price	+
417087	ATTRITOR - VHF & YETI Disilicate Sinter Diamond CadCam Bur 2.6mm FLAT END for d...	US\$ 48.00	<a href="#">+Info</a>
417088	ATTRITOR - VHF & YETI Disilicate Sinter Diamond CadCam Bur 2.4mm BALL END for d...	US\$ 48.00	<a href="#">+Info</a>
417081	ATTRITOR - VHF & YETI Disilicate Sinter Diamond CadCam Bur 1mm BALL END for dis...	US\$ 48.00	<a href="#">+Info</a>
417086	ATTRITOR - VHF & YETI Disilicate Sinter Diamond CadCam Bur 0.6mm BALL END for d...	US\$ 48.00	<a href="#">+Info</a>

**YETI & VHF GLASS CERAMIC MILLING BURs 4 & 5 Axes - YETI**

YETI - Milling Cutter Glass Ceramic CAD/CAM Burs

for 4 & 5 Axes machines



Code	Name	Price	+
417062	YETI - YETI digi MILL Milling Cutter M200-R4-32 (K/4) # 481-1121	US\$ 57.38	<a href="#">+Info</a>
417064	YETI - DIGI-Milling Cutters - Glass Ceramics (LiSi) CYLINDER BALL END ( Round )...	US\$ 62.70	<a href="#">+Info</a>
417063	YETI - DIGI-Milling Cutter Glass Ceramic ( LiSi) CONICAL - Radius 0.60 mm Len...	US\$ 62.70	<a href="#">+Info</a>
417091	YETI - DIGI-Milling Cutter Glass Ceramic ( LiSi) CYLINDER BALL END ( Flat ) - T...	US\$ 62.70	<a href="#">+Info</a>
417119	YETI - DIGI-Milling Cutters - Glass Ceramics (LiSi) CYLINDER FLAT END - 1.20 ...	US\$ 62.70	<a href="#">+Info</a>
417120	YETI - DIGI-Milling Cutters - Glass Ceramics (LiSi) CYLINDER FLAT END - 2.60 ...	US\$ 62.70	<a href="#">+Info</a>
417121	YETI - DIGI-Milling Cutter Glass Ceramic SET . 9 Pcs ( wet diamond sintered for...	US\$ 409.98	<a href="#">+Info</a>

## YETI DLC CAD/CAM MILLING BURS - YETI

YETI - DLC CAD/CAM MILLING BURS YETI



Code	Name	Price	+
417122	YETI - DIGI-Milling Cutter Zirconia PMMA Wax 0.60 mm Lenght: 40 mm # 481-0506...	US\$ 48.00	<a href="#">+Info</a>
417123	YETI - DIGI-Milling Cutter Zirconia PMMA Wax 2.0 mm Lenght: 40 mm # 481-0420 ...	US\$ 48.00	<a href="#">+Info</a>
417124	YETI - DIGI-Milling Cutter Zirconia PMMA Wax 1.0 mm Lenght: 40 mm # 481-0410 ...	US\$ 48.00	<a href="#">+Info</a>

### BESQUAL - 417130 Zircon Reader for 2mm di

US\$ 29.40

ECCO - Zircon Reader for 2mm diam x 12 mm length ( 1.6mm diam ) each - 3mm shank Total length 47mm ( compatible with Zircozhan- Amman girbach - Ceramill ] [ fresa para zirconia drill ] - Made in Germany

**BESQUAL - Zircon Reader for 2mm diam x 12 mm length**

**Burs**

**# 417130**

**Feature:**

- \* Compatible with Zircozhan
- \* Amman Girbach Ceramill Manual
- \* Manual System
- \* Every system with 3mm shank
- \* Made in Germany
- \* 3mm Shank
- \* 2mm Diameter
- \* 12mm length

## ZIRKONZAHN DLC M5 CAD/CAM MILLING BURS - SHARK

### ZIRKONZAHN DLC M5 CAD/CAM MILLING BURS - SHARK

M5 - TEC Milling - Bur 0.6mm

This M5 Zirkonzahn carbide burs are round end milling burs, and are designed to mill pre-sintered Zirconia with CAD/CAM Milling Unit M5 ( Precise Milling )

**Features:**

- \* Round End
- \* Carbide Milling Bur for M5
- \* Available in 3 Different Sizes 0.6mm - 1mm - 2mm

**ZIRKONZAHN is a trademark of ZIRKONZAHN GMBH**



Code	Name	Price	+
115581	ZirkonZahn M5 Cad/Cam DLC Milling Bur 0.5mm - each - Made in Germany	US\$ 20.00	<a href="#">+Info</a>
115576	ZirkonZahn M5 Cad/Cam DLC Milling Bur 1mm - each - Made in Germany	US\$ 20.00	<a href="#">+Info</a>
115570	ZirkonZahn M5 Cad/Cam DLC Milling Bur 2mm - each - Made in Germany	US\$ 20.00	<a href="#">+Info</a>

Zirconia abrasives

## Zirconia & Emax Grinders



## ABRASIVE DISCS - MPF

### ABRASIVE DISCS - MPF

**Features:**

Highly durable discs for ceramic and alloy restoration, without the risk of chipping or cracking the surface.

**Quantity:** 100pk

**Available in:** Small - Medium - Large



Code	Name	Price	+
123935	MPF BRUSH co - Aluminum Oxide Abrasive Disc small ( 22 x 0-2-0-25 x 1-8 ) 100 p...	US\$ 30.00	<a href="#">+Info</a>
123933	MPF BRUSH co - Aluminum Oxide Abrasive Disc medium (25 x 0-6 x 1-8 ) 100 pcs 123...	US\$ 30.00	<a href="#">+Info</a>
123934	MPF BRUSH co - Aluminum Oxide Abrasive Disc large ( 38 x 0-6 x 1-8 ) 100 pcs 123...	US\$ 44.00	<a href="#">+Info</a>

## COREMASTER - SHOFU

### COREMASTER - SHOFU

**Features:**

CoreMaster is designed for adjusting marginal areas and the thickness of aluminum oxide and zirconia copings.

- \* Diamond impregnated polishers eliminate microfractures
- \* Two grit sizes and three shapes
- \* Contouring and finishing
- \* No water cooling needed
- \* Use CoreMaster for: zirconia frameworks, alumina frameworks and porcelains



**US DENTAL DEPOT Inc is an Authorized SHOFU Dealer in the US.**

### Standard Grit



Code	Name	Price	+
106767	SHOFU - HP Coremaster KN7 HP - 2pk - ( #0139 ) ...	<b>US\$ 24.41</b>	<a href="#">+Info</a>
106773	SHOFU - HP Coremaster Standard Unmounted Wheel - 1pk - ( #0141 ) ...	<b>US\$ 34.42</b>	<a href="#">+Info</a>

## WHITE PEARLS ZIRCONIA GRINDER S - MPF Brush Co

### WHITE PEARLS ZIRCONIA GRINDERS - MPF Brush Co



**Features:**

Zirconia has become very popular in the dental industry in the past few years and its growing day by day.

Due to zirconias hardness and sensativity to heat, there was a great need for a rotary instrument that could cut and abrade this material with a minimum of heat generation as possible.

We have created a rotary instrument with a very high content of fine diamond particle which has been bonded with special binders that hold all the elements together, creating a superb Zirconia cutting material with exceptional life span, 2-3 times more life than any product on the market today. With our high content of diamond particle our product can cut and abrade with very little pressure needed from the handpiece, thus much less heat generation to the zirconia substructure.

Utilizing stainless iron mandrels for extra strength we guarantee a very high end looking and functioning product. We have laser printed onto our shanks for specific marking & branding for this exclusive product line.

### #UsddMpf

Code	Name	Price	+
123210	WHITEPEARL Invertd S Cone 6*8	US\$ 36.00	+Info
123209	WHITEPEARL Invertd Cone 12*6	US\$ 36.00	+Info
123208	WHITEPEARL Cone 3.5*11	US\$ 36.00	+Info
123207	WHITEPEARL Large Barrel 6.5*13	US\$ 36.00	+Info
123205	WHITEPEARL Thin Disc 22*0.7	US\$ 45.00	+Info
123206	WHITEPEARL Small Barrel 4.8*13	US\$ 18.00	+Info
123203	WHITEPEARL Small Disc 15*3.5	US\$ 45.00	+Info
123204	WHITEPEARL Disc 22*2	US\$ 45.00	+Info
123202	WHITEPEARL Large Disc 22*4.5	US\$ 52.00	+Info

## Z-GRINDERS ZIRCONIA GRINDERS - ATTRITOR

### Z-GRINDERS ZIRCONIA GRINDERS - ATTRITOR

Proud of this range of Sinter Diamond stones specially designed for working on all types of **Ceramics & Zirconia materials** .

With our unique combination of **Aluminum Oxide , Glass, Diamond Particle** and a range of ingredients that create this spectacular range of abrasive instruments we have created a product that is second to none in the market .

**Our unique product far exceeds the competitors in cutting ability and life span.**



## Z-GRINDERS

Diamond impregnated polishers for ceramics and zirconia materials. Superior in performance and much more economical due to its durability.

Code	Name	Price	+
112947	Z-Grindr Zirconia Wheel 15x3.5	US\$ 24.30	+Info
112228	Z-Grindr Zirconia Wheel 22x2	US\$ 24.30	+Info
100754	Z-Grindr Zirconia Wheel 22x4.5	US\$ 24.30	+Info
112948	Z-Grindr Zirconia Taper 3.5x11	US\$ 21.75	+Info
100756	Z-Grindr Zirconia Lg Barr 6x13	US\$ 17.91	+Info
100755	Z-Grindr Zirconia Med Cyl 5x13	US\$ 24.30	+Info



### Z-ZIRCON ZIRCONIA GRINDERS - AT TRITOR

#### Z-GRINDERS ZIRCONIA GRINDERS - ATTRITOR

Proud of this range of Sinter Diamond stones specially designed for working on all types of **Ceramics & Zirconia materials** .

With our unique combination of **Aluminum Oxide , Glass, Diamond Particle** and a range of ingredients that create this spectacular range of abrasive instruments we have created a product that is second to none in the market .

**Our unique product far exceeds the competitors in cutting ability and life span.**



### Z-ZIRCON

Code	Name	Price	+
115225	Zirconia Z-Zircon Med Inv Cone	US\$ 30.69	<a href="#">+Info</a>
114356	Zirconia Z-Zircon Disc 22x0.7	US\$ 28.69	<a href="#">+Info</a>
114177	Zirconia Z-Zircon Wheel 22x4.5	US\$ 26.57	<a href="#">+Info</a>
114176	Zirconia Z-Zircon Wheel 22x2	US\$ 30.69	<a href="#">+Info</a>
114175	Zirconia Z-ZirconWheel 15x3.5	US\$ 30.69	<a href="#">+Info</a>

### ZIRCOCERA ZIRCONIA GRINDERS - MPF Brush Co

#### ZIRCOCERA ZIRCONIA GRINDERS - MPF Brush Co

#### Features:

The GSS Series are a diamond impregnated rotary instrument of the highest quality, they boast a very high content of industrial diamond, bound together with a special ceramic binder. The GSS Series rotary instruments have extraordinary abrasive ability, they can cut through Zirconia with a minimal amount of pressure, and this alleviates the chances of micro fractures within the Zirconia substructure.

GSS Series have the ability to cut all types of ceramic material with great ease, and due to the fine grain. Industrial diamond and special binder, the ceramic surface is left with a smooth, clean finish.

Due to the composition of the diamond to binder ratio, the GSS Series boasts a very high service life, ergonomically it is by far the best solution for the ceramist today!



#UsddMpf

Code	Name	Price	+
123201	MPF BRUSH co ZIRCO CERA GSS Zirconia Grinder Wheel Inverted Small Cone 6mm x 8mm...	US\$ 34.00	+Info
123200	MPF BRUSH co ZIRCO CERA GSS Zirconia Grinder Wheel Inverted Cone 12mm x 6mm - 1p...	US\$ 34.00	+Info
123199	MPF BRUSH co ZIRCO CERA GSS Zirconia Grinder Wheel Cone 3.5mm x 11mm - 1pk- Sint...	US\$ 34.00	+Info
123198	MPF BRUSH co ZIRCO CERA GSS Zirconia Grinder Wheel Large Barrel 6.5mm x 13mm - 1...	US\$ 34.00	+Info
123197	MPF BRUSH co ZIRCO CERA GSS Zirconia Grinder Wheel Small Barrel 4.8mm x 13mm - 1...	US\$ 34.00	+Info
123196	MPF BRUSH co ZIRCO CERA GSS Zirconia Grinder Wheel Thin Disc 22mm x 0.7mm - 1pk-...	US\$ 44.00	+Info
123194	MPF BRUSH co ZIRCO CERA GSS Zirconia Grinder Wheel Small Disc 15mm x 3.5mm - 1pk...	US\$ 44.00	+Info
123195	MPF BRUSH co ZIRCO CERA GSS Zirconia Grinder Wheel Disc 22mm x 2mm - 1pk- Sinter...	US\$ 44.00	+Info
123192	MPF BRUSH co ZIRCO CERA GSS Zirconia Grinder Wheel Large Disc 22mm x 4.5mm - 1pk...	US\$ 44.00	+Info

**ATTRITOR - 416991 ZIRCONIA GRINDER KIT 8pcs** US\$ 115.00

ATTRITOR - Zirconia Grinder Kit Assorted - 8 pcs  
Zirconia & Emax Grinders  
# 416991

**Features:**

Proud of this range of Sinter Diamond stones specially designed for working on all types of Ceramics & Zirconia materials .

With our unique combination of Aluminum Oxide , Glass, Diamond Particle and a range of ingredients that create this spectacular range of abrasive instruments we have created a product that is second to none in the market . Our unique product far exceeds the competitors in cutting ability and life span.

**Quantity:** 8 pcs



**ZMAX ZIRCONIA GRINDERS - DIAMOND IMPREGNATED DISCS - MPF Brush Co**

ZMAX ZIRCONIA GRINDERS - DIAMOND IMPREGNATED DISCS - MPF Brush Co

**Features:**

Zirconia has become very popular in the dental industry



in the past few years and its growing day by day.

Due to zirconias hardness and sensitivity to heat, there was a great need for a rotary instrument that could cut and abrade this material with a minimum of heat generation as possible.

We have created a rotary instrument with a very high content of fine diamond particle which has been bonded with special binders that hold all the elements together, creating a superb Zirconia cutting material with exceptional life span, 2-3 times more life than any product on the market today. With our high content of diamond particle our product can cut and abrade with very little pressure needed from the handpiece, thus much less heat generation to the zirconia substructure.

Utilizing stainless iron mandrels for extra strength we guarantee a very high end looking and functioning product.

We have laser printed onto our shanks for specific marking & branding for this exclusive product line.

The ZMax come in 9 different shapes

### #UsddMpf

Code	Name	Price	+
123186	MPF BRUSH co ZMAX Zirconia Grinder Wheel Thin 22mm x 0.7mm - 1pk- Sintered Dia...	US\$ 48.00	+Info
123185	MPF BRUSH co ZMAX Zirconia Grinder Wheel 22mm x 2mm - 1pk- Sintered Diamond in...	US\$ 48.00	+Info
123184	MPF BRUSH co ZMAX Zirconia Grinder Wheel Small 15mm x 3.5mm - 1pk- Sintered Di...	US\$ 48.00	+Info
123183	MPF BRUSH co ZMAX Zirconia Grinder Wheel Large 22mm x 4.5mm - 1pk- Sintered Di...	US\$ 48.00	+Info
123190	MPF BRUSH co ZMAX Zirconia Grinder Wheel Inverted Cone 12mm x 6mm - 1pk- Sinter...	US\$ 44.00	+Info
123188	MPF BRUSH co ZMAX Zirconia Grinder Wheel Small Barrel 6.5mm x 13mm - 1pk- Sinte...	US\$ 44.00	+Info
123191	MPF BRUSH co ZMAX Zirconia Grinder Wheel Inverted Small Cone 6mm x 8mm - 1pk- S...	US\$ 44.00	+Info
123187	MPF BRUSH co ZMAX Zirconia Grinder Wheel Small Barrel 4.8mm x 13mm - 1pk- Sinte...	US\$ 44.00	+Info
123189	MPF BRUSH co ZMAX Zirconia Grinder Wheel Cone 3.5mm x 11mm - 1pk- Sintered Diam...	US\$ 44.00	+Info

### Zirconia abrasives

## Zirconia & Emax Polishers

### RAINBOW CERAMIC POLISHER - MPF Brush Co

#### RAINBOW CERAMIC POLISHER - MPF Brush Co

#### Features:

Premium diamond abrasive wheels for perfectly sealed ceramic or zirconia restoration.

For minimal reduction on zirconia marginal regions and all types of ceramic materials.

Available in color red, blue, gray and yellow



Code	Name	Price	+
123929	MPF BRUSH co MPF BRUSH co Rainbow- Ceramic Polisher- Coarse- Red 122-0001 - #USD...	US\$ 41.00	+Info
123930	MPF BRUSH co Rainbow- Ceramic Polisher- Medium- Blue 122-0002 - #USDDMPF BRUSH c...	US\$ 41.00	+Info
123931	MPF BRUSH co Rainbow- Ceramic Polisher- Grey- Fine 122-0003 - #USDDMPF BRUSH co ...	US\$ 41.00	+Info
123932	MPF BRUSH co Rainbow- Ceramic Polisher- Yellow- High Shine 122-0004 - #USDDMPF B...	US\$ 41.00	+Info

**ATTRITOR - 416558 ZIRCONIA POLISHER KIT US\$ 93.04**

**6pcs**

**ATTRITOR - Zirconia High Shine Bicolor Polisher Kit 6**

**pcs**

**Zirconia & Emax Polishers**

**# 416558**

**Features:**

Bicolor Diamond Impregnated Silicone Discs.

\* Excellent for minor surface scratches and prepares for a very high luster finish.

\* Excellent for Emax porcelain & Zirconia polishing and margin adjustments.

\* Grey Fine rubber disc a highly impregnated rubber disc for the removal of minor scratches on all types of ceramic materials with great ease, this range of rubber wheels has exceptional long

life properties with all ceramic materials

\* These wheels can be used to abrade and polish

ZIRCONIA AND ALL TYPES OF CERAMICS.

\* We advise that the RPM are maintained to our specifications ( 6000 rpm ) so as to obtain the optimum result and allows for a longer life span of the product.

**Description:** Zirconia polishers can replace EDENTA JOTA

**Color:** assorted

**Size:** 27 mm x 0.7 mm

**Contains Package:**

\* 1 Zirconia High Shine Yellow Super Fine

\* 1 Zirconia High Shine Grey Fine

\* 1 Zirconia High Shine Blue Medium

\* 1 Zirconia High Shine Red Coarse

\* 1 Dura Wheel - Super Fine White

\* 1 Dura Wheel Medium - BLUE



### ZIRCONIA POLISHERS POINT - ATTRITOR

#### ZIRCONIA POLISHERS POINT - ATTRITOR

- \* For fast shaping
- \* Scratch removal
- \* Leaves satin finish on all porcelains and final stage gloss finish on cercon
- \* Alumina and Zirconia copins



Code	Name	Price	+
100933	ATTRITOR - Dura Point Extra Fine - BLUE-Silicone Diamond point - for fast shapin...	US\$ 19.38	+Info
100932	ATTRITOR - Dura Point Fine PINK -Silicone Diamond point - for fast shaping -scr...	US\$ 19.38	+Info
100934	ATTRITOR - Dura Point Coarse - GREEN Silicone Diamond point - for fast shaping -...	US\$ 19.38	+Info

### ZIRCONIA POLISHERS WHEELS - ATTRITOR

#### ZIRCONIA POLISHERS WHEELS - ATTRITOR


- \* Diamond Impregnated Silicone Discs Two tone ( 22 x3mm )
- \* Diamond Impregnated Rubber Discs Mono Tone
- \* Diamond Impregnated Silicone Wheels



### HIGH SHINE SILICON WHEELS

These Mono tone color SILICONE WHEELS ( Square edge ) have a high content of refined diamond particle, which is bonded together with odorless Japanese Silicone material.

- \* EXTRA COARSE - BROWN
- \* COARSE - GRAY
- \* MEDIUM - BLUE
- \* FINE - PINK
- \* EXTRA FINE - WHITE



Code	Name	Price	+
100930	ATTR Dura Wheel- X-Coarse-Br	US\$ 21.75	+Info
100723	ATTR Dura Wheel -Coarse-Gray	US\$ 24.94	+Info
100722	ATTR IDura Wheel - Medium-Blue	US\$ 21.75	+Info
100763	ATTR Dura Wheel - Fine - Pink	US\$ 21.75	+Info
115563	ATTR Dura Wheel X Fine -White	US\$ 21.75	+Info

### HIGH SHINE SILICONE DISCS

These two tone color SILICONE DISCS have a very high content of refined diamond particle, which is bonded together with odorless Japanese Rubber material.



This product line is extremely popular with technicians around the world; the dimensions are 27mm x 0.7mm disc shape.

These wheels can be used to abrade and polish ZIRCONIA AND ALL TYPES OF CERAMICS  
We advise that the RPM are maintained to our specifications so as to obtain the optimum result and allows for a longer life span of the product.

Code	Name	Price	+
112956	High Shine Bicolor Yell S-Fine	US\$ 24.17	<a href="#">+Info</a>
112957	High Shine Bicolor Grey Fine	US\$ 24.17	<a href="#">+Info</a>
112958	High Shine Bicolor Blue Medium	US\$ 24.17	<a href="#">+Info</a>
112959	High Shine Bicolor Red Coarse	US\$ 24.17	<a href="#">+Info</a>

### POLISHING RUBBER DISCS

These Mono tone color rubber discs have a very high content of refined diamond particle, which is bonded together with odorless Japanese Rubber material.



This product line is extremely popular with technicians around the world; the dimensions are 27mm x 0.7mm disc shape.  
These wheels can be used to abrade and polish zirconia and all types of dental ceramic.

We advise that the RPM are maintained to our specifications so as to obtain the optimum result and allows for a longer life span of the product.

There abrasive grade range from Super coarse, coarse, medium, fine and superfine.

Code	Name	Price	+
112953	Zirconia Polish Lt Gray S-Fine	US\$ 21.10	<a href="#">+Info</a>
112952	Zirconia Polish Red Fine	US\$ 21.10	<a href="#">+Info</a>
112951	Zirconia Polish Blue Medium	US\$ 21.10	<a href="#">+Info</a>
112950	Zirconia Polish Green Coarse	US\$ 21.10	<a href="#">+Info</a>
112949	Zirconia Polish Gray Sup Corse	US\$ 21.10	<a href="#">+Info</a>



## Catalog Index / Indice

CAD/CAM Burs .....	1
Zirconia & Emax Grinders .....	12
Zirconia & Emax Polishers .....	17