

Autoclaves

Autoclave Units

AUTOCLAVE SERIE - ECCO

AUTOCLAVE SERIE - ECCO

Autoclave Euro B Standard

- * Adopting European with thermal vacuuming and drying function, it can be used for A type hollow instruments and solid instrument.
- * Thermodynamic vacuum guarantees elimination of cold air from cavity and makes sure of perfect saturation of water draining.
- * Three sterilizing water quantity selections are respectively different material, which embodies the individualize serve.
- * Automatic malfunction detecting system makes easy maintaining.



Autoclave New Euro B

This Autoclave is designed and made for places that need constant sterilization, such as clinics, medical departments, laboratories, etc. Being automatically controlled by a microprocessor and with a people-friendly interface, this sterilizer is safe, reliable and easy to operate. The different parameters and conditions in the operating process are displayed dynamically on a digital screen. Its automatic detection of malfunctions, protection of over-heating and excess pressure ensures the reliability of the sterilization. With a built-in condensed water collector dispelling the need to let off exhaust steam, this product is safe and environment-friendly.

Autoclave Emerald

The Emerald Autoclave sterilize surgical instruments, dressing, glassware, plastics. These Autoclave have: Dry function, Full Water Alarm, Auto Failure detection, Vacuum and two shelves.

Autoclave Manual

- * Compact design only 23"x15"x14"(60x39x36cm) 18 liter capacity
- * Temperature 251f or 269f (122c & 132c)
- * Sterilizing pressure 17.4-33.3ps (0.12mpa & 0.23mpa)
- * Time 25 or 6 minutes
- * Heater power 1000 watts/ 9 amps
- * Manual water supply
- * Overheating protection

AUTOCLAVE EURO B STANDARD

Code	Name	Price	+
106901	ECCO - Autoclave European B Standard TDK4220-18L - 220V1400 wats Size 445x590x40...	US\$ 2,309.77	+Info

AUTOCLAVE NEW EURO B

Code	Name	Price	+
106902	ECCO - Autoclave New Euro Class B TDK4330 -18L - 220V/1700 wats Size 480x590x455...	US\$ 2,519.75	+Info

AUTOCLAVE EMERALD

Code	Name	Price	+
112911	ECCO - Emerald Autoclave 23L - 110v Autoclave # EMERALD 23 - Dry function- Full...	US\$ 2,400.00	+Info

MANUAL

Code	Name	Price	+
115111	ECCO - Manual Autoclave 18Lts - 110v - 25cm*35cm chamber - 900w consumption - St...	US\$ 600.00	+Info
115493	ECCO - Manual Autoclave 18Lts - 220volts- 25cm*35cm chamber - 900w consumption - ...	US\$ 600.00	+Info

PRINTER FOR AUTOCLAVE

Code	Name	Price	+
106903	ECCO - Printer for Autoclave - New Euro Class B # 106902 Includes - Printer - Ad...	US\$ 230.98	+Info

ECCO - 106900 Autoclave Standard 18L

US\$ 1,100.00

ECCO - Autoclave Standard 18L
Autoclave Units
106900 - Mfg # TDK4210A 18L



Feature:

* Automatic malfunction detecting system makes easy maintaining.

Technical:

- * Voltage/Power: 110V/950W
- * Chamber Size: 170x320
- * Outside Size: 345x530x340MM
- * Net Weight: 29 kg

ECCO - 417320 Autoclave Standard 18L

US\$ 165.00

ECCO - Heating Element for Autoclave Standard 18L
Autoclave Units



www.usdentaldepot.com

SHARK - 123044 Ball Sterilizer 110v-220v

US\$ 120.00

SHARK Ball Sterilizer 110v-220v (sterilizador de bollilas de cuarzo)
SHARK -
Ball Sterilizer 110-220v
Autoclaves Units
123044



Features:

- * Best for metal tools
- * High temperature(250 Centigrade) to work
- * With Glass Beads in the inner pot to hold the metal tools being sterilized and to diffuse heat and keep hot
- * Convenient using, no pollution, energy saving
- CE Approval
- * UK standard plug

Technical Data

- * Volume of the Hole:
- * Voltage: 110V-220
- * Plug:UK plug
- * Power: 100W
- * Temperature: 250 Centigrade (Any temperature lower than 250 Centigrade will have no use for complete sterilization. To ensure the effect, we don't make it adjustable so it will automatically increase to 250 Centigrade as it is designed)



FREE SHIPPING
FOR ALL DOMESTIC ACCOUNTS!



US
DENTAL
DEPOT

+1 954 874 6325 
www.usdentaldepot.com 

STERIDENT - 225-201000 Steri-Dent Dry Heat # 200

US\$ 745.00

STERIDENT - Dry Heat Sterilizer Model #200 w/2 Trays 15.5x10.5x9in [Esterilizador]Steri-Dent 200 Dry Heat Sterilizer Certified #225-201000



The Steri-Dent 200 Dry Heat Sterilizer is a 2 tray model that kills bacteria by reaching 320 F (160 C) without moisture, and will not dull sharp instruments. Dry heat sterilization works without pressure, steam, or chemicals, which increases safety.

Warranty: 6 Months Parts and Labor. Extended Warranty available*.

The Steri-Dent 200 can also be used to dry wet packs or operate as a back up sterilizer in your practice - the perfect complement to an autoclave. The sterilizer is often used in dental offices or nail salons, and is not ideal for use as the single sterilizer in a large office due to the longer turnaround time

Features:

- * Operates on 115 V AC (single phase) or 230 VAC electrical power.
- * Maintenance-Free: Requires no routine cleaning.
- * Simple to Operate: Rapid-start heating element for fast operation.
- * Convenient: Utility trays can be sterilized with lids on and be stored without contamination.
- * Automatic timer with heat "hold" capability for continuous use in a busy practice.
- * High Quality: Stainless steel construction means no rust.
- * Value-Price: Low operating cost - no distilled water to buy.

Specifications:

- * Chamber Dimensions: 7.5" (L) x 11.25" (W) x 4.75" (H) (19 cm x 28.6 cm x 12 cm).
- * Chamber Volume: 6.08 L (1.6 Gallons).
- * Trays: 10.0" (L) x 6.25" (W).
- * Door: 11" x 4.5".
- * Dimensions: 15.5" x 10.5" x 9.5" (395 mm x 266 mm x 242 mm).
- * Power Current: 4 A.

* If you would like a custom extended warranty, please call or specify in the comments section during the checkout process. Additional cost may be applied. We will contact you before processing a custom order.



STERIDENT - 225-301000 Steri-Dent Dry Heat # 300

US\$ 925.00

STERIDENT- Dry Heat Sterilizer Model #300 W/3 Trays 19 1/2 X 14 1/2 X 11 [Esterilizador] **Steri-Dent 300 Dry Heat Sterilizer Certified #225-301000**



The Best Method for Sharp Instruments: Steri-Dent Dry Heat Sterilizers kill bacteria without moisture and will not dull sharp instruments. Plus, they can be used to dry wet packs or as a back-up sterilizer in your practice - the perfect complement to an autoclave.

Mercury Free Thermometer

Maintenance-Free: Unlike wet heat methods, Steri-Dent Dry Heat Sterilizers

require no routine cleaning.

Convenient: Utility trays can be sterilized with lids on and stored without contamination. Steri-Dent Sterilizers also have automatic timers with heat hold capability for continuous use in a busy practice.

Value-Priced: Low purchase cost and low operating cost no distilled water to buy! Steri-Dent Dry Heat Sterilizers Models 200 and 300 The Best Method for Sharp Instruments Steri-Dent Dry Heat Sterilizers kill bacteria without moisture and will not dull sharp instruments. Plus, they can be used to dry wet packs or as a back-up sterilizer in your practice the perfect compliment to an autoclave.

Safe and Simple to Operate Dry heat sterilization works without pressure, steam, or chemicals, which increases safety. Steri-Dent products feature a rapid-start heating element for fast operation. Thermal cutout protection prevents overheating.

Specifications

Sterilization Time: 60 minutes.

Unit Size: 19"W x 14.5"H x 11"D, 483 mm x 368 mm x 279 mm

Chamber Size: 13.5" x 6.75" x 8", 343 x 171 x 203 mm

Weight: 20 lbs, (9 Kg)

Shipping Weight/Dimensions: 45 lbs, 22" x 18" x 18"

Tray Size: (3 trays) 11.75"W x 7.25"D, 298 x 184 mm

Electrical: 220-240V, 750 watts

Sterilization: 320 Fahrenheit (160 Celsius) Temperature

Warranty: One year

High-Quality: Stainless steel construction means no rust ever!

Approvals

ETL FDA 510(k) K771070 UL 61010-1, IEC Safety Standard 61010-2-040

CAN/CSA C22.2 No. 61010-1



Catalog Index / Indice

Autoclave Units 1