

Instruments Endodontic



**ECCO- 786649 ENDODONTIC EXPLORER
DG-16**

US\$ 4.25



Ecco- Endodontic Explorer DG16 Double Side Instrument
#50-464 DG-16 Explorador Endodontico doble
**Ecco-
Endodontic Explorer DG16 Double Side Instrument
Instruments Endodontic
786649 - Mfg # 50-464**

Technical Specifications:

- Material:** Stainless Steel
- Rusting Prevention Procedure:** Passivated
- Ultrasonic Cleaned:** Yes
- Polished:** Yes
- Usage:** Left Hand or Right Hand
- Tests Performed:** Boil Test, Performance Test, Shape Test
- Packing:** Individually Packed
- QC Passed:** Yes

Common Uses:

Periodontal (gum) instruments are used in various gingivitis and periodontitis.

ENDODONTIC PLUGGERS - AMERICAN EAGLE
ENDODONTIC PLUGGERS - AMERICAN EAGLE

Features:

Endodontic Root Canal Pluggers are stainless steel instruments with a fine taper and blunt at the tip used for pressing or forcing a gutta percha cone into a root canal.

Handle material: Stainless Steel



| Code | Name | Price | + |
|--------|---|------------|-------|
| 119519 | AMERICAN EAGLE - Endo plugger 1-3 - # AEE1-3P [instruments] | US\$ 32.07 | +Info |
| 119520 | AMERICAN EAGLE - Endo plugger 3-5 - # AEE3-5P [instruments] | US\$ 32.07 | +Info |
| 119521 | AMERICAN EAGLE - Endo plugger 5-7 - # AEE5-7P [instruments] | US\$ 32.07 | +Info |
| 119522 | AMERICAN EAGLE - Endo plugger 9-11 - # AEE9-11P [instruments] | US\$ 32.07 | +Info |

ENDODONTIC SPREADERS - AMERICAN EAGLE

ENDODONTIC SPREADERS - AMERICAN EAGLE

Features:

Endodontic Spreaders are tapered metal instruments used to compress filling material against the sides of the canal to make room for additional cones during root canal procedures.

Description:

- * Material: Stainless Steel
- * Single-Ended
- * Standard



| Code | Name | Price | + |
|--------|---|-------------------|-----------------------|
| 119528 | AMERICAN EAGLE - Spreader d11t-25 - # AEED11T-25S [instruments] | US\$ 26.77 | +Info |
| 119527 | AMERICAN EAGLE - Spreader d11t-21 - # AEED11T-21S [instruments] | US\$ 26.77 | +Info |
| 119526 | AMERICAN EAGLE - Spreader d11-21 - # AEED11-21S [instruments] | US\$ 26.77 | +Info |

Instruments Endodontic

ENDO PLUGGER



AMERICAN EAGLE - 119525 Endo carver glick 1

US\$ 33.53

AMERICAN EAGLE - Endo Carver Glick 1
Instruments Endodontic
119525 - Mfg # AEEMG1CA



CARVER GLICK 1
AEEMG1CA
1/4" stainless steel handle

Features:

Endodontic Carver Glick # 1 is a double ended endodontic instrument

Handle Material: Stainless Steel

Size: 6.3 mm standard handle

FREE SHIPPING
FOR ALL DOMESTIC ACCOUNTS!



US
DENTAL
DEPOT

+1 954 874 6325 
www.usdentaldepot.com 

AMERICAN EAGLE - 119523 Endo excavator 31I US\$ 33.53

AMERICAN EAGLE - Endo Excavator 31I
Instruments Endodontic
119523 - Mfg # AEEM31LEXC



EXCAVATOR 31I
AE EM 31L EXC
1/4" stainless steel handle

Features:

Excavators are used for the removal of softened carious dentin.

Description:

- * Number of Ends: Double
- * Handle Material: Stainless Steel
- * Handle Type: Standard
- * Application: Endodontic, Restorative

Diameter: 1/4" - 6.3mm

AMERICAN EAGLE - 119524 Endo explorer dg-16 US\$ 24.51

AMERICAN EAGLE - Endo explorer dg-16 - # AEEMDG16
[instruments]AMERICAN EAGLE - Endo Explorer Dg-16
Instruments Endodontic
119524 - Mfg # AEEMDG16



EXPLORER DG 16
AE EM DG16
7/32" stainless steel handle

Features:

American Eagle Instruments Inc. Explorers are diagnostic instruments used in the detection of caries (decay) and calculus, exploration of pocket characteristics, furcations and restorations.

Description:

- * Handle Shape: Round
- * Number of Ends: Double
- * Handle Material: Stainless Steel
- * Handle Type: Standard
- * Autoclave Compatible

Diameter: 7/32" - 5.6 mm

Catalog Index / Indice

| | |
|------------------------------|---|
| Instruments Endodontic | 1 |
| ENDO PLUGGER | 2 |